



Wound Care

1-800-962-9888



COVIDIEN

Kendall™ AMD Antimicrobial Foam Dressings

Covidien Kendall™ AMD antimicrobial foam dressings have been developed to support your efforts in managing both moisture and bacterial balance. The dressing is an open-cell polyurethane foam and is impregnated with 0.5% Polyhexamethylene Biguanide (PHMB), a highly effective low-toxicity antiseptic. The dressings are effective against gram positive and gram negative bacteria. They are also effective against fungi and yeast. The Kendall AMD antimicrobial foam border dressings are designed to help prevent complications like friction and shear injuries.

Additional AMD antimicrobial products can be found on pages 9 and 10.

For more information, please visit: www.KendallAMD.com or www.KendallAMDFoam.com



Kendall AMD Antimicrobial Foam Dressings

- Contain PHMB (Polyhexamethylene Biguanide).
- Resist bacterial colonization within the dressing.
- Reduce bacterial penetration through the dressing.
- Provides protection against gram-, gram+ (including MRSA and VRE) and fungi/yeast microorganisms.
- Limits cross-contamination from and to the patient, clinician and environment.
- Are sterile and hypoallergenic.

Code	Description	Per Box	Ship Case
Kendall™ AMD Antimicrobial Foam Dressings			
55522AMD	2" x 2" (5 cm x 5 cm)	25	100
55535AMD	3.5" x 3" (8.8 cm x 7.5 cm) Fenestrated	10	50
55544AMD	4" x 4" (10 cm x 10 cm)	10	50
55548AMD	4" x 8" (10 cm x 20 cm)	10	50
55566AMD	6" x 6" (15 cm x 15 cm)	10	50
55588AMD	8" x 8" (20 cm x 20 cm)	10	50
55535PAMD	3.5" x 3" (8.8 cm x 7.5 cm) Fenestrated, with Topsheet	10	50
55544PAMD	4" x 4" (10 cm x 10 cm) with Topsheet	10	50
Kendall AMD Antimicrobial Foam Border Dressings			
55523BAMD	1.75" x 3.25" (4.4 cm x 8.3 cm) Pad: 1" x 1.75" (2.5 cm x 4.4 cm)	10	50
55544BAMD	3.5" x 3.5" (8.9 cm x 8.9 cm) Pad: 2" x 2" (5.0 cm x 5.0 cm)	10	50
55546BAMD	3.5" x 5.5" (8.9 cm x 14 cm) Pad: 2" x 4" (5.0 cm x 10 cm)	10	50
55566BAMD	5.5" x 5.5" (14 cm x 14 cm) Pad: 4" x 4" (10 cm x 10 cm)	10	50
55588BAMD	7.5" x 7.5" (19 cm x 19 cm) Pad: 6" x 6" (19 cm x 19 cm)	10	50
Kendall AMD Antimicrobial Foam Discs			
55511AMD	1" (2.54 cm) diameter, 4 mm hole	10	40
55512AMD	1" (2.54 cm) diameter, 7 mm hole	10	40

Do not use Covidien AMD antimicrobial dressings with Dakin's Solution™ cleanser or bleach solutions. These solutions deactivate the PHMB and will cause yellowing of the product.