

Medent Barr Guarantee

If You Are Not 100% Satisfied
With This BARR'S Branded
Product Return It for A 100%
Refund Used or Not

MEDENT
BARR
100%

Gloves Products



Product Catalogue

2015 | Edition 1



REACH FOR THE BEST
SUPERMAX

When only Medical Grade will do.



[Online Sales & Ordering](#)



HEAD OFFICE:
PH 04 527 7702 FAX 04 527 7499
PH OUTSIDE WELLINGTON: 0800 26 333 68
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)

About Medent Medical

After 30 years supplying the medical profession, Medent Medical remains 100% kiwi owned and 100% committed to delivering to you the best of medical, surgical, pharmaceutical and equipment supplies whenever you need them.

BISMOS - The innovative inventory management tool

We wanted to make it easier for you to order, store and manage your medical supplies. We came up with BISMOS. A complete stock management and ordering solution that takes the pain out of supply management.

With the help of BISMOS you can:

- ✓ Automatically generate your next order by scanning products in your storeroom
- ✓ Easily check stock levels
- ✓ Improve practice productivity
- ✓ Cut medical supply ordering in half
- ✓ Reduce out of stock items
- ✓ And much more

[Discover how BISMOS can work for you](#)

Providing you the best of medical, surgical & pharmaceutical supplies

With a constantly expanding inventory with thousands of products to choose from, including all of the major brands, when it comes to [medical supplies](#) think of Medent Medical as your 'go to' company. We have been supplying the medical profession with the best of medical, surgical, pharmaceutical and equipment supplies for more than 30 years. That is 30 years of working for you to simplify your ordering needs. You simply won't find a faster, more efficient way to manage your inventory.

Ask how we can help you:

- ✓ Optimise your ordering process
- ✓ Streamline inventory replacement
- ✓ Facilitate easy online ordering
- ✓ Enhance your productivity
- ✓ Improve ordering efficiency
- ✓ And save yourself both time and money

For more choice, better service and expert advice on your medical supplies, [contact us today](#)

BISMOS AUTOMATED ORDERING AND STORAGE SOLUTIONS



- ✓ Save your time spent ordering supplies
- ✓ Reduce order error rate
- ✓ Easy checking of what stock you have
- ✓ Flexible modular system




TABLE OF CONTENTS

CATEGORY: LATEX POWDER FREE NON STERILE GLOVES	1
Supermax Latex Powder Free Non Sterile Medical Grade Examination Glove.....	1
Ansell Latex Powder Free Micro-Touch Glove	2
Med-X Lite Latex Powder Free Textured Glove	2
HyTec Latex Powder Free Glove	2
CATEGORY: LATEX LOW POWDER NON STERILE GLOVES	3
Supermax Latex Low Powder Medical Grade Non Sterile Examination Glove.....	3
HyTec Clear Latex Low Powder Glove	4
CATEGORY: LATEX LOW POWDER STERILE GLOVES	5
Supermax Latex Low Powder Sterile Medical Grade Surgical Glove.....	5
Triflex Latex Powdered Sterile Surgical Gloves	6
Ansell Gammex Gloves Latex Powdered Surgical Sterile	6
CATEGORY: LATEX POWDER FREE STERILE GLOVES	7
BioGEL Latex Powder Free Sterile Gloves	7
Ansell Gammex Latex Powder Free Surgical Sterile Gloves	7
Ansell Nutex Dermashield Latex Powder Free Sterile Surgical Gloves	8
Intouch Latex Powder Free Sterile Gloves	8
Protexis Latex Powder Free Micro Sterile Gloves (Micro SMT)	9
Triumph Latex Powder Free Sterile Gloves	9
CATEGORY: LATEX FREE POWDER FREE NON STERILE GLOVES	14
Supermax Latex Free Powder Free Blue Nitrile Non Sterile Gloves.....	14
Elite Latex Free Powder Free Blue Nitrile Non Sterile Gloves	15
Kimclark Latex Free Powder Free Nitrile Non Sterile Exam Glove	15
CATEGORY: LATEX FREE POWDER FREE STERILE GLOVES	16
Protexis Latex Free Powder Free Sterile Surgical Gloves	16
Dermaprene Latex Free Powder Free Sterile Surgical Gloves	16
CATEGORY: VINYL CO-POLYMER POLYTHENE	17
Vinyl Low Powder Gloves	17
Polythene Blue Med-X	17
Vinyl Powder Free Gloves	17
Med-X Co-Polymer Polythene Glove	18
CATEGORY: OTHERS	18
Rubber Gloves General Purpose	18
Finger Cots	18
Glove Dispenser Single Barrs Flexi Perspex	19
Glove Dispenser Treble Barrs Flexi Perspex	19

CATEGORY: LATEX POWDER FREE NON STERILE GLOVES



SUPERMAX

Supermax Latex Powder Free Non Sterile Medical Grade Examination Glove



Features

- ✓ Ambidextrous
- ✓ Non Sterile Examination Gloves
- ✓ Natural Rubber Latex
- ✓ Low powder
- ✓ Low protein Contains 50µg/gm or less water extractable protein
- ✓ Beaded Cuff
- ✓ Micro Textured
- ✓ 100 gloves per Box
- ✓ Packed in ready to use dispenser box
- ✓ Easy on hands
- ✓ Disposable examination glove for medical, industrial & food industries
- ✓ Chlorinated, white glove



510k registration



Standard Malaysian Glove Certificate Approval 2021



ISO 11193-1 Listing 71106

- Australian TGAL 63625
- Directive 89/686/EEC Article 11B-PPE from SGS, UK
- Conformity Identification Seal from OCP, Brasil

Supermax Gloves also meet the following standards and testing requirements:

- Water tightness: ASTM D5151:99
- Dimension: ASTM D3767:01
- Tensile: ASTM D412:98 & ASTM D573:99
- Protein Testing: ASTM D5712:99
- Residual Powder Testing: ASTM D6124:01
- Biocompatibility Testing (Samples sent to PTC, USA for testing)
- Bioburden Testing (Samples sent to microbiological lab for testing)

Medent-Code	Mfr Code	Description	Unit of Sale
GLLPF1	SM08225	Gloves Supermax Latex Powder Free Medical Grade Extra Small	Box 100
GLLPF2	SM08226	Gloves Supermax Latex Powder Free Medical Grade Small	Box 100
GLLPF3	SM08227	Gloves Supermax Latex Powder Free Medical Grade Medium	Box 100
GLLPF4	SM08228	Gloves Supermax Latex Powder Free Medical Grade Large	Box 100
GLLPF5	SM08229	Gloves Supermax Latex Powder Free Medical Grade Extra Large	Box 100

[Online Sales & Ordering](#)



HEAD OFFICE:
PH 04 527 7702 FAX 04 527 7499
PH OUTSIDE WELLINGTON: 0800 26 333 68
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)

Ansell Latex Powder Free Micro-Touch Glove

Features and Benefits

- ✓ Micro-Touch examination gloves combine excellent tactility with high levels of comfort
- ✓ Designed for caring for frail skin while delivering a positive, non-aggressive grip
- ✓ Glove of choice for those wishing to double glove



Medent-Code	Mfr Code	Description	Unit of Sale
GLNP1	4560	Ansell SensiClean Powder Free Latex Examination Glove Small	Box 100
GLNP2	4562	Ansell SensiClean Powder Free Latex Examination Glove Medium	Box 100
GLNP3	4564	Ansell SensiClean Powder Free Latex Examination Glove Large	Box 100
GLANDE	4574	Ansell DermaClean Powder Free Latex Examination Glove Medium	Box 100

Med-X Lite Latex Powder Free Textured Glove

Features and Benefits

- ✓ General purpose, single use, form fitting, ambidextrous, good protection
- ✓ For glove wearers with sensitive skin



Medent-Code	Mfr Code	Description	Unit of Sale
GLLPFS	05081S	Med-X Lite Latex Powder Free Textured Glove Small	Box 100
GLMDXPFL	05036 L	Med-X Latex Powder Free Textured Gloves Large	Box 100

HyTec Latex Powder Free Glove

Features and Benefits

- ✓ Recommended for use in healthcare, laboratories, automotive and printing
- ✓ Latex offers ultimate flexibility, dexterity and sensitivity
- ✓ Close fitting protection, non-sterile and ambidextrous



Medent-Code	Mfr Code	Description	Unit of Sale
GLCL2	CLPFL	HyTec Latex Powder Free Large Glove	Box 100
GLCL3	NLPFM	HyTec Latex Powder Free Medium Glove	Box 100
GLCL4	NLPFS	HyTec Latex Powder Free Small Glove	Box 100
GLCL2A	NLPFXL	HyTec Latex Powder Free Extra Large Glove	Box 100
GLLPF1A	CLPFXS	HyTec Latex Powder Free Extra Small Glove	Box 100

CATEGORY: LATEX LOW POWDER NON STERILE GLOVES



SUPERMAX

Supermax Latex Low Powder Medical Grade Non Sterile Examination Glove



Features

- ✓ Ambidextrous
- ✓ Non Sterile Examination Gloves
- ✓ Natural Rubber Latex
- ✓ Low powder
- ✓ Low protein
- ✓ Beaded Cuff
- ✓ Fine touch surface
- ✓ Textured surface with fine grip
- ✓ 100 gloves per Box
- ✓ Packed in ready to use dispenser box



510k registration



Standard Malaysian Glove Certificate Approval 2021



ISO 11193-1 Listing 71106

- Australian TGA L 63625
- Directive 89/686/EEC Article 11B-PPE from SGS, UK
- Conformity Identification Seal from OCP, Brasil

Supermax Gloves also meet the following standards and testing requirements:

- Water tightness: ASTM D5151:99
- Dimension : ASTMD3767:01
- Tensile : ASTMD412:98 & ASTMD573:99
- Protein Testing : ASTMD5712:99
- Residual Powder Testing : ASTMD6124:01
- Biocompatibility Testing (Samplesent to PTC, USA for testing)
- Bioburden Testing (Samplesent to microbiological lab for testing)

Medent-Code	Mfr Code	Description	Unit of Sale
GLL2	SM28826	Gloves Supermax Latex Powdered Medical Grade Small	Box 100
GLL3	SM28827	Gloves Supermax Latex Powdered Medical Grade Medium	Box 100
GLL4	SM28828	Gloves Supermax Latex Powdered Medical Grade Large	Box 100

HyTec Clear Latex Low Powder Glove

Features and Benefits

- ✓ Recommended for use in healthcare, laboratories, automotive and printing
- ✓ Latex offers ultimate flexibility, dexterity and sensitivity
- ✓ Close fitting protection, non-sterile and ambidextrous



Medent-Code	Mfr Code	Description	Unit of Sale
GLCL1	CLM	HyTec Clear Latex Low Powder Glove Medium	Box 100
GLCL5	CLXS	HyTec Clear Latex Low Powder Glove Extra Small	Box 100



CATEGORY: LATEX LOW POWDER STERILE GLOVES



SUPERMAX

Supermax Latex Low Powder Sterile Medical Grade Surgical Glove



Features

- ✓ Sterile surgical glove
- ✓ Natural rubber latex
- ✓ Low powder
- ✓ Low protein for ease of donning
- ✓ Sterilized by Gamma Irradiation
- ✓ 100% inspected
- ✓ Smooth surface with beaded cuff
- ✓ 50 pairs of gloves per Box



510k registration



Standard Malaysian Glove Certificate Approval 2021



ISO 11193-1 Listing 71106

- Australian TGA L 63625
- Directive 89/686/EEC Article 11B-PPE from SGS, UK
- Conformity Identification Seal from OCP, Brasil

Supermax Gloves also meet the following standards and testing requirements:

- Water tightness : ASTM D5151:99
- Dimension : ASTM D3767:01
- Tensile : ASTM D412:98 & ASTM D573:99
- Protein Testing : ASTM D5712:99
- Residual Powder Testing : ASTM D6124:01
- Biocompatibility Testing (Samples sent to PTC, USA for testing)
- Bioburden Testing (Samples sent to microbiological lab for testing)

Medent-Code	Mfr Code	Description	Unit of Sale
GLS2	SS6.5P	Gloves Supermax 6.5 (Sterile) Surgical Latex Low Powder	Pair
GLS3	SS7.0P	Gloves Supermax 7.0 (Sterile) Surgical Latex Low Powder	Pair
GLS4	SS7.5P	Gloves Supermax 7.5 (Sterile) Surgical Latex Low Powder	Pair
GLS5	SS8.0P	Gloves Supermax 8.0 (Sterile) Surgical Latex Low Powder	Pair

[Online Sales & Ordering](#)




HEAD OFFICE:
PH 04 527 7702 FAX 04 527 7499
PH OUTSIDE WELLINGTON: 0800 26 333 68
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)

Triflex Latex Powdered Sterile Surgical Gloves

Features and Benefits

- ✓ Natural rubber latex formulation
- ✓ Cream color, powdered, standard thickness
- ✓ Smooth, anti-slip finish for a firm grip
- ✓ Extensive rinse cycle reduces latex allergens and proteins
- ✓ Interlocking locking cuff helps prevent roll down



Medent-Code	Mfr Code	Description	Unit of Sale
GLT0	2D7250I	Triflex® 5.5 Sterile Latex Powdered Surgical Gloves	Pair
GLT01	2D7250I	Triflex® 5.5 Sterile Latex Powdered Surgical Gloves	Box 40
GLT1	2D7251	Triflex® 6.0 Sterile Latex Powdered Surgical Gloves	Pair
GLT11	2D7251	Triflex® 6.0 Sterile Latex Powdered Surgical Gloves	Box 40
GLT2	2D7252I	Triflex® 6.5 Sterile Latex Powdered Surgical Gloves	Pair
GLT22	2D7252I	Triflex® 6.5 Sterile Latex Powdered Surgical Gloves	Box 40
GLT3	2D7253	Triflex® 7.0 Sterile Latex Powdered Surgical Gloves	Pair
GLT33	2D7253	Triflex® 7.0 Sterile Latex Powdered Surgical Gloves	Box 40
GLT4	2D7254	Triflex® 7.5 Sterile Latex Powdered Surgical Gloves	Pair
GLT44	2D7254	Triflex® 7.5 Sterile Latex Powdered Surgical Gloves	Box 40
GLT5	2D7255I	Triflex® 8.0 Sterile Latex Powdered Surgical Gloves	Pair
GLT55	2D7255I	Triflex® 8.0 Sterile Latex Powdered Surgical Gloves	Box 40
GLT6	2D7256I	Triflex® 8.5 Sterile Latex Powdered Surgical Gloves	Pair
GLT66	2D7256I	Triflex® 8.5 Sterile Latex Powdered Surgical Gloves	Box 40
GLT7	2D7257	Triflex® 9.0 Sterile Latex Powdered Surgical Gloves	Pair
GLT77	2D7257	Triflex® 9.0 Sterile Latex Powdered Surgical Gloves	Box 40



Ansell Gammex Gloves Latex Powdered Surgical Sterile

Features and Benefits

Ansell Gammex Latex Powdered Surgical Glove

- ✓ Powdered glove offering strength, fit and feel
- ✓ Exceptional comfort and reduced hand fatigue
- ✓ Soft and malleable
- ✓ Maximum tactile sensitivity in fingers and palm
- ✓ Finely controlled powder levels
- ✓ Low allergenicity glove, low residual chemical levels
- ✓ Highest barrier protection standards
- ✓ Excellent moisture control



Medent-Code	Mfr Code	Description	Unit of Sale
GLGAM0	4151N	Gloves Gammex Ansell # 6.0 Latex Powdered Surgical Sterile	Box 40
GLGAM1	4154N	Gloves Gammex Ansell # 7.5 Latex Powdered Surgical Sterile	Box 40
GLGAM2	4155N	Gloves Gammex Ansell # 8.0 Latex Powdered Surgical Sterile	Box 40
GLGAM3	4153N	Gloves Gammex Ansell # 7.0 Latex Powdered Surgical Sterile	Box 40
GLGAM4	4152N	Gloves Gammex Ansell # 6.5 Latex Powdered Surgical Sterile	Box 40



CATEGORY: LATEX POWDER FREE STERILE GLOVES

BioGEL Latex Powder Free Sterile Gloves

Features and Benefits

- ✓ Excellent sterile, powder free natural rubber latex glove
- ✓ Great barrier protection in general procedures
- ✓ Recommended for all general procedures
- ✓ Unique Biogel coating on the inner surface, makes donning easier even with damp hands
- ✓ Soothes the skin to help prevent moisture loss
- ✓ Beaded cuff for extra security
- ✓ Every glove is air inflation tested, and visually inspected for quality and user confidence
- ✓ Powder free to eliminate starch powder related complications
- ✓ Low dermatitis potential.

Medent-Code	Mfr Code	Description	Unit of Sale
GLBG65	GEN	BioGEL 6.5 Latex Powder Free Sterile Gloves	Box 50
GLBG75	GEN465	BioGEL 7.5 Latex Powder Free Sterile Gloves	Box 50
GLBG8	GEN466	BioGEL 8.0 Powder Free Latex Sterile Gloves	Box 50



Ansell Gammex Latex Powder Free Surgical Sterile Gloves

Features and Benefits

- ✓ Provide excellent donning and comfort through the use of a hydrophilic/hydrophobic polymer lining allowing for both wet and dry donning
- ✓ Wide thumb ball area combined with Ansell's soft latex formulation ensures greater comfort and less hand fatigue
- ✓ Excellent latex moisture control maintains shape and fit with hand
- ✓ Thinner and shaped finger pads for improved finger sensitivity
- ✓ Tapered straight cuff for improved gown grip
- ✓ Sticky band for improved gown sleeve adherence
- ✓ Cuff printing of brand, size and L and R for easy identification



Gammex® PF



Powder Free Latex Surgical Gloves
The evolution in powder free surgical gloves



Ansell

Medent-Code	Mfr Code	Description	Unit of Sale
GLGAM1PF	4154PF	Gammex Ansell 7.5 Latex Powder Free Surgical Sterile Glove	Box 40
GLGAM2PF	4155PF	Gammex Ansell Latex Powder Free Surgical Sterile 8.0 Glove	Box 40

Ansell Nutex Dermashield Latex Powder Free Sterile Surgical Gloves

Features and Benefits

- ✓ Nutex DermaShield provides enhanced donning ability through its polymer lining
- ✓ DermaShield finish enables you to don with a damp hand, such as intra-operatively or dry hand
- ✓ Nutex DermaShield specific grip characteristics enable you to positively grip instruments
- ✓ Smooth upper surface allows for effective manipulation of delicate items such as suture material
- ✓ Cream coloured natural rubber latex to reduce glare from operating lights when using magnification instruments



Powder Free Sterile Surgical Gloves
The first choice in polymer coated powder free gloves

Ansell

Medent-Code	Mfr Code	Description	Unit of Sale
GLN1	4454D	Ansell Nutex 7.5 DermaShield Latex Powder Free Sterile Surgical Glove	Box 40
GLN3	4452D	Ansell Nutex 6.5 DermaShield Latex Powder Free Sterile Surgical Glove	Box 40

Intouch Latex Powder Free Sterile Gloves

Features and Benefits

INTOUCH Surgical Powder-free Latex gloves are an alternative choice over powdered gloves as these do not contain airborne powder particles. Gloves are polymer coated to promote ease of dry hand donning. They also have a stronger grip on the working surface




Medent-Code	Mfr Code	Description	Unit of Sale
GLPINT2	W05ST6.0	INTOUCH 6.0 Latex Powder Free Sterile Surgical Gloves	Pair
GLPINT3	W05ST6.5	INTOUCH 6.5 Latex Powder Free Sterile Surgical Gloves	Pair
GLPINT4	W05ST7	INTOUCH 7.0 Latex Powder Free Sterile Surgical Gloves	Pair
GLPINT5	W05ST7.5	INTOUCH 7.5 Latex Powder Free Sterile Surgical Gloves	Pair
GLPINT6	W05ST8	INTOUCH 8.0 Latex Powder Free Sterile Surgical Gloves	Pair
GLPINT7	W05ST8.5	INTOUCH 8.5 Latex Powder Free Sterile Surgical Gloves	Pair

- ✓ Individual Sterile Pouches
- ✓ Left and Right
- ✓ Surgical Sterile Grade
- ✓ Sealed Pouches
- ✓ Polymer Coated

Protexis Latex Powder Free Micro Sterile Gloves (Micro SMT)

Features and Benefits

- ✓ Natural rubber latex formulation
- ✓ Powder-free, helps to reduce powder-related skin irritation and post-operative complications associated with granulomas
- ✓ Exclusive glove mold with independent thumb design that allows for an anatomical fit and more natural movement in the fingers, thumb and palm
- ✓ Thinner (20 percent) than Protexis™ Latex Surgical Gloves for enhanced tactile sensitivity
- ✓ Triple-dip coating technology integrates nitrile and latex for effortless donning and glove strength
- ✓ Excellent choice as an outer glove when double-gloving
- ✓ Proprietary interlocking, beaded cuff design reduces roll-down



Medent-Code	Mfr Code	Description	Unit of Sale
GLP0	2D72NT65	Protexis™ 6.5 Latex Powder-Free Micro Sterile Surgical Gloves with Nitrile Coating	Box 50
GLP0E	2D72NT65	Protexis™ 6.5 Latex Powder-Free Micro Sterile Surgical Gloves with Nitrile Coating	Pair
GLP4	2D72NT70	Protexis™ 7.0 Latex Powder-Free Micro Sterile Surgical Gloves with Nitrile Coating	Box 50
GLP5	2D72NT75X	Protexis™ 7.5 Latex Powder-Free Micro Sterile Surgical Gloves with Nitrile Coating	Box 50
GLP6	2D72NT80	Protexis™ 8.0 Latex Powder-Free Micro Sterile Surgical Gloves with Nitrile Coating	Box 50
GLP8	2D72NT85X	Protexis™ 8.5 Latex Powder-Free Micro Sterile Surgical Gloves with Nitrile Coating	Box 50



Triumph Latex Powder Free Sterile Gloves

Features and Benefits

- ✓ Designed for neurosurgery, cardiovascular surgery, plastic surgery, and ophthalmology procedures
- ✓ Can be used for double-gloving with Triumph Green
- ✓ Can be used as an overglove when used for double-gloving purposes
- ✓ Low protein natural rubber latex material
- ✓ 8 mils thickness
- ✓ Smooth texture grip
- ✓ White colour
- ✓ Synthetic polymer coating donning agent

Medent-Code	Mfr Code	Description	Unit of Sale
GLP3	MSG2365	Triumph 6.5 Micro Latex Powder-Free Sterile Surgical Gloves	Pair
GLP7	MSG2370	Triumph 7.0 Micro Latex Powder-Free Sterile Surgical Gloves	Pair



CATEGORY: LATEX FREE POWDER FREE NON STERILE GLOVES



SUPERMAX

Supermax Latex Free Powder Free Blue Nitrile Non Sterile Gloves



Features

- ✓ Ambidextrous
- ✓ Non Sterile Exam Gloves
- ✓ 6-Mil Nitrile Synthetic
- ✓ Online Powder Free
- ✓ No NRL Protein
- ✓ High tensile to meet toughest demands
- ✓ Powder Free and 100% Latex free
- ✓ Beaded Cuff
- ✓ Finger Textured
- ✓ Blue Color
- ✓ 100 gloves per Box



510k registration



Standard Malaysian Glove
Certificate Approval 2021



ISO 11193-1 Listing
71106

- Australian TGA L 63625
- Directive 89/686/EEC Article 11B-PPE from SGS, UK
- Conformity Identification Seal from OCP, Brasil

Supermax Gloves also meet the following standards and testing requirements:

- Water tightness : ASTM D5151:99
- Dimension : ASTM D3767:01
- Tensile : ASTM D412:98 & ASTM D573:99
- Protein Testing : ASTM D5712:99
- Residual Powder Testing : ASTM D6124:01
- Biocompatibility Testing (Samples sent to PTC, USA for testing)
- Bioburden Testing (Samples sent to microbiological lab for testing)

Medent-Code	Mfr Code	Description	Unit of Sale
GLNPFS2	SM98896	Gloves Supermax Blue Nitrile Latex Free Powder Free Small	Box 200
GLNPFS3	SM98897	Gloves Supermax Blue Nitrile Latex Free Powder Free Medium	Box 200
GLNPFS4	SM98898	Gloves Supermax Blue Nitrile Latex Free Powder Free Large	Box 200

[Online Sales & Ordering](#)



HEAD OFFICE:
PH 04 527 7702 FAX 04 527 7499
PH OUTSIDE WELLINGTON: 0800 26 333 68
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)

Elite Latex Free Powder Free Blue Nitrile Non Sterile Gloves

Features and Benefits

- ✓ Non sterile, ambidextrous, single use
- ✓ Sky Blue with beaded cuff
- ✓ Latex protein free
- ✓ Textured fingers
- ✓ Thicker glove 4mil (0.10mm)
- ✓ Smooth finish



Medent-Code	Mfr Code	Description	Unit of Sale
GLNPFE2	C_GNPFEPS	Gloves Elite Latex Free Powder Free 4mil Nitrile Small	Box 100
GLNPFE3	C_GNPFEPM	Gloves Elite Latex Free Powder Free 4mil Nitrile Medium	Box 100
GLNPFE4	C_GNPFEP L	Gloves Elite Latex Free Powder Free 4mil Nitrile Large	Box 100
GLNPFE5	C_GNPFEPXL	Gloves Elite Latex Free Powder Free 4mil Nitrile Extra Large	Box 100

Kimclark Latex Free Powder Free Nitrile Non Sterile Exam Glove

Features and Benefits

- ✓ Non sterile, ambidextrous, single use
- ✓ Sky Blue with beaded cuff
- ✓ Latex protein free
- ✓ Textured fingers
- ✓ Thicker glove (0.10mm)
- ✓ Smooth finish





Medent-Code	Mfr Code	Description	Unit of Sale
GLNPF2	13940	Kimclark Latex Free Powder Free Nitrile Exam Glove Small	Box 200
GLNPF3	13941	Kimclark Latex Free Powder Free Nitrile Exam Glove Medium	Box 200
GLNPF4	13942	Kimclark Latex Free Powder Free Nitrile Exam Glove Large	Box 200
GLNPF5	13943	Kimclark Latex Free Powder Free Nitrile Exam Glove Extra Large	Box 170

CATEGORY: LATEX FREE POWDER FREE STERILE GLOVES

Protexis Latex Free Powder Free Sterile Surgical Gloves

Features and Benefits

- ✓ Best-selling unique synthetic polyisoprene surgical glove
- ✓ Chemical structure identical to natural rubber latex, without the proteins
- ✓ Offers healthcare workers an alternative solution that performs like latex in terms of comfort, fit and feel
- ✓ Powder-free, helps to reduce powder-related skin irritation and post-operative complications associated with granulomas
- ✓ Proprietary interlocking, beaded cuff design reduces roll-down
- ✓ Exclusive glove mold with independent thumb design that allows for an anatomical fit and more natural movement in the fingers, thumb and palm



Medent-Code	Mfr Code	Description	Unit of Sale
GLDU1	2D72PT70X	Gloves Protexis PI # 7.0 Surgical Latex Free	Box 50
GLDU2	2D72PT75X	Gloves Protexis PI # 7.5 Surgical Latex Free	Box 50
GLDU3	2D72PT65X	Gloves Protexis PI # 6.5 Surgical Latex Free	Box 50
GLDUM1	2D73PM65	Gloves Protexis PI Micro #6.5 Surgical Latex Free	Box 50
GLDUM2	2D73PM70	Gloves Protexis PI Micro #7.0 Surgical Latex Free	Box 50
GLDUM3	2D73PM80	Gloves Protexis PI Micro #8.0 Surgical Latex Free	Box 50

Dermaprene Latex Free Powder Free Sterile Surgical Gloves

- ✓ Gammex powder free neoprene surgical gloves
- ✓ Straight, thinner cuff to improve softness and comfort
- ✓ Wide 'sticky band' positioned at the cuff end to secure gown and prevent roll down
- ✓ Accelerator free (no chemical residuals)
- ✓ Puncture resistant, offering strength, softness, comfort and grip
- ✓ Provides secure handling of instruments in wet and dry conditions
- ✓ Thinner fingers and palm for increased tactile sensitivity
- ✓ Green colour for instant identification as a synthetic glove, reducing confusion and risk to both staff and patients with sensitive skin



Ansell

Medent-Code	Mfr Code	Description	Unit of Sale
GLD6	4212N	Gloves Dermaprene Sterile Latex Free Powder Free 6.5	Box 25
GLD8	4213N	Gloves Dermaprene Sterile Latex Free Powder Free 7.0	Box 25
GLD7	4214N	Gloves Dermaprene Sterile Latex Free Powder Free 7.5	Box 25
GLD9	4216N	Gloves Dermaprene Sterile Latex Free Powder Free 8.5	Box 25



Ansell

Ansell

CATEGORY: VINYL CO-POLYMER POLYTHENE

Vinyl Low Powder Gloves

Features and Benefits

- ✓ Non sterile, ambidextrous, single use, lightly powdered
- ✓ With beaded cuff, latex protein free, smooth finish

Medent-Code	Mfr Code	Description	Unit of Sale
GLCV1	-	Gloves Clear Vinyl Low Powder Large	Box 100
GLCV2	-	Gloves Clear Vinyl Low Powder Medium	Box 100
GLCV3	-	Gloves Clear Vinyl Powderfree Medium	Box 100
GLCV4	-	Gloves Clear Vinyl Powderfree Small	Box 100
GLCV5	-	Gloves Clear Vinyl Low Powder Small	Box 100
GLCV6	-	Gloves Clear Vinyl Low Powder Small	Box 100
GLCV7	-	Gloves Clear Vinyl Powderfree Large	Box 100
GLV3	05019L	Gloves Clear Vinyl Powdered Large	Box 100
GLBV1	-	Gloves Blue Vinyl Low Powder Medium	Box 100
GLBV4	-	Gloves Blue Vinyl Low Powder Large	Box 100



Polythene Blue Med-X

Features and Benefits

- ✓ X Light Weight Polythene Wrist Length
- ✓ Ambidextrous, lightweight protection, disposable

Medent-Code	Mfr Code	Description	Unit of Sale
GLPB2	05017M	FABRI-CELL FAB-X Gloves Polythene Blue Med-X Medium Blue	Pack 100



Vinyl Powder Free Gloves

Features and Benefits

- ✓ Non sterile, ambidextrous, single use, powder free
- ✓ With beaded cuff, latex protein free, smooth finish, thicker glove 4mil (0.10mm)

Medent-Code	Mfr Code	Description	Unit of Sale
GLCV3	-	Gloves Clear Vinyl Powder Free Medium	Box 100
GLCV4	-	Gloves Clear Vinyl Powder Free Small	Box 100
GLCV7	-	Gloves Clear Vinyl Powder Free Large	Box 100
GLBVPF2	C-GVPFBS	Gloves Blue Vinyl Powder Free Small	Box 100
GLBVPF3	C-GVPFBM	Gloves Blue Vinyl Powder Free Medium	Box 100
GLBVPF4	C-GVPFBL	Gloves Blue Vinyl Powder Free Large	Box 100



Med-X Co-Polymer Polythene Glove

Features and Benefits

- ✓ General purpose, single use
- ✓ Ambidextrous, lightweight protection



Medent-Code	Mfr Code	Description	Unit of Sale
GLCO2	01025M-	Med-X Co-Polymer Polythene Ster Glove (Medium)	Box 100
GLCO4	05006L	Med-X Co-Polymer Polythene Glove (Large)	Box 100
GLCO5	5006	Med-X Co-Polymer Polythene Glove (Medium)	Box 100

CATEGORY: OTHERS

Rubber Gloves General Purpose

Non sterile kitchen gloves

Medent-Code	Mfr Code	Description	Unit of Sale
GLRBL	05024-	Rubber Gloves General Purpose Blue Large	Pack 12
GLRYM	05022 (MED)	Rubber Gloves General Purpose Yellow Medium	Pack 12
GLRYL	05022 (LGE)	Rubber Gloves General Purpose Yellow Large	Pack 12



Finger Cots

Features and Benefits

- ✓ Designed to cover just the fingers and finger tips with a layer of latex material
- ✓ Useful for wound protection or industrial uses, but are not intended for exam use
- ✓ Cover and protect bandaged finger wounds, allowing employees with minor wounds to return to work while maintaining sanitary conditions
- ✓ Can be used for work with sensitive equipment, such as clock parts or computer hardware, to prevent fingerprints or body oil transfer
- ✓ Easy to grip and put on



Medent-Code	Mfr Code	Description	Unit of Sale
GLFCOTL	196L	Finger Cots Large Natural White Latex	Box 100

Glove Dispenser Single Barrs Flexi Perspex

Barr's plastic portable and single glove dispenser

Medent-Code	Mfr Code	Description	Unit of Sale
GLDISP	GLDISP	Glove Dispenser Single Barrs Flexi Perspex	Each



Glove Dispenser Treble Barrs Flexi Perspex

Barr's white plastic portable and treble glove dispenser

Medent-Code	Mfr Code	Description	Unit of Sale
GLDISP3	10336	Glove Dispenser Treble Barrs Flexi Perspex	Each



To Order Contact Us at:

Medent Medical

Unit 4, 110 Whakatiki St
PO Box 40356
Upper Hutt
New Zealand

Free Phone: 0800 263 3368

Phone: 64 4 527 7702

Free Fax: 0800 527 7499

Fax: 64 4 5277499

Email Enquires: enquiry@medent.co.nz

Finding a product and ordering it – you choose an option that suits you.
Either talk with one of our team or order right here online!

- Order from our sales representative when they visit you.
To arrange a visit please phone us on 0800 263 3368
- You can phone us on 0800 263 3368
or Fax us on 0800 527 7499

All products and more available at www.medent.co.nz

Medent Barr Guarantee

If You Are Not 100% Satisfied With This BARR'S Branded Product Return It for A 100% Refund Used or Not

MEDENT BARR 100%



BARR'S
The Practitioners First Choice




CardinalHealth
Essential to care™



Kimberly-Clark
Trusted Clinical Solutions*




MEDLINE
Sterile*
Triumph



Ansell

REACH FOR THE BEST
SUPERMAX

When only Medical Grade will do.




iNtouch



HYTEC



FINE TOUCH
DISPOSABLES



EliteHealthcare

[Online Sales & Ordering](#)



medent
MEDICAL

HEAD OFFICE:
PH 04 527 7702 FAX 04 527 7499
PH OUTSIDE WELLINGTON: 0800 26 333 68
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)