

Medent Barr Guarantee

If You Are Not 100% Satisfied
With This BARR'S Branded
Product Return It for A 100%
Refund Used or Not

MEDENT
BARR
100%

Disability Products



Product Catalogue

2015 | Edition 15



[Online Sales & Ordering](#)



HEAD OFFICE:
PH 04 527 7702 FAX 04 527 7499
PH OUTSIDE WELLINGTON: 0800 26 333 68
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)

About Medent Medical

After 30 years supplying the medical profession, Medent Medical remains 100% kiwi owned and 100% committed to delivering to you the best of medical, surgical, pharmaceutical and equipment supplies whenever you need them.

BISMOS - The innovative inventory management tool

We wanted to make it easier for you to order, store and manage your medical supplies. We came up with BISMOS. A complete stock management and ordering solution that takes the pain out of supply management.

With the help of BISMOS you can:

- ✓ Automatically generate your next order by scanning products in your storeroom
- ✓ Easily check stock levels
- ✓ Improve practice productivity
- ✓ Cut medical supply ordering in half
- ✓ Reduce out of stock items
- ✓ And much more

[Discover how BISMOS can work for you](#)

Providing you the best of medical, surgical & pharmaceutical supplies

With a constantly expanding inventory with thousands of products to choose from, including all of the major brands, when it comes to [medical supplies](#) think of Medent Medical as your 'go to' company. We have been supplying the medical profession with the best of medical, surgical, pharmaceutical and equipment supplies for more than 30 years. That is 30 years of working for you to simplify your ordering needs. You simply won't find a faster, more efficient way to manage your inventory.

Ask how we can help you:

- ✓ Optimise your ordering process
- ✓ Streamline inventory replacement
- ✓ Facilitate easy online ordering
- ✓ Enhance your productivity
- ✓ Improve ordering efficiency
- ✓ And save yourself both time and money

For more choice, better service and expert advice on your medical supplies, [contact us today](#)

BISMOS AUTOMATED ORDERING AND STORAGE SOLUTIONS



- ✓ Save your time spent ordering supplies
- ✓ Reduce order error rate
- ✓ Easy checking of what stock you have
- ✓ Flexible modular system




TABLE OF CONTENTS

CATEGORY: BANDAGES-CAST	1
Smith & Nephew Delta Cast Elite	1
Smith & Nephew Delta Cast Soft	2
Smith & Nephew Delta Dry Cast	2
Smith & Nephew Dynacast Prelude	3
Smith & Nephew Gypsona* BP Plaster of Paris Bandages	4
Prime Cast Polyester	4
Primenet Polyester Stockinette	5
3M Scotchcast Plus Casting Tape	5
3M Scotchcast Softcast Tape	5
CATEGORY: ORTHOPAEDIC PADDING	5
Smith & Nephew Soffban Orthopaedic Padding	5
Smith & Nephew Soffban Procell Cast Liner (Goretex)	5
Primepad Polyester Undercast Padding (Soffban Alt)	6
Reston Foam Roll Adhesive Light Support	6
CATEGORY: CASTING PROTECTOR	7
Aquashield Waterproof Cast Protector	7
Protector Shin/Calf Wrap	7
CATEGORY: BRACE/SPLINT/SUPPORT	8
Walker Moon Boot	8
Ankle Stabiliser	8
AFO-Deluxe-Allcare	8
Procure Beach Boy Cast Boot	9
BSN Cast Shoe Adult	9
Cast Boot Shoe	9
Procure Wrap Patella Knee	9
Elastic Knee Support/Stabiliser	10
Arm Sling - Immobiliser With Waist Strap	10
Sling Arm With Strap Support	10
Mesh Arm Sling	10
Universal Wrist Splint	11
OPPO Wrist Splint With Strap	11
Neoprene Wrist Stabiliser	11
Thumb brace	11
Finger Splint Zimmer (Aluminium)	12
Finger Stall Leather	12
Finger Splints Stax	12

Thumb Stall Leather 12

Collar Neck Hard Adjustable Adult..... 13

Cardboard Splint 13

Cervical Collars 13

CATEGORY:MOBILITY 14

Crutches Adult Double Adjustable Forearm (Pair)..... 14

Stoppers For Crutches (Pair) 14

CATEGORY: BANDAGES-CAST

Smith & Nephew Delta Cast Elite

Description

Smith & Nephew Delta Cast Elite is a polyester cast tape designed with ease of application and patient comfort in mind. The patented fibreglass-free casting material combines strength and rigidity with soft edges for enhanced patient comfort.

Provides rigid casting options for primary and secondary immobilisation, where strength and radiolucency are required.

Patented polyester substrate provides exceptional moulding, and the unique resin system ensures maximum lamination and end laydown for long lasting casts.



Features and Benefits

- ✓ Smooth cast with soft edges
- ✓ Fibreglass free
- ✓ Clean, easy to apply and remove
- ✓ Highly conformable
- ✓ Strong yet lightweight
- ✓ Durable to reduce cast breakdown
- ✓ Smooth finish
- ✓ Consistent set times
- ✓ Enhanced radiograph viewing

Medent Code	Mfr Code	Description	Unit of Sale
BACAD1	72279-01	Smith & Nephew Delta Cast Elite White 7.5cm x 3.6m	Box 10
BACAD4	72278-02	Smith & Nephew Delta Cast Elite Black 10cm x 3.6m	Box 10
BACAD12	72276-00	Smith & Nephew Delta Cast Elite Blue 5.0cm x 3.6m	Box 10
BACAD13	72276-01	Smith & Nephew Delta Cast Elite Blue 7.5cm x 3.6m	Box 10
BACAD14	72275-01	Smith & Nephew Delta Cast Elite Red 7.5cm x 3.6m	Box 10
BACAD16	72275-02	Smith & Nephew Delta Cast Elite Red 10.0cm x 3.6m	Box 10
BACAD18	72277-01	Smith & Nephew Delta Cast Elite Purple 7.5cm x 3.6m	Box 10
BACAD19	72275-00	Smith & Nephew Delta Cast Elite Red 5cm x 3.6m	Box 10
BACAD20	72278-00	Smith & Nephew Delta Cast Elite Black 5cm x 3.6m	Box 10
BACAD21	72278-01	Smith & Nephew Delta Cast Elite Black 7.5cm x 3.6m	Box 10
BACAD22	72276-02	Smith & Nephew Delta Cast Elite Blue 10cm x 3.6m	Box 10
BACAD23	72277-00	Smith & Nephew Delta Cast Elite Purple 5cm x 3.6m	Box 10
BACAD24	72277-02	Smith & Nephew Delta Cast Elite Purple 10cm x 3.6m	Box 10
BACAD25	72279-03	Smith & Nephew Delta Cast Elite White 5cm x 3.6m	Box 10
BACAD26	72279-02	Smith & Nephew Delta Cast Elite White 10cm x 3.6m	Box 10
BACAD27	72279-03	Smith & Nephew Delta Cast Elite White 12.5cm x 3.6m	Box 10
BACAD29	73455-00	Smith & Nephew Delta Cast Elite Pink 5cm x 3.6m	Box 10
BACAD30	73455-01	Smith & Nephew Delta Cast Elite Pink 7.5cm x 3.6m	Box 10
BACAD31	73455-02	Smith & Nephew Delta Cast Elite Pink 10cm x 3.6m	Box 10
BACAD32	73454-00	Smith & Nephew Delta Cast Elite Green 5cm x 3.6m	Box 10
BACAD33	73454-01	Smith & Nephew Delta Cast Elite Green 7.5cm x 3.6m	Box 10
BACAD34	73454-02	Smith & Nephew Delta Cast Elite Green 10cm x 3.6m	Box 10
BACAD35	72273-22	Smith & Nephew Delta Cast Elite Camo 5cm x 3.6m	Box 10
BACAD36	72273-06	Smith & Nephew Delta Cast Elite Camo 7.5cm x 3.6m	Box 10
BACAD37	72273-17	Smith & Nephew Delta Cast Elite Pastel 5cm x 3.6m	Box 10
BACAD38	72273-07	Smith & Nephew Delta Cast Elite Pastel 7.5cm x 3.6m	Box 10



Smith & Nephew Delta Cast Soft

Description

Smith & Nephew Delta-Cast Soft is fibreglass-free and was designed to provide variable support to fractures and soft tissue injuries. Its semi-rigidity allows for more flexibility and mobility, thus helping to reduce immobilisation related problems such as atrophy or joint stiffness. As a result, the overall healing process tends to be shortened.



Features and Benefits

- ✓ Earlier patient mobility as cast will allow a limited range of movement
- ✓ Improved circulation, faster healing, shorter rehabilitation period
- ✓ Help avoid or minimise muscle atrophy and other immobilisation related issues
- ✓ Lightweight, breathable material offers high wearing comfort
- ✓ Smooth surface and edges for a snag-free finish
- ✓ Available in a wide range of attractive colours

Medent-Code	Mfr Code	Description	Unit of Sale
BACAD8	72708-08	Smith & Nephew Delta Cast Soft Red 5.0cm x 3.6m	Box 10
BACAD9	72708-05	Smith & Nephew Delta Cast Soft Blue 5.0cm x 3.6m	Box 10
BACAD10	72708-18	Smith & Nephew Delta Cast Soft Red 7.5.0cm x 3.6m	Box 10
BACAD11	72708-25	Smith & Nephew Delta Cast Soft Blue 10.0cm x 3.6m	Box 10
BACAD17	72708-02	Smith & Nephew Delta Cast Soft White 2.5cm x 3.6m	Box 10
BACAD28	72708-12	Smith & Nephew Delta Cast Soft White 5cm x 3.6m	Box 10



Smith & Nephew Delta Dry Cast

Description

The unique patented technology in Smith & Nephew Delta-Dry® has created a water resistant cast padding which allows easy application by the clinician and affords the patient all the lifestyle benefits of a cast they can get wet.

Delta-Dry® is a water resistant padding while Delta-Dry® Softliner is a water resistant padded stockinette and can be used without additional padding.



Features and Benefits

- ✓ Better long term water resistance
- ✓ Offers more under cast protection
- ✓ Eliminates the need for stockinette
- ✓ Allows patients to shower, swim and bathe like normal
- ✓ No special tools are needed for application or removal
- ✓ Conformable gives patient increased comfort and protection



Medent-Code	Mfr Code	Description	Unit of Sale
BACDRY	73443-04	Smith & Nephew Delta Dry Cast Padding 7.5cm x 2.4m	Box 12
BACDRY1	73443-03	Smith & Nephew Delta Dry Cast Padding 5.0cm x 2.4m	Box 12
BACDRYS1	72779-04	Smith & Nephew Delta Dry Cast Softliner 10.0cm x 10m	Each
BACDRYS5	72779-02	Smith & Nephew Delta Dry Cast Softliner 5.0cm x 10m	Each
BACDRYS7	72779-03	Smith & Nephew Delta Dry Cast Softliner 7.5cm x 10m	Each

Smith & Nephew Dynacast Prelude

Description

Smith & Nephew Dynacast Prelude is a synthetic splinting system incorporating polyurethane coated fibre glass covered by polypropylene padding, to offer easy application for the support and immobilisation of fresh fractures and soft tissue injuries. It is presented in a foil covered roll format for flexibility in use and application.

Features and Benefits

- ✓ A single interlock knit layer of polyurethane resin-impregnated fibreglass
- ✓ The latest innovation in splinting
- ✓ Smoother
- ✓ More conformable
- ✓ Lamination free



Medent-Code	Mfr Code	Description	Unit of Sale
BACAD2	71415-00	Smith & Nephew Dynacast Prelude 10.0cm x 4.6m	Each
BACAD5	71421-00	Smith & Nephew Dynacast Prelude 7.5cm x 4.6m	Each
BACAD6	71416-00	Smith & Nephew Dynacast Prelude 12.5cm x 4.6m	Each
BACAD7	71419-00	Smith & Nephew Dynacast Prelude 20.0cm x 2.3m	Each
BACAD15	71417-00	Smith & Nephew Dynacast Prelude 15.0cm x 4.6m	Each

Smith & Nephew Gypsona* BP Plaster of Paris Bandages

Description

Gypsona bandages combine fast setting, high quality plaster of Paris with the conformable and stable interlocked weave of leno cloth. The bandage demonstrates minimum plaster loss while the rich creamy texture allows easy handling.



Medent-Code	Mfr Code	Description	Unit of Sale
BACAG1	87060527	Smith & Nephew Gypsona Cast White 5.0cm x 2.75m	Box 12
BACAG1E	87060527	Smith & Nephew Gypsona Cast White 5.0cm x 2.75m	Each
BACAG2	87061535	Smith & Nephew Gypsona Cast White 15.0cm x 3.5m	Box 12
BACAG2E	87061535	Smith & Nephew Gypsona Cast White 15.0cm x 3.5m	Each
BACAG3	87067535	Smith & Nephew Gypsona Cast White 7.5cm x 3.5m	Box 24
BACAG3E	87067535	Smith & Nephew Gypsona Cast White 7.5cm x 3.5m	Each
BACAG4	87061035	Smith & Nephew Gypsona Cast White 10.0cm x 3.5m	Box 24
BACAG4E	87061035	Smith & Nephew Gypsona Cast White 10.0cm x 3.5m	Each
BACAGS1	87004915	Smith & Nephew Gypsona Slab 15.0cm x 20m	Each
BACAGS2	87004913	Smith & Nephew Gypsona Slab 10.0cm x 20m	Each
BACAGS3	87004916	Smith & Nephew Gypsona Slab 20.0cm x 20m	Each



Features and Benefits

- ✓ Rapid immobilisation of fractures
- ✓ Leno cloth for good conformability
- ✓ Round cores for improved bandage handling
- ✓ Low plaster loss
- ✓ Moisture proof packaging

Prime Cast Polyester

Description

PrimeCast™ is a synthetic casting tape impregnated with polyurethane resin which is activated by water before use. Quality of PrimeCast™ is well proven through years of successful sales worldwide and it is certainly one of the top quality products currently available in the market.

Features and Benefits

- ✓ A single interlock knit layer of polyurethane resin-impregnated fibreglass
- ✓ The latest innovation in splinting, smoother, more conformable
- ✓ Lamination free



Medent-Code	Mfr Code	Description	Unit of Sale
BACAP102B	CP04-BLACK	Prime Cast Polyester Black 10.2cm x 3.6m	Each
BACAP102BL	CPO4-BLUE	Prime Cast Polyester Blue 10.2cm x 3.6m	Each
BACAP102PI	CP04-PINK	Prime Cast Polyester Pink 10.2cm x 3.6m	Each
BACAP102P	CP04-PURPLE	Prime Cast Polyester Purple 10.2cm x 3.6m	Each
BACAP127BL	CP05-BLUE	Prime Cast Polyester Blue 12.7cm x 3.6m	Each
BACAP76B	CP03-BLACK	Prime Cast Polyester Black 7.6cm x 3.6m	Each
BACAP76PI	CP03-PINK	Prime Cast Polyester Pink 7.6cm x 3.6m	Each
BACAP76P	CP03-PURPLE	Prime Cast Polyester Purple 7.6cm x 3.6m	Each
BACAP76BL	CP03-BLUE	Prime Cast Polyester Blue 7.6cm x 3.6m	Each

Primenet Polyester Stockinette

Features and Benefits

- ✓ Cast stockinette made of skin friendly, finest quality synthetic material
- ✓ Use under synthetic cast, Plaster of Paris bandage and protection on amputation stump
- ✓ Protect patients' skin from maceration and irritation during casting
- ✓ Packed in easy to use dispenser carton



Medent-Code	Mfr Code	Description	Unit of Sale
BACAPN1	NET05P	Primenet 05 Polyester 5cm x 25m Stockinette	Each
BACAPN2	NET07P	Primenet 7.5 Polyester 7.5cm x 25m Stockinette	Each
BACAPN3	NET10P	Primenet 10 Polyester 10cm x 25m Stockinette	Each

3M Scotchcast Plus Casting Tape

Features and Benefits

- ✓ Lightweight for easy mobility, strong and durable casting tape
- ✓ Water-activated resin lets you control working time
- ✓ Allows for weight bearing in 20 minutes
- ✓ Combines the benefits of a fiberglass casting tape with the handling ease of plaster
- ✓ Stretches to conform well to contours
- ✓ Porous weave allows the skin to breathe
- ✓ Contains a lubricating agent in the resin
- ✓ Available in a variety of colors and sizes to maximize patient preferences



Medent-Code*	Mfr Code	Description	Unit of Sale
BACAS100B	82004A	Scotchcast Plus 3M 100mm x 3.6m Black	Box 10
BACAS125R	82005R	Scotchcast Plus 3M 125mm x 3.6m Red	Box 10
BACAS125W	82005	Scotchcast Plus 3M 125mm x 3.6m White	Each
BACAS50R	82002W	Scotchcast Plus 3M 50mm x 3.6m Red	Box 10
BACAS75G	82003G	Scotchcast Plus 3M 75mm x 3.6m Green	Box 10

*For illustration purposes only, other colours available for each size.

3M Scotchcast Softcast Tape

Features and Benefits

- ✓ Fiberglass casting tape that does not set to a completely rigid state, but remains flexible and soft even when it is completely cured
- ✓ Consists of a knitted fiberglass fabric impregnated with a unique polyurethane resin
- ✓ Finished bandage will retain its shape while remaining resilient



Medent-Code*	Mfr Code	Description	Unit of Sale
BACASC100	82104	Scotchcast Softcast 3M 100mm x 3.6m White	Pack 10
BACASC50	82102	Scotchcast Softcast 3M 50mm x 3.6m White	Pack 10
BACASC50B	82102B	Scotchcast Softcast 3M 100mm x 3.6m Blue	Pack 10
BACASC75	82103	Scotchcast Softcast 3M 75mm x 3.6m White	Pack 10

*For illustration purposes only, other colours available for each size.

CATEGORY: ORTHOPAEDIC PADDING

Smith & Nephew Soffban Orthopaedic Padding

Description

Synthetic, protective orthopaedic padding, that blends excellent conformability and good cushioning, with a water-shedding ability that will keep the patients skin dry.

The padding tears quickly and easily around awkward contours, minimising patient trauma and enhancing ease of application for the clinical user

Soffban Natural tears and teases apart easily, aiding rapid application and reducing the risk of stressing the fracture. It is also autoclavable, making it suitable for use in theatre or sterile procedures.

Made entirely from rayon fleece, Soffban# Natural provides soft, comfortable cushioning and protection for the patient, yet still conforms well for easy shaping to all body contours and joints. Ideal for use under Gypsona plaster of Paris bandages, the SoffbanNatural fleece absorbs sweat and exudate, reducing the risk of skin maceration and irritation.



Medent-Code	Mfr Code	Description	Unit of Sale
BASOF1	71464-01	Smith & Nephew Soffban 5CM x 2.7M	Box 12
BASOF2	71486-05	Smith & Nephew Soffban 10CM x 2.7M	Box 12
BASOF3	71467-01	Smith & Nephew Soffban 7.5CM x 2.7M	Box 12
BASOF4	71465-01	Smith & Nephew Soffban 15CM x 2.7M	Box 12
BASOF2N	71469-00	Smith & Nephew Soffban Natural 10CM x 2.7M	Box 12
BASOF3N	71473-00	Smith & Nephew Soffban Natural 7.5CM x 2.7M	Box 12
BASOF5	1090040	Smith & Nephew Soffban Sterile 7.5CM x 2.7	Each

Features and Benefits

- ✓ Highly conformable
- ✓ Easily tearable
- ✓ High loft
- ✓ Feathers out
- ✓ Low moisture

Smith & Nephew Soffban Procell Cast Liner (Goretex)

Description

Waterproof, washable, highly breathable quilted undercast liner developed to maintain patient quality of life. The innovative cast liner design will permit patients to bathe, shower or swim in a synthetic cast, allowing normal everyday activities to continue unconstrained during the treatment period.

Features and Benefits

- ✓ Reduced inconvenience
- ✓ Improved hygiene
- ✓ Accelerated recovery
- ✓ Simplified removal



Medent-Code	Mfr Code	Description	Unit of Sale
BACASO1	67002181	Smith & Nephew Soffban Procel Waterproof Cast Liner (Goretex) 10CM x 1.7M	Box 12
BACASO2	NCL2	Smith & Nephew Soffban Procel Waterproof Cast Liner (Goretex) 5CM x 1.7M	Box 12



Primepad Polyester Undercast Padding (Soffban Alt)

Features and Benefits

- ✓ Undercast padding for cushioning under synthetic cast and
- ✓ Plaster of Paris bandage
- ✓ Made of 100% finest quality polyester
- ✓ Keeps the patient's skin dry and comfortable during long-term wear, easy to tear
- ✓ Smooth surface and good air permeability

Medent-Code	Mfr Code	Description	Unit of Sale
BACAP01	PAD05P	Primepad Polyester 5cm x 2.7m	Box 12
BACAP02	PAD07P	Primepad Polyester 7.5cm x 2.7m	Box 12
BACAP03	PAD10P	Primepad Polyester 10cm x 2.7m	Box 12
BACAP04	PAD15P	Primepad Polyester 15cm x 2.7m	Box 12



Reston Foam Roll Adhesive Light Support

Features and Benefits

Self-adhering open cell foam designed to protect skin from splints, casts and prostheses. Product is not intended to be used as a pressure reducing device. It can also be used as a padding for beds and chairs.

Medent-Code	Mfr Code	Description	Unit of Sale
DPFRRO	1563L	Reston Foam Roll Adhesive 10cmx5m Light Support	Box 5



CATEGORY:CASTING PROTECTOR

Aquashield Waterproof Cast Protector

Features and Benefits

- ✓ Keeps a cast or bandage dry which is critical to the safe healing process
- ✓ Made from heavy-gauge polyurethane, a material renowned for its strength and elasticity
- ✓ Features a pre-formed opening that can be stretched to 5 times its size yet return to within 3% of its original dimension
- ✓ Slides easily over a cast or bandage to create a comfortable & watertight seal



Medent-Code	Mfr Code	Description	Unit of Sale
DPAQ1	A17	Aquashield Waterproof Cast Protector (Sml Half Arm)	Each
DPAQ2	A21	Aquashield Waterproof Protector (Reg Half Arm)	Each
DPAQ3	A26	Aquashield Waterproof Protector (Sml Full Arm)	Each
DPAQ4	A31	Aquashield Waterproof Protector (Reg Full Arm)	Each
DPAQ5	L20	Aquashield Waterproof Protector (Sml Halfleg)	Each
DPAQ6	L25	Aquashield Waterproof Protector (Reg Half Leg)	Each
DPAQ7	L33	Aquashield Waterproof Protector (Sml Full Leg)	Each
DPAQ8	L44	Aquashield Waterproof Protector (Reg Full Leg)	Each



Protector Shin/Calf Wrap

Features and Benefits

- ✓ Calf Leg Protector Velcro Fastening. Inner Fabric Brushed Polyester for Skin Protection Core Fabric Quilted Tetron Fibre Outer Fabric Micro Fleece for Softness and Protection
- ✓ This item will give higher protection than the tube type due o the construction of the product.
- ✓ Can be used as a Heavy Duty Arm Protector.
- ✓ Can be washed on a wool cycle, and gentle cycle dried



Medent-Code	Mfr Code	Description	Unit of Sale
DPRPRO2	CLV10	Protector Shin/Calf Wrap Small	Pair
DPRPRO3	CLV20	Protector Shin/Calf Wrap Medium	Pair
DPRPRO4	CLV30	Protector Shin/Calf Wrap Large	Pair

CATEGORY:BRACE/SPLINT/SUPPORT

Walker Moon Boot

Features and Benefits

- ✓ Contoured strut design, ensuring a better fit
- ✓ Unique shock absorbing sole reduces impact from the heel strike to improve comfort during ambulation
- ✓ The rocker bottom is specially engineered to be low and wide to help promote a natural gait and increase stability

Medent-Code	Mfr Code	Description	Unit of Sale
DPBLMBSS	ALL20103830	Walker Moon Boot 11 Inch High Small	Each
DPBLMBMS	ALL20103840	Walker Moon Boot 11 Inch High Medium	Each
DPBLMBLS	ALL20103850	Walker Moon Boot 11 Inch High Large	Each
DPBLMBS	ALL21101630	Walker Moon Boot 17 Inch High Small	Each
DPBLMBM	ALL21101640	Walker Moon Boot 17 Inch High Medium	Each
DPBLMBL	ALL21101650	Walker Moon Boot 17 Inch High Large	Each
DPBLMBLXS	ALL20103820	Walker Moon Boot 11 Inch High Xsmall	Each
DPBLMBXS	ALL21101620	Walker Moon Boot 17 Inch High Xsmall	Each



Ankle Stabiliser

Features and Benefits

- ✓ Helps support & stabilise the ankle
- ✓ Lace up design for a snug and customised fit
- ✓ Non-elastic straps lock the heel to resist inversion & eversion
- ✓ Elastic closure strap 'locks in' tied laces & stabilising straps
- ✓ Fits left or right ankle
- ✓ For mild to moderate ankle sprain, ankle instability, recurrent ankle sprain



Medent-Code	Mfr Code	Description	Unit of Sale
DPASL	AOA1830	Ankle Stabiliser (Large)	Each
DPASM	AOA1830	Ankle Stabiliser (Medium)	Each
DPASS	AOA1830	Ankle Stabiliser (Small)	Each

AFO-Deluxe-Allcare

Features and Benefits

- ✓ Designed to support flacid drop foot
- ✓ Semi rigid calf section for toe clearance and support
- ✓ Thin and flexible foot plate for easy trimming a full length sole



Medent-Code	Mfr Code	Description	Unit of Sale
DPADALL	ALLAFODLX-LL	Afo - Deluxe (L-Large)	Each
DPADARL	ALLAFODLX-RL	Afo - Deluxe (R-Large)	Each
DPADALM	ALLAFODLX-LM	Afo - Deluxe (L-Medium)	Each
DPADARM	ALLAFODLX-RM	Afo - Deluxe (R-Medium)	Each
DPADALS	ALLAFODLX-LS	Afo - Deluxe (L-Small)	Each
DPADARS	ALLAFODLX-RS	Afo - Deluxe (R-Small)	Each

Procare Beach Boy Cast Boot

Features and Benefits

- ✓ Sandal style boot with E.V.A rocker sole
- ✓ With sturdy elastic strap
- ✓ Velcro closure



Medent-Code	Mfr Code	Description	Unit of Sale
DPRBCBBS	PR81073	Procure Beach Boy Cast Boot Small	Each
DPRBCBBM	PR81075	Procure Beach Boy Cast Boot Medium	Each
DPRBCBBL	PR81077	Procure Beach Boy Cast Boot Large	Each

Cast Boot Shoe

Features and Benefits

Cast shoe is a premium shoe, constructed of only quality material to aid in mobilisation and to protect the cast from breakdown during periods of weight bearing.



Medent-Code	Mfr Code	Description	Unit of Sale
DPRBDCM	PR81015	Cast Boot Delux Medium	Each
DPRBDCL	PR81017	Cast Boot Delux Large	Each
DPRBOOTRS	PR81113	Cast Boot Rocker Small	Each
DPRBOOTRM	PR81115	Cast Boot Rocker Medium	Each
DPRBOOTCRL	PR81117	Cast Boot Rocker Large	Each

BSN Cast Shoe Adult

Features and Benefits

- ✓ BSN Cast shoe is a premium shoe, constructed of only quality material to aid in mobilisation and to protect the cast from breakdown during periods of weight bearing.



Medent-Code	Mfr Code	Description	Unit of Sale
DPSECS	58-35010	BSN Cast Shoe Adult Small	Each
DPSECM	58-35011	BSN Cast Shoe Adult Medium	Each
DPSECL	58-35012	BSN Cast Shoe Adult Large	Each

* Small (5.5 - 7.5 size), Medium (8 - 10 size), Large (9 - 11.5 size)

Procure Wrap Patella Knee

Features and Benefits

- ✓ Neoprene wrap with open patella (1/8 inches)
- ✓ Sleeve can accommodate buttress pad for additional patellar support
- ✓ Choice of nylon one side or nylon two sides



Medent-Code	Mfr Code	Description	Unit of Sale
DPRWKPXXL	PR82461	Procure Wrap Patella Knee XXL	Each

Elastic Knee Support/Stabiliser

Features and Benefits

- ✓ Reduces swelling from minor injuries
- ✓ Flexible compression provides knee support
- ✓ Made of seamless, four-way stretch elastic breathable material for a comfortable fit
- ✓ Open patella
- ✓ For knee pain, articular effusion & swelling, soft tissue injuries, chronic inflammation after operation or trauma



Medent-Code	Mfr Code	Description	Unit of Sale
DPEKSS	AOK1430	Elastic Knee Support Small	Each
DPEKSM	AOK1440	Elastic Knee Support Medium	Each
DPEKSL	AOK1450	Elastic Knee Support Large	Each

Arm Sling - Immobiliser With Waist Strap

Features and Benefits

Has the same features and benefits as that of the sling arm with strap support but has an additional padded cotton waist strap for greater stability.



Medent-Code	Mfr Code	Description	Unit of Sale
DPSIM	AOS1640	Arm Sling - Immobiliser With Waist Strap Medium	Each
DPSIL	AOS1650	Arm Sling - Immobiliser With Waist Strap Large	Each
DPSIS	AOS1630	Arm Sling - Immobiliser With Waist Strap Small	Each

Sling Arm With Strap Support

Features and Benefits

- ✓ Breathable, lightweight and durable material
- ✓ Padded cotton shoulder strap for greater comfort
- ✓ Envelope style with thumb support
- ✓ Universal design fits right or left arm
- ✓ Supports and immobilises the arm following soft tissue injuries, fracture or surgery



Medent-Code	Mfr Code	Description	Unit of Sale
DPSS	AOS1730	Arm Sling With Strap (Small)t	Each
DPSM	AOS1740	Arm Sling With Strap (Medium)	Each
DPSSL	AOS1750	Arm Sling With Strap (Large)	Each

Mesh Arm Sling

Features and Benefits

- ✓ Supports and immobilises the arm following soft tissue injuries, fracture or surgery
- ✓ Vented nylon mesh with shoulder strap
- ✓ Envelope style with thumb support
- ✓ Universal design fits right or left arm
- ✓ For shoulder, dislocation/subluxation, upper arm fracture, post surgery



Medent-Code	Mfr Code	Description	Unit of Sale
DPMASS	AOS2330	Mesh Arm Sling (Small)t	Each
DPMASM	AOS2340	Mesh Arm Sling (Medium)	Each
DPMASL	AOS2350	Mesh Arm Sling (Large)	Each

Universal Wrist Splint

Features and Benefits

- ✓ Metal palmar splint and elastic wrist strap to stabilise the wrist
- ✓ Universal design - fit the left or right hand
- ✓ Provides compression and support to injured and weak wrist
- ✓ For carpal tunnel syndrome, soft tissue injuries, wrist pain, wrist strains, wrist sprains



Medent-Code	Mfr Code	Description	Unit of Sale
DPUWSL	AOW11BE50	Universal Wrist Splint L	Each
DPUWSM	AOW11BE40	Universal Wrist Splint M	Each
DPUWSS	AOW11BE30	Universal Wrist Splint S	Each
DPUWSEL	AOW11BE60	Universal Wrist Splint L	Each

Neoprene Wrist Stabiliser

Features and Benefits

- ✓ Palmer and dorsal metal splints stabilise the hand and wrist
- ✓ Provides compression and heat retention
- ✓ Supports injured and weak wrist
- ✓ Hook and loop closure for easy adjustment and secure fit
- ✓ Wrist strap with D-ring for additional support and stability
- ✓ Made of 3mm neoprene
- ✓ Available in left or right hand



Medent-Code	Mfr Code	Description	Unit of Sale
DPNWSS	DPNWSS	Neoprene Wrist Stabiliser SL	Each
DPNWSSR	DPNWSSR	Neoprene Wrist Stabiliser SR	Each
DPNWSM	DPNWSM	Neoprene Wrist Stabiliser ML	Each
DPNWSMR	DPNWSMR	Neoprene Wrist Stabiliser MR	Each
DPNWSL	DPNWSL	Neoprene Wrist Stabiliser LL	Each
DPNWSLR	DPNWSLR	Neoprene Wrist Stabiliser LR	Each

Small 14 - 16cm, Medium Left 16 - 19cm, Large Left 19 - 21.5cm

OPPO Wrist Splint With Strap

Features and Benefits

- ✓ Cotton / elastic construction with a circumferential strap for improved wrist stabilisation
- ✓ Preformed aluminium stay, velcro closure



Medent-Code	Mfr Code	Description	Unit of Sale
DPACES	OPP2182S	Wrist Splint With Strap, S	Each
DPACEM	OPP2182M	Wrist Splint With Strap M	Each
DPACEL	OPP2182L	Wrist Splint With Strap L	Each

Small 14 - 16cm, Medium Left 16 - 19cm, Large Left 19 - 21.5cm

Thumb brace

Features and Benefits

- ✓ Stabilises the joints of thumb, relieves pain
- ✓ Provides thermal support and compression
- ✓ Anatomically contoured metal thumb splint to stabilise the thumb
- ✓ Hook and loop closure for easy adjustment
- ✓ 3mm thick neoprene
- ✓ For joint injury, rheumatoid arthritis, thumb pain, sprain/strain



Medent-Code	Mfr Code	Description	Unit of Sale
DPTBS	AOW1230	Thumb Brace Small	Each
DPTBM	AOW1240	Thumb Brace Medium	Each
DPTBL	AOW1250	Thumb Brace Large	Each

Finger Splint Zimmer (Aluminium)

Features and Benefits

- ✓ For immobilisation of fingers after injury or as a splint during treatment
- ✓ Also used as an arm immobiliser in neonatal treatment



Medent-Code	Mfr Code	Description	Unit of Sale
DPRSZ1	11001089	Finger Splint Zimmer 1/2 Inch x 18 Inch (Aluminium)	Each
DPRSZ2	11003089	Finger Splint Zimmer 1 Inch x 18 Inch (Aluminium)	Each

Finger Splints Stax

Features and Benefits

- ✓ Designed to protect and support the distal joint of the finger
- ✓ Sold individually and come either clear or skin tone
- ✓ Also known as Mallet finger splints



Medent-Code	Mfr Code	Description	Unit of Sale
DPRSF1	11001010006	Finger Splints Stax 5cm	Each
DPRSF2	11001020006	Finger Splints Stax 5.5cm	Each
DPRSF3	11001030006	Finger Splints Stax 5.8cm	Each
DPRSF8	11001040006	Finger Splints Stax 6.1cm	Each
DPRSF4	11001050006	Finger Splints Stax 6.3cm	Each
DPRSF6	11001055006	Finger Splints Stax 6.8cm	Each
DPRSF5	11001060006	Finger Splints Stax 7.3cm	Each
DPRSF7	11001070006	Finger Splints Stax 8.4cm	Each

Finger Stall Leather

Features and Benefits

- ✓ Soft leather fingerstalls help protect dressings
- ✓ Helps reduce the sweating that is associated with PVC
- ✓ Made from soft, flexible, durable leather



Medent-Code	Mfr Code	Description	Unit of Sale
DPRSFLM	782661	Finger Stall Leather Medium	Each
DPRSFL	782688	Finger Stall Leather Large	Each
DPRSFLXL	782696	Finger Stall Leather Extra Large	Each

Thumb Stall Leather

Features and Benefits

- ✓ Soft leather thumbstalls help protect dressings
- ✓ Helps reduce the sweating that is associated with PVC
- ✓ Made from soft, flexible, durable leather



Medent-Code	Mfr Code	Description	Unit of Sale
DPRSTLL	783250	Thumb Stall Leather (Large)	Each
DPRSTLXL	407283	Thumb Stall Leather (Extra Large)	Each

Collar Neck Hard Adjustable Adult

Features and Benefits

- ✓ Ambu Adjustable Neck Collar - Adult (Multifit) - Single



Medent-Code	Mfr Code	Description	Unit of Sale
DPCOLMF	COLNECADJB	Collar Neck Hard Adjustable Adult	Each

Cervical Collars

Features and Benefits

- ✓ Soft foam collars designed to limit cervical flexion
- ✓ Made from high density, non latex foam
- ✓ Covered in natural 100% cotton stockinette to provide a comfortable and contoured fit



Medent-Code	Mfr Code	Description	Unit of Sale
DPCOL1	SM402	Cervical Collar 100 x 480mm. (Large)	Each
DPCOL2	SM401	Cervical Collar 90 x 480mm. (Medium)	Each
DPCOL3	SM400	Cervical Collar 75 x 480mm. (Small)t	Each

Cardboard Splint

Features and Benefits

The Cardboard Splint is a low cost, compact and convenient splint with a capacity to immobilise a broad range of injuries



Medent-Code	Mfr Code	Description	Unit of Sale
DPCARL	FRHE-535	Cardboard Splint Large (L77cm, W33cm)	Each
DPCARS	FRHE-536	Cardboard Splint Small (L41cm, W25cm)	Each

CATEGORY:MOBILITY

Crutches Adult Double Adjustable Forearm (Pair)

Features and Benefits

Ergonomic Crutches. Ergonomic shaped handle to help distribute the pressure more evenly across the hand. The cuff provides a secure fitting to the forearm and they are adjustable between cuff and handle & handle and ground.

Medent-Code	Mfr Code	Description	Unit of Sale
DPRCRUT2	CRU691800	Crutches Adult Double Adjustable Forearm Large (890mm -1370mm)	Pair
DPRCRUT1	CRU691400	Crutches Adult Double Adjustable Forearm Medium (890 X 1200mm)	Pair
DPRCRUT5	CRU691900	Crutches Adult Double Adj Forearm Ergonomic Tall	Pair
DPRCRUT4	CRU691500	Crutches Adult Double Adj Forearm Ergonomic Med	Pair



Ergonomic Crutches

Stoppers For Crutches (Pair)

- ✓ Standard Crutch Tip
- ✓ Grey non-marking rubber tips
- ✓ High quality rubber tips reinforced with steel to prolong their life
- ✓ Provide good non-slip traction
- ✓ Suit most canes, walking sticks, crutches and anything requiring a non-slipper rubber stopper



Medent-Code	Mfr Code	Description	Unit of Sale
DPRCRUTS	16MM	Stoppers For Crutches Pair (Internal 16mm)	Pair
DPRCRUTS1	CRU686320	Stoppers For Crutches Pair (Internal 20mm)	Pair

To Order Contact Us at:

Medent Medical

Unit 4, 110 Whakatiki St
PO Box 40356
Upper Hutt
New Zealand

Free Phone: 0800 263 3368

Phone: 64 4 527 7702

Free Fax: 0800 527 7499

Fax: 64 4 5277499

Email Enquires: enquiry@medent.co.nz

Finding a product and ordering it – you choose an option that suits you.
Either talk with one of our team or order right here online!

- Order from our sales representative when they visit you.
To arrange a visit please phone us on 0800 263 3368
- You can phone us on 0800 263 3368
or Fax us on 0800 527 7499

All products and more available at www.medent.co.nz

Medent Barr Guarantee

If You Are Not 100% Satisfied
With This BARR'S Branded
Product Return It for A 100%
Refund Used or Not

MEDENT
BARR
100%

