

ViroVac[®]

Surgical Smoke Evacuator



Compact Design

Allows for easy placement on shelf or counter and accommodates stackability.

Patented 3-Port Filter Design

Versatile 3-port filter design features covered, hidden male ports and accommodates a variety of tubing sizes providing versatility in procedure and accessory setup.

Variable Filter Life (up to 35 hours)

Delivers unsurpassed value and reduces per procedure costs.

Filter Tracking Technology

Automatically tracks and records filter life through a patented technology.

Simple LED Display

Two button design (motor speed and on/standby) makes this system extremely easy to operate.

Quiet Operation

Addresses noise concerns.

GPCEM Conference Offer

Save \$500 on list price Plus SAVE another \$500 with our free starter pack including: Filter + Autoclavable Blade + Ball + Needle

Was ~~\$4000~~ Now \$2999

Offer ex GST ends July 31st 2019

Interested? Contact Cooper at Medent Medical 0275302671



Auckland Ph 09 574 6832 Fax 09 574 6476
Wellington Ph 04 527 7702 Fax 04 527 7499
National Ph 0800 363 3368 Fax 0800 527 7499



ViroVac®

Surgical Smoke Evacuator

ViroVac surgical smoke plume evacuator provides ideal flow for physician offices and outpatient surgery centres. It combines customizable (variable) filter life, a compact shelf design, 3-port filter with filter life tracking technology, while exceeding expectations for value and efficacy.



Technical Information

Noise Levels
(Nominal Values) @
maximum power < 55 dBA

Dimensions
(H x W x D)
6 in x 11 in x 15.5 in
15.2 cm x 27.9 cm x 39.4 cm

Flow Rates
25 cfm (708 lpm) using 7/8 in
(22 mm) tubing

Weight
With filter installed: 12.0 lbs
(5.0 kg)
Without filter installed: 9.5 lbs
(4.0 kg)

Warranty
One (1) Year Full Warranty

Product Ordering Information

Description	Catalogue Number	Quantity
ViroVac® (100/120 Volt Systems) Includes one filter and one footswitch	VV120	1/Each
ViroVac® (200/240 Volt Systems) Includes one filter and one footswitch	VV220	1/Each
ViroVac® (100/120 Volt Systems) Includes one footswitch	VV120 NF	1/Each
V (200/240 Volt Systems) Includes one footswitch	VV220 NF	1/Each
EZLink® Automatic activation device can be used to remotely turn the evacuator on and off, saving filter life	EZLINK01	1/Case
ViroSafe® Fluid Trap Designed to capture fluids to help maintain the integrity of the filtration system and assist with maintaining filter life	VSFT10	10/Case
ViroSafe® ViroVac® Filter Four layers of proven filtration gross particulate, ULPA media, activated carbon, and post filter	VS353 VS35302 VS35304	1/Case 2/Case 4/Case

Accessories
(each sold separately)



EZLink®

**ViroSafe®
Fluid Trap**



Auckland Ph 09 574 6832 Fax 09 574 6476
Wellington Ph 04 527 7702 Fax 04 527 7499
National Ph 0800 363 3368 Fax 0800 527 7499



Filtration Group®

BUFFALO FILTER