



3M™ Red Dot™ Paediatric Solid Gel Electrode 2248

NPC Code	FDK223	
Indication for Use	Long Term Monitoring	
Dimensions	Electrode Size:	4.44cm diameter
	Skin Contact Size:	4.44cm diameter
	Adhesive Area:	10.45 sq cm
	Height excluding connector:	0.13cm
Electrode Materials	Backing Material:	3M™ Micropore™ Tape
	Backing Material Adhesive:	Micropore tape adhesive
	Connector:	Stainless Steel Snap
	Release Liner:	Si coated paper
	Sensor Material:	Silver/silver-chloride coated plastic
Sensor	Gel system:	Crossed linked guar gum (solid gel)
	Gel area:	5.03 sq cm
	Sensor area:	2.45 sq cm
Lifetime	Recommended max application time:	3 days
	Sealed pouch:	2 years
X-Ray and MRI	X-Ray:	NO
	MRI:	NO
Environmental Issues	PVC-free electrode:	YES
	Latex-free electrode:	YES
	PVC-free packaging	YES
Packaging Quantities	Pouch:	50

3M Health Care Limited

3M House
Morley Street
Loughborough
Leicestershire LE11 1EP
Telephone (01509) 611611
Facsimile (01509) 237288
www.3mhealthcare.co.uk

3M Ireland

3M House, Adelphi Centre
Dun Laoghaire
Co Dublin
Ireland
Telephone (01) 280 3555
Facsimile (01) 280 3509

3M and Red Dot are trademarks
of the 3M Company.
© 3M Health Care 2008.

