

*Medent Barr Guarantee*

If You Are Not 100% Satisfied  
With This BARR'S Branded  
Product Return It for A 100%  
Refund Used or Not

MEDE  
NT  
BARR  
100%

**medent** - **BARR'S**  
**MEDICAL**

# PRODUCT CATALOGUE

2015 | Edition 21



 **BARR'S**  
The Practitioners First Choice

[Online Sales & Ordering](#)

**medent**  
**MEDICAL**

HEAD OFFICE:  
PH 04 527 7702 FAX 04 527 7499  
PH OUTSIDE WELLINGTON: 0800 26 333 68  
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)



## About Us

Medent-Barr Agencies has been trading in the New Zealand market for more than **30 years** and are well recognised as an efficient supplier to many medical organisations including public hospitals and retirement homes.

**Our expertise is appreciated** by busy healthcare professionals and we work to provide you with a personalised service to meet your exact needs.

We are **100% New Zealand** owned and well positioned with our head office and distribution centre in Wellington.

## Just like you, we believe that the human touch is important

We are proud of the efforts we go to, add a **human touch** to our services and are one of the few suppliers still maintaining a team of sales representatives.

Our Wellington based service team includes a pharmacist and nursing/paramedic qualified team member to **better assist you** with complex queries.

Specialist staff consistently work to maintain and develop our systems and ordering processes to ensure the most efficient and professional service for you.

With Medent-Barr you will work alongside our knowledgeable and efficient team to access the best medical supplies for your clients and patients. You can be sure you will get the right information and the right product for the right application.

We have sole agency agreements and partnerships with internationally recognised medical manufacturers. This means you have an assurance of both top quality products from world leaders in medical manufacturing and a wealth of professional knowledge regarding clinical information and product applications.

We also conduct regular training with manufacturers to ensure that we deliver the best value to customers by providing the latest advice on product selection and use.

You have my personal guarantee of quality product and quality service.

**Martin Barr**  
Managing Director



### Medent Barr Guarantee

If You Are Not 100% Satisfied  
With This BARR'S Branded  
Product Return It For A 100%  
Refund Used or Not

MEDENT  
BARR  
100%

## About Medent Medical

After 30 years supplying the medical profession, Medent Medical remains 100% kiwi owned and 100% committed to delivering to you the best of medical, surgical, pharmaceutical and equipment supplies whenever you need them.

## BISMOS - The innovative inventory management tool

We wanted to make it easier for you to order, store and manage your medical supplies. We came up with BISMOS. A complete stock management and ordering solution that takes the pain out of supply management.

With the help of BISMOS you can:

- ✓ Automatically generate your next order by scanning products in your storeroom
- ✓ Easily check stock levels
- ✓ Improve practice productivity
- ✓ Cut medical supply ordering in half
- ✓ Reduce out of stock items
- ✓ And much more

[Discover how BISMOS can work for you](#)

## Providing you the best of medical, surgical & pharmaceutical supplies

With a constantly expanding inventory with thousands of products to choose from, including all of the major brands, when it comes to [medical supplies](#) think of Medent Medical as your 'go to' company. We have been supplying the medical profession with the best of medical, surgical, pharmaceutical and equipment supplies for more than 30 years. That is 30 years of working for you to simplify your ordering needs. You simply won't find a faster, more efficient way to manage your inventory.

Ask how we can help you:


- ✓ Optimise your ordering process
- ✓ Streamline inventory replacement
- ✓ Facilitate easy online ordering
- ✓ Enhance your productivity
- ✓ Improve ordering efficiency
- ✓ And save yourself both time and money

For more choice, better service and expert advice on your medical supplies, [contact us today](#)

## BISMOS AUTOMATED ORDERING AND STORAGE SOLUTIONS



- ✓ Save your time spent ordering supplies
- ✓ Reduce order error rate
- ✓ Easy checking of what stock you have
- ✓ Flexible modular system



## TABLE OF CONTENTS


<b>CATEGORY: BANDAGES – COHESIVE</b> .....	1
Barr's Elastic Cohesive Bandage.....	1
Barr's Cohesive Kids Bandage Colour Combo.....	1
<b>CATEGORY: BANDAGES – COMPRESSION/CREPE</b> .....	2
Barr's Conforming Elastic Gauze Bandage.....	2
Barr's Softicrepe Bandage.....	2
<b>CATEGORY: BANDAGES – SLING/TUBULAR</b> .....	3
Barr's Sling Triangular Bandage Disposable.....	3
Barr's Tubicare Tubular Bandage.....	3
<b>CATEGORY: CATHETERS / URINARY</b> .....	4
Barr's Catheter Drainage Bag 2 Litre Push/Pull Outlet, 90cm.....	4
Barr's Catheter Drainage Bag 2 Litre Push/Pull Outlet, 120cm.....	4
Barr's Catheter Drainage Leg Bag 750ml 'T-Tap' Outlet.....	5
Barr's Catheter Drainage Leg Bag 2L 'T-Tap' Outlet.....	5
Barr's Catheter Drainage Leg Bag 750ml 'T-Tap' Outlet Adjustable Tube.....	6
Barr's 2 Litre Catheter Bag No Bottom Outlet.....	6
<b>CATEGORY: DRESSINGS – DRY – GAUZE/NON-ADHERENT</b> .....	7
Barr's Sterile Non-Woven Gauze Swabs.....	7
Barr's 4ply Gauze Swabs Non Sterile, Non Woven.....	7
Barr's Soft Cotton Gauze Swabs.....	8
Barr's Help-It Non-Adherent Dressing Pad.....	8
<b>CATEGORY: DRESSINGS - DRY - COTTON WOOL</b> .....	9
Barr's Cotton Wool Roll Basic Grade 500g.....	9
Barr's Cotton Wool Zig-Zag Fold.....	9
<b>CATEGORY: DRESSINGS - WOUND MANAGEMENT</b> .....	10
Surgin Pore Surgical Dressing.....	10
Bosui Pad.....	10
Barr's Foam Bandage 10cm x 3mm x 3m Non Sterile.....	11
<b>CATEGORY: EQUIPMENT- MEDICAL</b> .....	12
Digital Blood Pressure Monitor (Manual) LCD Screen (Alpk2).....	12
Sphygmomanometer Pocket Aneroid Double Tube Alpk2 Set (500-V-EU).....	12
Sphygmomanometer Pocket Aneroid Head Only for Double Tube Set (Alpk2).....	13
Sphygmomanometer One Hand Pocket Aneroid Single Tube Alpk2 (Set).....	13
Sphygmomanometer Mercurial Desk 300mm Alpk2 (300-V-EU).....	14
Sphygmomanometer Bulb With Control Valves (Alpk2) .....	14
Stethoscope Coloured (Alpk2).....	15
Sphygmomanometer Cuff with Bag (Alpk2).....	15

<b>CATEGORY: GLOVES</b> .....	16
Supermax Powder – Free Gloves .....	16
Supermax Lightly Powdered Gloves .....	16
Supermax Latex Free Powder Free Bluple Nitrile Non Sterile Gloves .....	17
Glove Dispenser Single Barr’s Flexi Perspex .....	18
Glove Dispenser Treble Barr’s Flexi Perspex .....	18
<b>CATEGORY: INFECTION CONTROL</b> .....	19
Barr’s Smoothlube Lubricating Gel 500 ml Pump Pot .....	19
Barr’s Alcohol Prep Pads .....	19
<b>CATEGORY: SUTURES</b> .....	20
DemeTech Nylon Mono Sutures .....	20
DemeTECH’s PGA 910 Absorbable, Colourless (Polyglycolic Acid) Sutures .....	21
Demetech Plain Catgut Beige Sutures .....	22
<b>CATEGORY: CONSUMABLES</b> .....	23
Barr’s Paper Tape (White) .....	23
Barr’s Autoclave Water 3L .....	23
Barr’s Cotton Applicator 15cm Single Ended (100’s) .....	24
Barr’s Non-Sterile Cotton Applicator 7.5cm Double Ended (100’s) .....	24
Barr’s Jumbo Cotton Applicator 20’S .....	24



# Bandages

- ✓ [Barr's Elastic Cohesive Bandages](#)
- ✓ [Barr's Cohesive Kids Bandage Colour Combo](#)
- ✓ [Barr's Conforming Elastic Gauze Bandage](#)
- ✓ [Barr's Softicrepe Bandages](#)
- ✓ [Barr's Sling Triangular Bandage Disposable](#)
- ✓ [Barr's Tubicare Tubular Bandage](#)



## CATEGORY: BANDAGES – COHESIVE

### Barr's Elastic Cohesive Bandage

#### Features & Benefits

- It sticks to itself and not to skin, hair or clothing
- Individually wrapped
- Air permeable
- Cool and comfortable to wear



✓ White cohesive bandage

Similar Products : Lastoban, Rinkilastic, Handygauze

Product Code	Size	Unit of Sale	Price
BAELR1B	2.5cm x 4.5m	Each	
BAELR2B	5cm x 4.5m	Each	
BAELR3B	7.5cm x 4.5m	Each	
BAELR10B	10cm x 4.5m	Each	
BAACOB1	10cm x 4.5m	Each	



✓ Blue cohesive bandage

### Barr's Cohesive Kids Bandage Colour Combo

- ✓ Barr's Help-It Cohesive Kids Colour Combo Bandage
- ✓ Contains Latex

Similar Products : Handy Gauze

Product Code	Size	Unit of Sale	Price
BAELRCC1	2.5cm x 4.5cm	Pack of 10	
BAELRCC2	5.0cm x 4.5cm	Pack of 10	





## CATEGORY: BANDAGES – COMPRESSION/CREPE

### Barr's Conforming Elastic Gauze Bandage

- ✓ A light conforming elastic gauze bandage
- ✓ Made from high quality porous cotton/nylon blend

#### Features & Benefits

- Each bandage has soft non-fray edges is fully washable, re-usable and autoclavable
- Ideally suited for light support and bandage conforms to all body contours, so ideal for covering joints and bone structures
- Slit-strip easy opening pack
- 1.8m unstretched
- Approximately 4m stretched



Product Code	Size	Unit of Sale	Price
BACFB5	7.5cm x 1.8m	Packet of 12	
BACFB1	10cm x 1.8m	Packet of 12	

**Similar Products :** Hospilux, Crinx, Handyband

### Barr's Softicrepe Bandage

- ✓ High quality 100% cotton bandage with feathered edges

#### Features & Benefits

- Hospital use light weight crepe bandage
- Use for dressing retention and light support
- Individually wrapped bandages
- Each bandage comes with two dressing clips
- Recommended for single use
- Available in four different widths; 5cm, 7.5cm, 10cm, and 15cm
- 1.5m in length un-stretched, 4m stretched



Product Code	Size	Unit of Sale	Price
BACRELW1	5cm x 1.5m	Each	
BACRELW2	7.5cm x 1.5m	Each	
BACRELW3	10cm x 1.5m	Each	
BACRB4	15cm x 1.5m	Each	

**Similar Products :** Hospicrepe, Idealcrepe, Elastolite

## CATEGORY: BANDAGES – SLING/TUBULAR

### Barr's Sling Triangular Bandage Disposable

- ✓ Help-It Disposable Triangular Bandage

#### Features & Benefits

- Provides support and protection for injured arms, wrists and hands
- Comes with two safety pins



Similar Products : Smith & Nephew Triangular Bandage


Product Code	Size	Unit of Sale	Price
BASLT2B	110 x 110 x 15 cm	Each	
BASLT3B	96 x 96 x 127 cm	Each	

### Barr's Tubicare Tubular Bandage

- ✓ Tubicare is a tubular bandage

#### Features & Benefits

- Designed to provide lasting effective support with complete freedom of movement
- Made of cotton with covered elastic threads laid into the fabric
- Design creates spirals that enable free range of movement
- Adjust to the contours of the body and distribute pressure evenly over the surface
- Easy to apply and reapply
- No pins or tapes needed to secure
- Box designed with dispensing opening



Similar Products : Tubigrip

Product Code	Size	Unit of Sale	Price
BATBTCC	'C', 6.75cm x 10m, medium arm, small ankles	Each	
BATBTCD	'D', 7.5cm x 10m, large arms, medium ankles, small knees	Each	
BATBTCE	'E', 8.75cm x 10m, large ankles, medium knees, small thighs	Each	
BATBTCF	'F', 10cm x 10m, large knees, medium thighs	Each	



## Catheter/Urinary Drainage Bag

- ✓ [Barr's Catheter Drainage Bag 2 Litre Push/Pull Outlet, 90cm](#)
- ✓ [Barr's Catheter Drainage Bag 2 Litre Push/Pull Outlet, 120cm](#)
- ✓ [Barr's Catheter Drainage Leg Bag 750ml 'T-Tap' Outlet](#)
- ✓ [Barr's Catheter Drainage Leg Bag 2L 'T-Tap' Outlet](#)
- ✓ [Barr's Catheter Drainage Leg Bag 750ml 'T-Tap,' Adjustable, 65cm](#)
- ✓ [Barr's 2 Litre Catheter Bag No Bottom Outlet, 90cm](#)



## CATEGORY: CATHETERS / URINARY

### Barr's Catheter Drainage Bag 2 Litre Push/Pull Outlet, 90cm

- ✓ Catheter drainage bag with a 90cm or 120 long tubing system
- ✓ 2 litre capacity with a push/pull outlet
- ✓ Available in sterile and non-sterile packaging

#### Features & Benefits

- Graduated approximate measurements along side of bag
- Bag has opaque backing
- Holes on top of bag to accommodate hanger
- Anti-reflux valve on inlet tube, stopping backflow from bag



Similar Products : Coloplast, Conveen, Simpla, Platinum

Product Code	Size	Unit of Sale	Price
CABNPPB	2 Litre (90cm tube)	Packet of 10	
CABNPPBE	2 Litre (90cm tube)	Each	

### Barr's Catheter Drainage Bag 2 Litre Push/Pull Outlet, 120cm

- ✓ Catheter drainage bag with a 120cm long tubing system
- ✓ 2 litre capacity with a push/pull drainage outlet
- ✓ Available in sterile and non-sterile packaging

#### Features & Benefits

- Double sealed for extra protection from leakage
- Graduated approximate measurements along side of bag
- Easy see-through bag
- Fast draining
- Anti-reflux valve on inlet tube, stopping backflow from bag



Similar Products : Coloplast, Conveen, Simpla, Platinum

Product Code	Size	Unit of Sale	Price
CABPP120	2L (120cm tube)	Packet of 10	
CABPP120E	2L (120cm tube)	Each	


### Barr's Catheter Drainage Leg Bag 750ml 'T-Tap' Outlet

- ✓ Catheter drainage bag, non-sterile with a 10cm - 60cm long tubing system
- ✓ 750ml capacity with a 'T-Tap' outlet
- ✓ Available in sterile and non-sterile packaging

#### Features & Benefits

- Elastic leg straps included
- Connector attached to bottom of 'T-Tap' enabling 2nd bag to be connected without disconnecting drainage system
- Anti-reflux valve on inlet tube, stopping backflow from bag
- Non-woven backing for improved comfort

**Similar Products :** Coloplast, Conveen, Simpla, Platinum



Product Code	Size	Unit of Sale	Price
CABTT10	750ml, 10cm tube	Packet of 10	
CABTT10E	750ml, 10cm tube	Each	
CABTT45	750ml, 45cm tube	Packet of 10	
CABTT45E	750ml, 45cm tube	Each	
CABTT60	750ml, 60cm tube	Packet of 10	
CABTT60E	750ml, 60cm tube	Each	

### Barr's Catheter Drainage Leg Bag 2L 'T-Tap' Outlet

- ✓ Catheter drainage bag, non-sterile with a 90cm - 120cm long tubing system
- ✓ 2L capacity with a 'T-Tap' outlet
- ✓ Available in sterile and non-sterile packaging

#### Features & Benefits

- Double sealed for extra protection from leakage
- Graduated approximate measurements along side of bag
- Easy see-through bag
- Fast draining
- Anti-reflux valve on inlet tube, stopping backflow from bag

**Similar Products :** Coloplast, Conveen, Simpla, Platinum




Product Code	Size	Unit of Sale	Price
CABTT90	2L, 90cm tube	Packet of 10	
CABTT90E	2L, 90cm tube	Each	
CABTT120	2L, 120cm tube	Packet of 10	
CABTT120E	2L, 120cm tube	Each	

## Barr's Catheter Drainage Leg Bag 750ml 'T-Tap' Outlet Adjustable Tube

- ✓ Catheter drainage leg bag, 750ml 'T-Tap' outlet, 65cm adjustable
- ✓ Available in sterile and non-sterile packaging

### Features & Benefits

- Elastic leg straps included
- Adjustable tube length, can cut tubing down to wanted length and attach supplied connector
- Connector attached to bottom of 'T-Tap' enabling 2nd bag to be connected without disconnecting drainage system
- Non-woven backing for improved comfort
- Anti-reflux valve on inlet tube, stopping backflow from bag



Product Code	Size	Unit of Sale	Price
CABD4	750ml, adjustable tube (65cm)	Packet of 5	
CABD4E	750ml, adjustable tube (65cm)	Each	

**Similar Products :** Coloplast, Conveen, Simpla, Platinum

## Barr's 2 Litre Catheter Bag No Bottom Outlet

- ✓ 2 Litre (no bottom outlet) catheter bag
- ✓ Available in sterile and non-sterile packaging

### Features & Benefits

- Catheter drainage bag with a 90cm long tubing system but no bottom outlet
- Double sealed for extra protection from leakage
- Graduated approximate measurements along side of bag
- Easy see-through bag, fast draining
- Anti-reflux valve on inlet tube, stopping backflow from bag



**Similar Products :** Coloplast, Conveen, Simpla, Platinum

Product Code	Size	Unit of Sale	Price
CABNNB	2L, 90cm tube, no bottom outlet	Packet of 10	
CABNNBE	2L, 90cm tube, no bottom outlet	Each	



## Dressing-Dry/Wound Management

- ✓ Barr's Sterile Non-Woven Gauze Swabs
- ✓ Barr's Gauze Swabs Non Sterile, Non Woven 4ply
- ✓ Barr's Soft Cotton Gauze Swabs
- ✓ Barr's Help-It Non-Adherent Dressing Pad
- ✓ Barr's Cotton Wool Roll Basic Grade 500g
- ✓ Barr's Cotton Wool Zig-Zag Fold
- ✓ Surgin Pore Surgical Dressing
- ✓ Bosui Pad
- ✓ Barr's Foam Bandage, Non Sterile



## CATEGORY: DRESSINGS – DRY – GAUZE/NON-ADHERENT

### Barr's Sterile Non-Woven Gauze Swabs

- ✓ Individually packed sterile non-woven gauze swabs

#### Features & Benefits

- High absorbency 4ply squares
- Either in two or five per individual peelable pouches



Product Code	Size	Unit of Sale	Price
DDGSB3	7.5 x 7.5cm (5's)	Box of 20	
DDGSB3E	7.5 x 7.5cm (5's)	Each	
DDGSB4	7.5 x 7.5cm (2's)	Box of 50	
DDGSB4E	7.5 x 7.5cm (2's)	Each	

Similar Products : Propax Gauze

### Barr's 4ply Gauze Swabs Non Sterile, Non Woven

- ✓ Non-Sterile Non-Woven Gauze Swabs

#### Features & Benefits

- Non-sterile non-woven 4ply gauze swabs
- Comes in two sizes, either 5x5cm or 7.5x7.5cm
- In sealed packets of 100 swabs



Product Code	Size	Unit of Sale	Price
DDGSB1	5x5cm (5's)	Pack of 100	
DDGSB2	7.5x7.5cm (2's)	Pack of 100	

Similar Products : Mesoft, Sofisorb, Multisorb



## Barr's Soft Cotton Gauze Swabs

- ✓ Non-sterile woven soft cotton gauze swab


### Features & Benefits

- 100% natural cotton fibres
- High absorbency and strength
- Widely useable general wound dressing
- 100 swabs per packets



Similar Products : Mesoft, Sofisorb, Multisorb

Product Code	Size	Unit of Sale	Price
DDGS6B	5x5cm 4ply	Packet of 100	
DDGS3B	7.5x7.5cm 4ply	Packet of 100	



## Barr's Help-It Non-Adherent Dressing Pad

- ✓ Non-adherent dressing pad that consists of a thin layer of absorbent cotton fibres
- ✓ Covered on both sides with a poly-net non-adherent film

### Features & Benefits

- Cost effective non-adherent pad
- Can be used either way up
- Effectual in reducing the adherence to a wide variety of wounds
- Available in two different sizes
- Box designed with dispensing opening

Similar Products : Telfa non-adherent pad, Melolin



Product Code	Size	Unit of Sale	Price
DDNAT1B	7.5x10cm	Box of 100	
DDNAT2B	5x7.5cm	Box of 100	

## CATEGORY: DRESSINGS - DRY - COTTON WOOL

### Barr's Cotton Wool Roll Basic Grade 500g

- ✓ Cotton Wool Roll Basic Grade

#### Features & Benefits

- Cotton wool roll basic grade
- Large 500g pack
- High grade cotton wool

Similar Products : Cutisoft



Product Code	Size	Unit of Sale	Price
DDCWRB	500g	Each	

### Barr's Cotton Wool Zig-Zag Fold

- ✓ Cotton wool with a zigzag fold

#### Features & Benefits

- Large 500gm pack high-grade cotton wool, approximately 15cm width
- Comes in convenient resalable pack

Similar Products : Cutisoft



Product Code	Size	Unit of Sale	Price
DDCWZIGB	500gm	Each	

## CATEGORY: DRESSINGS - WOUND MANAGEMENT

### Surgin Pore Surgical Dressing

- ✓ Soft cloth dressing
- ✓ Can be used to cover all wounds
- ✓ Particularly useful for covering surgical wound sites



#### Features & Benefits

- Surgin Pore consists of a non-woven fabric backing, which ensures superior air permeability
- With an island Telfa type dressing pad, which has a non-adherent net covering to reduce the dressing sticking to the wound while still being absorbent
- Surgin Pore is water resistant
- Uses a hypoallergenic light adhesive so it is easy on the skin
- Comes in six different convenient sizes

Product Code	Size	Unit of Sale	Price
DWICSP3	5x8cm, (#50)	Box of 50	
DWICSP4	6x10cm, (#60)	Box of 50	
DWICSP5	8x13cm, (#80)	Box of 30	


**Similar Products :** Primapor, Cutiplast, Cosmopor

### Bosui Pad

- ✓ Waterproof transparent hypoallergenic adhesive film dressing
- ✓ With a non-adherent skin coloured absorbent pad wound sites

#### Features & Benefits

- Bosui Pad protects the wound from bacteria, so reducing the risk of secondary infections occurring
- As a waterproof dressing, it allows for bathing, showering and swimming while keeping the wound dry
- Has moisture vapour permeability, allowing skin to breathe and prevents maceration of wound site
- With an absorbent non-adherent polyethylene mesh covering the wound pad to reduce sticking
- Comes in six different convenient sizes



**Similar Products :** Opsite Post-Op, Tegaderm with pad, Molnlycke Mepore

Product Code	Size	Unit of Sale	Price
DWIDTPB5	5x7cm	Box of 50	
DWIDTPB6	6x9cm	Box of 50	
DWIDTPB8	8x12cm	Box of 30	

## Barr's Foam Bandage 10cm x 3mm x 3m Non Sterile

- ✓ Non-sterile Foam Bandage

### Features & Benefits

- White foam bandage manufactured in New Zealand
- Use for non-sterile light compression applications
- Also used to protect slow healing wounds from knocks, scrapes, and rubbing
- Apply over normal dressing
- Not suitable for direct wound contact
- Standard length is 3m but can be cut to order



**Similar Products :** Molnlycke Mepilex

Product Code	Size	Unit of Sale	Price
BAFOAM2	10cm x 3mm x 3m	Each	



## Equipment-Medical

- ✓ Digital Blood Pressure Monitor (Man) LCD Screen (Alpk2)
- ✓ Sphygmomanometer Pocket Aneroid Double Tube (Alpk2)
- ✓ Sphygmomanometer Pocket Aneroid Head Only (Alpk2)
- ✓ Sphygmomanometer One Hand Pocket Aneroid ST (Alpk2)
- ✓ Sphygmomanometer Mercurial Desk 300mm (Alpk2)
- ✓ Sphygmomanometer Bulb With Control Valves (Alpk2)
- ✓ Stethoscope Coloured (Alpk2)
- ✓ Sphygmomanometer Cuff with Bag (Alpk2)



**CATEGORY: EQUIPMENT- MEDICAL****Digital Blood Pressure Monitor (Manual) LCD Screen (Alpk2)**

- ✓ Manual digital blood pressure monitor with LCD screen

**Features & Benefits**

- Large display screen makes easy to read
- Manual inflate, auto deflate
- Can store up to 30 measurements and their averages, when the number exceeds 30 the oldest data will be replaced with the newest data
- Cuff size 22-32cm circumference when in position
- Pressure cuff indication range 0-300mmHg
- Carry bag and 4 x AA batteries included
- Clinically proven accuracy

**Similar Products :** Accosan Greenlight, Omron, Panasonic



Product Code	Size	Unit of Sale	Price
EMSGDIA1		Each	

**Sphygmomanometer Pocket Aneroid Double Tube Alpk2 Set (500-V-EU)**

- ✓ Pocket aneroid sphygmomanometer, 300mmHg no pin stop manometer
- ✓ Standard adult cuff; 23-36cm circumference
- ✓ Durable tubing, bulb and cuff
- ✓ Durable vinyl carry bag, very portable

**Features & Benefits**

- One-handed use, with thumb and fingertip control
- Lightweight diecast case
- Supplied in pouch case with brass zip
- Range 0-300mmHg
- Latex free inflation system
- Supplied with adult size velcro cuff with coiled tubing

**Similar Products :** Welch Allyn, Riester, Erka, Kobold



Product Code	Size	Unit of Sale	Price
EMSGPA	Pocket Aneroid Double Tube	Each	

## Sphygmomanometer Pocket Aneroid Head Only for Double Tube Set (Alpk2)

- ✓ Head Only
- ✓ For Double Tube Set

For more information, please click [HERE](#)

Similar Products : Welch Allyn, Riester

Product Code	Size	Unit of Sale	Price
RSSGPH		Each	



## Sphygmomanometer One Hand Pocket Aneroid Single Tube Alpk2 (Set)

- ✓ One handed trigger release, inflate and deflate, blood pressure monitor (K2-700)
- ✓ Supplied with standard adult cuff
- ✓ Durable plastic design
- ✓ Durable bulb and cuff, vinyl carry bag, very portable

Similar Products : Accoson, Welch Allyn, Riester, Erka, Kobold



Product Code	Size	Unit of Sale	Price
EMSGPAST	One handed pocket size	Each	

## Sphygmomanometer Mercurial Desk 300mm Alpk2 (300-V-EU)

- ✓ Mercury blood pressure monitor
- ✓ Reliable sphygmomanometers widely used in New Zealand Primary Care

### Features & Benefits

- Japanese quality manufacturing
- Heavy duty nylon zipper case
- Range of parts
- Four cuff sizes available for most models
- Single, double tube models
- Models range from accurate mercury through to fully electronic



Similar Products :Riester, ADC

Product Code	Size	Unit of Sale	Price
EMSGDM	300mmHg	Each	

## Sphygmomanometer Bulb With Control Valves (Alpk2)

- ✓ Alpk2 Sphygmomanometer Control Valve with Bulb (No. 90).



For more information, please click [HERE](#)

Similar Products : Baum, Welch Allyn

Product Code	Size	Unit of Sale	Price
RSSGPBC	#90	Each	





## Stethoscope Coloured (Alpk2)

- ✓ Alpk2 brand with box, professional nursing stethoscope

### Features & Benefits

- Nurse aid stethoscope
- Dual head (diaphragm and bell)
- Comes in a range of bright colours
- Chestpiece: Bell & diaphragm
- Two way type, aluminum
- Binaural: Brass binaural
- Y-Tube: PVC Y tube
- Weight: 130g



**Similar Products :** 3M, Riester, Liberty, Erka Kobold, Heine, Welch Allyn

Product Code	Size	Unit of Sale	Price
EMSTNA		Each	

## Sphygmomanometer Cuff with Bag (Alpk2)

For more information, please click [HERE](#)

**Similar Products :** Accoson, GIMA, Riester, Erka Kobold, Heine, Welch Allyn



- ✓ Double Tube



- ✓ Single Tube

Product Code	Size	Unit of Sale	Price
RSSGPCB1	Cuff Adult Nylon Double Tube	Each	
RSSGPCB2	Cuff Child Nylon Double Tube	Each	
RSSGPCB3	Cuff Thigh Nylon Double Tube (40.6 - 66cm)	Each	
RSSGPCB4	Cuff Adult Large Nylon, Double Tube	Each	
RSSGPCB9	Cuff Adult Large Nylon Single Tube	Each	



## Examination Gloves

- ✓ Supermax Powder – Free Gloves
- ✓ Supermax Lightly Powdered Gloves
- ✓ Supermax Latex Free Powder Free Blue Nitrile Non Sterile Gloves
- ✓ Glove Dispenser Single Barr's Flexi Perspex
- ✓ Glove Dispenser Treble Barr's Flexi Perspex



## CATEGORY: GLOVES

### Supermax Gloves


- Made of natural rubber latex, ambidextrous, latex, medical grade, low -protien, non-sterile, hypoallergenic examination gloves
- Provide an excellent biological barrier, but not intended for use as a chemical barrier
- Every Supermax Glove comes with multiple quality standards and international registrations
- Similar Products include Gammex, Triflex, Protexis, Triumph, InTouch, Elite



ISO 11193-1 : ARTG 76912 : AUST STD 4011 BRITISH STANDARD : ASTM / USA : TUV GERMANY

### Supermax Powder – Free Gloves

- ✓ Easy on hands, chlorinated powder-free
- ✓ Micro textured surface with beaded cuff



Product Code	Size	Unit of Sale	Price
GLLPF1	X-Small	Box 100	
GLLPF2	Small	Box 100	
GLLPF3	Medium	Box 100	
GLLPF4	Large	Box 100	
GLLPF5	X-Large	Box 100	

### Supermax Lightly Powdered Gloves

- ✓ Lightly powdered for ease of donning
- ✓ Smooth surface with beaded cuff



Product Code	Size	Unit of Sale	Price
GLL2	Small	Box 100	
GLL3	Medium	Box 100	
GLL4	Large	Box 100	
GLL5	X-Large	Box 100	



# SUPERMAX

## Supermax Latex Free Powder Free Blue Nitrile Non Sterile Gloves



### Features

- Ambidextrous
- Non Sterile Exam Gloves
- 6-Mil Nitrile Synthetic
- Online Powder Free
- No NRL Protein
- High tensile to meet toughest demands
- Powder Free and 100% Latex free
- Beaded Cuff
- Finger Textured
- Blue Color
- 100 gloves per Box



510k registration



Standard Malaysian Glove  
Certificate Approval 2021



ISO 11193-1 Listing  
71106

- Australian TGA L 63625
- Directive 89/686/EEC Article 11B-PPE from SGS, UK
- Conformity Identification Seal from

### Supermax Gloves also meet the following standards and testing requirements:

- ✓ Water tightness : ASTM D5151:99
- ✓ Dimension : ASTM D3767:01
- ✓ Tensile : ASTM D412:98 & ASTM D573:99
- ✓ Protein Testing : ASTM D5712:99
- ✓ Residual Powder Testing : ASTM D6124:01
- ✓ Biocompatibility Testing (Samples sent to PTC, USA for testing)
- ✓ Bioburden Testing (Samples sent to microbiological lab for testing)

Product Code	Description	Unit of Sale
GLNPFS2	Gloves Supermax Blue Nitrile Latex Free Powder Free Small	Box 200
GLNPFS3	Gloves Supermax Blue Nitrile Latex Free Powder Free Medium	Box 200
GLNPFS4	Gloves Supermax Blue Nitrile Latex Free Powder Free Large	Box 200

[Online Sales & Ordering](#)




HEAD OFFICE:  
PH 04 527 7702 FAX 04 527 7499  
PH OUTSIDE WELLINGTON: 0800 26 333 68  
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)

## Glove Dispenser Single Barr's Flexi Perspex

- ✓ Barr's plastic portable and single glove dispenser

Product Code	Size	Unit of Sale	Price
GLDISP	Single	Each	



## Glove Dispenser Treble Barr's Flexi Perspex

- ✓ Barr's white plastic portable and treble glove dispenser

Product Code	Size	Unit of Sale	Price
GLDISP3	Treble	Each	





## Infection Control

- ✓ [Barr's Smoothlube Lubricating Gel 500 ml Pump Pot](#)
- ✓ [Barr's Alcohol Prep Pads](#)



## CATEGORY: INFECTION CONTROL

### Barr's Smoothlube Lubricating Gel 500 ml Pump Pot

- ✓ Clinically Formulated Lubricating Gel

#### Features & Benefits

- Formulated for use with disposable devices Barr's or Welch Allyn Vaginal Speculums Anoscopes and Sigmoidoscopes
- Clinically Formulated
- Hypoallergenic
- Odourless
- Microbiologically tested
- Water Based
- Non Staining
- Easy Application
- Handy refill pack available



Product Code	Size	Unit of Sale	Price
EMLUB4B	500G PUMP POT	Each	

**Similar Products :** Welch Allyn, Optilube, Surgigel

### Barr's Alcohol Prep Pads

- ✓ Alcohol prep w ipes

#### Features & Benefits

- Contain one pad (2ply) saturated with 70% Isopropyl alcohol
- Rapid evaporation to reduce sting sensation

**Similar Products :** Webcol alcohol w ipes, Midi-Swab



Product Code	Size	Unit of Sale	Price
ICPALCB	35x30mm 2ply	Box of 200	



# Sutures

- ✓ DemeTech Nylon Mono Sutures
- ✓ DemeTECH's Polyglycolic Acid (PGA) 910 Sutures
- ✓ Demetech Plain Catgut Beige Sutures





## CATEGORY: SUTURES

### DemeTech Nylon Mono Sutures

- ✓ Nylon monofilament suture

#### Features & Benefits

- Non-absorbable, sterile surgical suture composed of the long chain aliphatic polymers Nylon 6 and Nylon 6.6
- Blue, dyed with FDA listed dye(s)
- Elicit minimal inflammatory reaction in tissue, usually followed by gradual encapsulation of the suture by fibrous connective tissue
- Indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures
- Easy passage through tissues due to its low frictional coefficient
- Minimal tissue reactivity
- Excellent elasticity
- Exceptional for skin closures
- Small USP sizes are designed specifically for microsurgical and ophthalmic procedures



**Similar Products :** Tyco B Braun and J & J Nylon Sutures


Product Code	Size	Unit of Sale	Price
SUDETN1	45cm 3/0 Blue Ds19 3/8 Reverse Cut	Box of 12	
SUDETN2	45cm 4/0 Blue Ds19 3/8 Reverse Cut	Box of 12	
SUDETN3	45cm 5/0 Blue Ds19 3/8 Reverse Cut	Box of 12	
SUDETN4	45cm 6/0 Blue Ds12 3/8	Box of 12	

## DemeTECH's PGA 910 Absorbable, Colourless (Polyglycolic Acid) Sutures

- ✓ Synthetic, absorbable, sterile, surgical suture composed of 100% glycolide
- ✓ Chemical formula of DemeTECH's PGA is  $(CaHaO_2)_n$

### Features & Benefits

- Composed of polymers of glycolide acids
- Coated to allow for easier passage through tissues with only minimal drag
- Facilitate easy tissue passage, precise knot placement
- Smooth tie down, a decreased tendency to irritate tissue
- Greater tensile strength, a longer absorption period
- Clinical trials have shown that, approximately after two weeks 84% of DemeTECH's PGA initial strength remained
- After 3 weeks 51% of their tensile strength remained
- At four weeks, approximately 23% of DemeTECH's PGA linear tensile strength remained
- Nonantigenic, nonpyrogenic and elicit only mild tissue reactivity during the absorption process



**Similar Products:** Tyco B Braun and J&J Nylon Sutures

Product Code	Size	Unit of Sale	Price
SUDETPGA1	3.0 DS19 45CM	Box of 12	
SUDETPGA2	4.0 DS19 45CM	Box of 12	
SUDETPGA3	5.0 DS16 45CM	Box of 12	

## Demetech Plain Catgut Beige Sutures

- ✓ Absorbable, sterile, surgical suture
- ✓ Composed of highly purified connective tissue (mostly collagen) derived from either beef or sheep intestines

### Features & Benefits

- Covered in chromic salt to prolong integrity, minimize tissue reactivity, and minimize the chance of adverse reactions
- Tensile strength is retained for up to 14 days, and sometimes as long as 21 days
- Allows for smooth passage through tissue
- Easy knot tying characteristics
- Chromicized strands minimize the possibility of fraying and maximize the potential and integrity of its Chromic Catgut



**Similar Products :** Tyco B Braun and J & J Nylon Sutures

Product Code	Size	Unit of Sale	Price
SUDETPGA5	3.0 75CM 19MM <b>3/8 Circle</b> Reverse Cutting	Box of 12	



## Consumables

- ✓ [Barr's Paper Tape \(White\)](#)
- ✓ [Barr's Autoclave Water 3L](#)
- ✓ [Barr's Cotton Applicator 15cm Single Ended](#)
- ✓ [Barr's Cotton Applicator 7.5cm Double Ended](#)
- ✓ [Barr's Jumbo Cotton Applicator](#)




## CATEGORY: CONSUMABLES

### Barr's Paper Tape (White)

- ✓ Lightweight paper tape

#### Features & Benefits

- Gentle removal from skin
- Porous for superior air permeability
- Latex free
- Easy to tear by hand
- Individually wrapped
- Comes in two different widths, 12mm and 25mm



Similar Products : Micropore

Product Code	Size	Unit of Sale	Price
TAPAB1	12mm	Each	
TAPAB1BX	12mm	BOX/24	
TAPAB2	25mm	Each	
TAPAB2BX	25mm	BOX/12	

### Barr's Autoclave Water 3L

- ✓ Distilled water comes in 3L bottles, for use with autoclaves

#### Features & Benefits

- Distilled water using evaporation technique guarantees the quality for this product
- Safe for drinking if your autoclave is already full



Product Code	Size	Unit of Sale	Price
MEWATD4	3L	Each	

### Barr's Cotton Applicator 15cm Single Ended (100's)

- ✓ Cotton tipped w ooden applicator
- ✓ Single ended, non-sterile 15cm in length

Product Code	Size	Unit of Sale	Price
INAPPB2	15cm	Packet of 100	



### Barr's Non-Sterile Cotton Applicator 7.5cm Double Ended (100's)

- ✓ Cotton tipped w ooden applicator
- ✓ Non-sterile 7.5cm double ended

Product Code	Size	Unit of Sale	Price
INAPPB1	7.5cm	Packet of 100	



### Barr's Jumbo Cotton Applicator 20'S

- ✓ Jumbo cotton tipped applicator
- ✓ Single ended, 18cm in length, w ith a plastic shaft

Product Code	Size	Unit of Sale	Price
INAPPBJ	18cm S/E	Packet of 20	



## To Order Contact Us at:

### Medent Medical

Unit 4, 110 Whakatiki St  
PO Box 40356  
Upper Hutt  
New Zealand

**Free Phone: 0800 263 3368**

Phone: 64 4 527 7702

**Free Fax: 0800 527 7499**

Fax: 64 4 5277499

**Email Enquires:** [enquiry@medent.co.nz](mailto:enquiry@medent.co.nz)

Finding a product and ordering it – you choose an option that suits you.  
Either talk with one of our team or order right here online!

- Order from our sales representative when they visit you.  
To arrange a visit please phone us on 0800 263 3368
- You can phone us on 0800 263 3368  
or Fax us on 0800 527 7499

All products and more available at [www.medent.co.nz](http://www.medent.co.nz)

