



**NEW MODEL**



Introducing the new improved  
 WH121  
**WALMAY  
 CRYOSPRAY**

The well proven Walmay Cryospray has given over two decades of proven reliable service to Medical Practitioners. It is the instrument of choice for Dermatologists and Practitioners.

- ▶ Australian designed & manufactured
- ▶ ARTG Number 170073
- ▶ European approval CE 0805



- 1 Thicker sheathing around the Probe, giving it greater tensile strength.
- 2 Improved metallurgical principles in construction of the flask collar.
- 3 The centre of gravity of the Walmay Cryospray WH121 is towards the top of the instrument, making it easy for the doctor to balance and steady the instrument in operation.
- 4 Fine gauze Strainer on the end of the Siphon Tube to reduce contamination accessing the Regulator mechanism.
- 5 The Carry Bag containing five spray tips, Instructions for Use and FIVE YEAR Warranty card.
- 6 The base Stand ring provides a greater surface area for stability of the instrument when not in use.



Distributed by  
**MESENT MEDICAL**  
 www.medent.co.nz  
 Ph: 0800 263 3368

**FIVE YEAR WARRANTY**