

KENDALL



The First Choice of Vet is the First Choice for Pets™

tyco
Healthcare

Animal Health

Medical

- Needles & Syringes
- Biopsy Products
- Sharps Disposal
- Wound Care Dressings & Tape
- Kennel Pads & Supplies
- Catheters & Tubes
- Blood Collection

Suture

- MAXON™
- BIOSYN™
- POLYSORB™, DEXON™ II, DEXON S
- Chromic Gut
- NOVAFIL™
- MONOSOF™, DERMALON™, SURGILON™
- SURGIPRO™/SURGIPRO MESH
- SOFSILK™
- TI·CRON™

Operating Room Supplies

- Laparotomy Sponges
- X-Ray Detectable Sponges
- Hand Scrub
- SURGIDRAPE™

Surgical

- Open Surgery
 - Surgical Staplers
 - Skin Staplers
 - Staple Removers
- Minimally Invasive Surgery
 - Trocars & Access Ports
 - Hand Instruments
 - Endoscopic Suturing
 - Endoscopic Staplers
 - Clip Applicators
 - Energy Products


Chemotherapy Protection & Safety

- Protective Gloves
- Protective Apparel
- Chemotherapy Spill Kits
- Absorbent Products
- Vial Venting System
- IVA™ Security Seals
- Chemotherapy Waste Bags

Programs

1-800-962-9888, #4
www.kendallhq.com

ARGYLE Red Rubber Robinson Catheters




ARGYLE Red Rubber Robinson Catheters

- Red radiopaque rubber.
- Tapered funnel connector.
- Smooth rounded tip.

Code	French Size	Qty/Case
<i>Sterile</i>		
8887660085	8 Fr. (2.7mm)	100
8887660101	10 Fr. (3.3mm)	100
8887660127	12 Fr. (4.0mm)	100
8887660143	14 Fr. (4.7mm)	100
8887660168	16 Fr. (5.3mm)	100
8887660184	18 Fr. (6.0mm)	100
8887660200	20 Fr. (6.7mm)	100

ARGYLE Penrose Tubing



ARGYLE Penrose Tubing

- Radiopaque Latex.
- Individually packaged sterile.

Code	Flat Size	Diameter	Length	Qty/Case
8888513002	7/16" (11mm)	1/4" (6mm)	12" (30cm)	50
8888513200	1/2" (13mm)	5/16" (8mm)	12" (30cm)	50
8888513408	5/8" (16mm)	3/8" (10mm)	12" (30cm)	50
8888513606	3/4" (19mm)	1/2" (13mm)	12" (30cm)	50
8888513804	1" (25mm)	5/8" (16mm)	12" (30cm)	50
8888514000	1 1/4" (32mm)	3/4" (19mm)	12" (30cm)	50
8888514208	1 3/8" (35mm)	7/8" (22mm)	12" (30cm)	50
8888514406	1 1/2" (38mm)	1" (25mm)	12" (30cm)	50
8888514604	7/16" (11mm)	1/4" (6mm)	18" (46cm)	50
8888514802	1/2" (13mm)	5/16" (8mm)	18" (46cm)	50
8888515007	5/8" (16mm)	3/8" (10mm)	18" (46cm)	50
8888515205	3/4" (19mm)	1/2" (13mm)	18" (46cm)	50
8888515403	1" (25mm)	5/8" (16mm)	18" (46cm)	50
8888515601	1 1/4" (32mm)	3/4" (19mm)	18" (46cm)	50
8888515809	1 1/2" (38mm)	1" (25mm)	18" (46cm)	50