



A Wide Selection of Skin Staplers
to Meet a Wide Variety of Procedural Needs



COVIDIEN

positive results for life™

Single Use Skin Staplers

Covidien is a leader in skin staplers and offers a wide variety of skin staplers to meet most procedural needs. Each Covidien stapler will place one staple each time the instrument handles are activated.



MultiFire Premium™ Single Use Skin Stapler

The MultiFire Premium™ single use skin stapler is available with different sizes of staples, each contained in its own single use loading unit.

- Reloadable stapler can reload a maximum of 2 times for a total of 105 (35 x 3) staples
- 360° rotating head maximizes access
- Each stapler contains stainless steel staples of 0.51 Mm diameter wire. The staple span is 10.2 Mm before closure and approximately 4.8 Mm wide and 3.4 Mm high when closed. The 3.4 Mm staple height enhances cosmesis.
- The Premium™ extractor has spring-loaded handle which resets the starting position automatically



Appose™ Single Use Skin Stapler

The Appose™ single use skin stapler can be used for the closure of skin in a wide variety of surgical procedures.

- Two staple width options provide wound closure flexibility
- Two staple height options provide enhanced cosmesis
- Clear plastic cartridge enhance staple count visibility
- Ergonomically designed staple remover handle provides comfort and ease of use

Single Use Skin Staplers



Signet™ Single Use Skin Stapler

- Device comes preloaded with 35 staples which can reduce intraoperative instrument set-up time
- Heavy gauge staple for optimal staple strength



DFS™ Single Use Fascia Stapler with Stainless Steel Staples

- 20 Each wide staples in each DFS™ stapler provides an adequate supply for complete approximation
- 13.8mm wide staple width post-closure enables staples to gain sufficient tissue purchase

Order Information

Order Code	Description	Units per box
8886803512	Appose™ single use skin stapler	12
8886803712	Appose™ single use skin stapler	12
070614	DFS™ - 20W single use fascia stapler	6
059035	MultiFire Premium™ - 35 single use skin stapler	6
059037	MultiFire Premium™ - 35W single use skin stapler	6
150462	Premium™ skin staple remover	12
059036	MultiFire Premium™ - 35 SULU	12
059038	MultiFire Premium™ - 35W SULU	12
054006	Signet™ - 35W single use skin stapler	6

Always refer to Instructions for Use provided with the product for complete indications, contraindications, warnings, precautions, and instructions.



COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. © 2011 Covidien. 4.11 M110257