



Dressing Wound



PRODUCT CATALOGUE

2015 | Edition 1



[Online Sales & Ordering](#)



HEAD OFFICE:
PH 04 527 7702 FAX 04 527 7499
PH OUTSIDE WELLINGTON: 0800 26 333 68
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)

About Medent Medical

After 30 years supplying the medical profession, Medent Medical remains 100% kiwi owned and 100% committed to delivering to you the best of medical, surgical, pharmaceutical and equipment supplies whenever you need them.

BISMOS - The innovative inventory management tool

We wanted to make it easier for you to order, store and manage your medical supplies. We came up with BISMOS. A complete stock management and ordering solution that takes the pain out of supply management.

With the help of BISMOS you can:

- ✓ Automatically generate your next order by scanning products in your storeroom
- ✓ Easily check stock levels
- ✓ Improve practice productivity
- ✓ Cut medical supply ordering in half
- ✓ Reduce out of stock items
- ✓ And much more

[Discover how BISMOS can work for you](#)

Providing you the best of medical, surgical & pharmaceutical supplies

With a constantly expanding inventory with thousands of products to choose from, including all of the major brands, when it comes to [medical supplies](#) think of Medent Medical as your 'go to' company. We have been supplying the medical profession with the best of medical, surgical, pharmaceutical and equipment supplies for more than 30 years. That is 30 years of working for you to simplify your ordering needs. You simply won't find a faster, more efficient way to manage your inventory.

Ask how we can help you:

- ✓ Optimise your ordering process
- ✓ Streamline inventory replacement
- ✓ Facilitate easy online ordering
- ✓ Enhance your productivity
- ✓ Improve ordering efficiency
- ✓ And save yourself both time and money

For more choice, better service and expert advice on your medical supplies, [contact us today](#)

BISMOS AUTOMATED ORDERING AND STORAGE SOLUTIONS



- ✓ Save your time spent ordering supplies
- ✓ Reduce order error rate
- ✓ Easy checking of what stock you have
- ✓ Flexible modular system




TABLE OF CONTENTS

CATEGORY: DRESSING WOUND – ACRYLIC	1
3M™ Tegaderm™ Absorbent Clear Acrylic Dressing.....	1
CATEGORY: DRESSING WOUND – ALGINATE	2
Algisite M Calcium Alginate Dressing Smith & Nephew.....	2
Algivon Alginate Dressing With Manuka Honey.....	2
Convatec Aquacel Dressing.....	3
Convatec Aquacel Foam Dressing.....	3
Convatec Kaltostat Calcium Alginate Dressing.....	4
Seasorb Calcium Alginate Dressing.....	4
3M Tegaderm Alginate Rope Dressing.....	4
CATEGORY: DRESSING WOUND – ANTIMICROBIAL	5
Actilite Manuka Honey Dressings With Activon.....	5
Alevyn Ag Silver Smith & Nephew.....	5
Aquacel Convatec Silver (Ag) Dressing.....	6
AMD Kendall Antimicrobial Foam Dressing.....	7
Inadine Adaptic Tulle Gras.....	7
AMD Kerlix Antimicrobial Gauze Roll	7
Iodosorb Dressing (Iodine).....	7
Silver Dressing Acticoat Smith & Nephew.....	8
Silver Dressing Atrauman Ag.....	8
3M Tegaderm Ag Mesh Silver Dressing	9
CATEGORY: DRESSING WOUND – CLOSURE	10
Leukosan Skinlink Smith & Nephew.....	10
Premium Skin Extractors (Staple Remover).....	10
Royal Skin Staplers 35w.....	10
Wound Closure System Steristrip 3M Nexcare.....	11
Wound Closure Leukostrip Smith & Nephew.....	12
CATEGORY: DRESSING WOUND – CLOTH WITH PAD	12
Covermed Cloth Dressing With Pad Roll Smith & Nephew.....	12
Cutiplast (With Pad) Island Cloth Smith & Nephew.....	13
Medipore Island Dressing 3M.....	13
MoInlycke Mepore Dressing.....	14
Primapore (With Pad) Island Cloth Smith & Nephew.....	14
Surgin Pore Island Cloth.....	15
CATEGORY: DRESSING WOUND – CLOTH WITHOUT PAD	15
Medipore Soft Cloth Surgical Tape 3M.....	15
Fixomull Stretch Soft Cloth Smith & Nephew	16
Hypafix Soft Cloth Retent Sheet Smith & Nephew.....	16
MoInlycke Mefix Soft Cloth Strip.....	17

CATEGORY: DRESSING WOUND – FILM WITH PAD	17
Bosui Pad Film Island Dressing.....	17
Cutifilm Plus Island Dressing Smith & Nephew.....	18
MoInlycke Mepore Film & Pad	18
Opsite Post-Op Island Dressing Smith & Nephew.....	19
Tegaderm (Pad) Island Dressing 3M.....	19
CATEGORY: DRESSING WOUND – FILM WITHOUT PAD	20
Opsite Flexifix Film Dressing Smith & Nephew.....	20
Opsite Flexigrid Film Dressing Smith & Nephew.....	20
Opsite IV 3000 Dressing Smith & Nephew.....	21
Tegaderm Film Dressing 3M.....	21
Tegaderm IV Film Dressing 3M.....	22
Tegaderm Film Roll 3M.....	22
CATEGORY: DRESSING WOUND – FOAM	23
Allevyn Hydro Cellular Dressing Smith & Nephew	23
Allevyn Adhesive Hydro Cellular Dressing Smith & Nephew	23
Allevyn Cavity Circular Foam Hydro Dressing Smith & Nephew	24
Allevyn Gentle Border Adhesive Smith & Nephew	24
Allevyn Thin Hydro Colloid Dressing Smith & Nephew.....	25
Convatec Versiva XC Adhesive Dressing.....	25
Copa Hydrophilic Foam Dressings.....	26
MoInlycke Lyofoam.....	26
MoInlycke Mepilex Foam	26
MoInlycke Mepilex Foam Border.....	27
Tegaderm Foam Adhesive Dressing Oval 3M.....	27
Tegaderm Foam Non Adhesive Dressing 3M.....	28
MoInlycke Mepilex Foam Border.....	28
CATEGORY: DRESSING WOUND – GEL S/HYDROGELS/SILICONE/POWDER.....	29
Activon Manuka Honey Tube.....	29
Tricotex Wound Layering Dressing Smith & Nephew	29
Aqueous Cream (BP).....	30
Cica Care Silicone Gel Sheet Smith & Nephew	30
Convatec Duoderm Gel Tube	30
Convatec Saf Gel Tube.....	30
Durafiber Absorbent Non Woven Gelling Fibre Smith & Nephew	31
Duratouch Dressings Silicon Smith & Nephew.....	31
Flamazine Cream (Silver Sulfadiazine 1%) Pom Smith & Nephew.....	31
Gel Burnaid Tube.....	32
Gel Intrasite Wound Smith & Nephew.....	32
Iodosorb Ointment Tubes Smith & Nephew.....	32
MoInlycke Mepiform.....	33
MoInlycke Mepitel Soft Silicone.....	33
MoInlycke Mesalt Debriders.....	33
Opsite Flexifix Gentle Roll Film Dressing with Silicone Gel Smith & Nephew	34

Purilon Hydrogel Tube	34
Silflex Dressings Silicon.....	35
Solosite Gel Tube Smith & Nephew.....	35
CATEGORY: DRESSING WOUND – HYDROCOLLOIDS.....	36
Collodion Flex BP Liquid To Dry.....	36
Comfeel Mini Hydro Colloid Dressing.....	36
Convatec Duoderm Hydrocolloid Dressing	37
Cutinova Hydro Dressing Smith & Nephew.....	38
Replicare Thin Hydro Colloid Dressing Smith & Nephew	38
Tegaderm Hydrocolloid Thin Dressing 3M.....	38
CATEGORY: DRESSING WOUND – MEDICATED	39
Alevyn Life Dressing Medium Smith & Nephew	39
Gelocast Zinc Impregnated Bandage (Stretched) Roll Smith & Nephew.....	39
Ichthopaste Bandage Smith & Nephew.....	40
Viscopaste Zinc Paste Bandage Smith & Nephew.....	40
Zipzoc Zinc Oxide Stockings Smith & Nephew	40
CATEGORY: DRESSING WOUND – ODOUR CONTROL	41
Carbonet Odour Control Dressing Smith & Nephew.....	41
Convatec Carboflex Odour Control.....	41
CATEGORY: DRESSING WOUND – PLASTERS (BANDAIDS)	42
Bandaid Plasters	42
Covermed Injection Plasters Smith & Nephew.....	42
Coverplast First Aid Dressing Smith & Nephew.....	43
Elastoplast Hypo Allergy Transparent Spots Smith & Nephew	44
Help-it Fabric Plaster.....	44
Healband Plastic Spots/Plasters.....	44
Injection Pad IV Pressure Pad Defries	45
Nexcare Plasters 3M	45
Nexcare Opticlude Regular Plasters (3M 1539).....	45
Plasters Blue Metal Detect Waterproof	46
Wound Dressing.....	46
CATEGORY: DRESSING WOUND – PROCEDURE PACKS.....	46
Dialysis Procedure Tray Smith & Nephew.....	46
Propax Lumps and Bumps Smith and Nephew.....	46
Theatre Pack Sterile Dressing Pack.....	47
Wound Dressing Pack 14 Sagepack Gauze.....	47
Wound Dressing Pack Propax 6 Swabs Smith & Nephew.....	47
Wound Dressing Pack Propax Gauze Balls Smith & Nephew.....	47
CATEGORY: DRESSING WOUND – RIBBON (PACKING) GAUZE.....	48
Intrasite Conformable Ribbon (Packing) Gauze Smith & Nephew.....	48
Ribbon Gauze Roll.....	48
CATEGORY: DRESSING WOUND – SPRAYS/CREAMS/OINTMENTS	48
Barrier Cream Dimethicone 5%.....	48

Cadesorb Pro Modulating Tube.....	49
Cavilon Barrier Spray No Sting 3M.....	49
Cavilon Durable Barrier Cream 3M.....	49
Cavilon Extra Dry Skin Cream 3M.....	50
Cetaphil Moisturizing Lotion (Pump).....	50
Convatec Duoderm Paste.....	50
Cream Anthisan.....	51
Cream Anti-flamme Tube.....	51
Cream Aquacare.....	51
Cream Arnica Tube Weleda.....	52
Cream Arnica Natures Kiss.....	52
Cream Calendula Healing Tube Weleda.....	52
Cream Nutraplus 10% Tube.....	53
Cream Secura (Extra Protection) Tube Smith & Nephew.....	53
Emulsifying Ointment.....	53
Flamazine Cream (Silver Sulfadiazine 1%) PM Smith & Nephew.....	54
Health E Fatty Cream.....	54
Pro Block Sunscreen, SPF.....	54
Hirudoid Cream.....	55
Liquid Paraffin Sterile.....	55
Lotion Alpha Keri.....	55
DP Lotion-HC 1%.....	56
Lyderm Cream.....	56
Mebo Anti Itch Ointment.....	56
Mebo Burn Repair Ointment.....	57
Mebo Wound Repair Ointment.....	57
NU-GEL® Hydrogel with Alginate.....	57
Oil Alpha Keri.....	58
Ointment Bepanthen.....	58
Ointment Ungvita.....	58
Ointment Vitamin A.....	59
Opsite Spray Dressing.....	59
Paste Magnesium Sulphate.....	59
Sunblock SPF 50 Pro Bloc.....	60
Vaseline Jelly Pot.....	60
Waxsol Ear Drop.....	60
Zinc & Castor Ointment.....	61
CATEGORY: DRESSING WOUND – TULLE GRAS	61
Actisorb Plus 25 Tulle Gras.....	61
Activon Manuka Honey Tulle Dressings.....	61
Adaptic Digit Finger Dressing.....	62
Adaptic Non Adherent Tulle Gras.....	62
Atrauman Dg Tulle Gras.....	62
Bactigras (Antiseptic) Tulle Gras Smith & Nephew.....	63

CuticerinTulle Gras Smith & Nephew	63
Jelonet Tulle Gras Sterile Smith & Nephew	63
Paranet Tulle Gras.....	64
CATEGORY: DRESSING WOUND – WIPES	64
Bzk Antiseptic Towelette Alcohol Free Wipes.....	64
Cavilon No Sting Barrier Film Wipes 3M.....	64
Convatec Remove Adhesive Wipes Convacare.....	65
Wipes Remove Smith & Nephew.....	65
CATEGORY: DRESSING WOUND – OTHERS.....	65
Kendall Argyle Penrose Tubing.....	65
Puritan Wound Depth Indicator.....	66
Seroma Wound Drainage System & Suc Res.....	66
Wound Mapping Grid 3M.....	66
Wound 3M Measuring Guide Pad.....	67

CATEGORY: DRESSING WOUND – ACRYLIC


3M™ Tegaderm™ Absorbent Clear Acrylic Dressing

Clear acrylic absorbent pad allows wound monitoring without disturbing the dressing. Enables clinicians to monitor pressure ulcers, skin tears and other moderately exuding wounds without changing the dressing. Saves clinical time.



Features & Benefits

- ✓ Effective barrier to contaminants such as fluid, bacteria and pathogens
- ✓ Hypoallergenic and latex-free
- ✓ Dressing has no upper wear time, fewer dressing changes aid healing
- ✓ Patients can shower with dressing in place
- ✓ Provides excellent healing rates while increasing patient comfort
- ✓ Reduces the risk of infection due to protection against bacterial and viral contamination




Medent Code	Mfr Code	Description	Unit of Sale
DWHTDA	90800	3M™ Tegaderm™ Absorbent Clear Acrylic Dressing, Oval 7.6 x 9.5 cm	Box 5
DWHTDA1	90801	3M™ Tegaderm™ Absorbent Clear Acrylic Dressing Oval 11.1 x 12.7 cm	Box 5
DWHTDA2	90802	3M™ Tegaderm™ Absorbent Clear Acrylic Dressing Square 14.9 x 15.2 cm	Box 5
DWHTDA3	90803	3M™ Tegaderm™ Absorbent Clear Acrylic Dressing Oval 14.2 x 15.8 cm	Box 5

CATEGORY: DRESSING WOUND – ALGINATE

Algisite M Calcium Alginate Dressing Smith & Nephew

Smith & Nephew Algisite M is a calcium-alginate dressing which forms a soft, integral gel when it comes into contact with wound exudate. It utilises the proven benefits of moist wound management, which results in faster healing of wounds.



Features & Benefits

- ✓ Easy to remove
- ✓ Conforms to wound contours
- ✓ Highly absorbent
- ✓ Fast gelling, high mannuronic acid fibres
- ✓ Minimal fibre shed construction
- ✓ High integrity when wet
- ✓ Soft and conformable



Medent Code	Mfr Code	Description	Unit of Sale
DWCAAL1	66000520	Smith & Nephew Algisite M Calcium-Alginate Dressing 10 x 10 cm	Box 10
DWCAAL2	66000519	Smith & Nephew Algisite M Calcium-Alginate Dressing 5 x 5 cm	Box 5
DWCARAL	66000522	Smith & Nephew Algisite M Calcium-Alginate Rope Dressing 2G (40CM)	Box 5

Algivon Alginate Dressing With Manuka Honey


Algivon is a soft alginate dressing impregnated with 100% medical grade Manuka honey. The alginate fibres enable a sustained, slower release of honey. Algivon is an ideal choice for wetter wounds as the alginate has a small capacity to absorb, meaning the honey isn't washed away with exudate therefore staying at the wound site for longer.

The dressing is very soft and conformable, ideal for cavities and debriding and de-sloughing large areas of necrotic and sloughy tissue.

Features & Benefits

Manuka honey brings the following benefits:

- ✓ Anti-bacterial, kills harmful bacteria
- ✓ Anti-inflammatory
- ✓ Eliminates odours without masking them
- ✓ Osmotic effect, drawing harmful tissue away from the wound bed
- ✓ Maintains the ideal moist wound healing environment




Medent Code	Mfr Code	Description	Unit of Sale
DWCAAL3	CR4231	Algivon Plus Ribbon Manuka Honey Alginate 2.5cm x 20cm	Box 5
DWCAAL4	CR3659	Algivon Manuka Honey Alginate Dressings 10 x 10 cm	Box 5

Convatec Aquacel Dressing

AQUACEL Dressing is a sterile, white, fibrous dressing derived from 100% sodium carboxymethylcellulose and incorporates the unique Hydrofiber Technology. It is a primary, versatile dressing indicated for use on moderately and highly exuding chronic and acute wounds.

Medent Code	Mfr Code	Description	Unit of Sale
DWCAA1	420671	Convatec Aquacel Dressing Extra 5 x 5 cm	Box 10
DWCAA2	420672	Convatec Aquacel Dressing Extra 10 x 10 cm	Box 10
DWCARA	403770	Convatec Aquacel Rope/Ribbon 2 x 45cm	Box 5



Convatec Aquacel Foam Dressing

AQUACEL® Foam Dressing has a soft absorbent Foam pad, an AQUACEL® interface, a gentle silicone adhesive, and a waterproof/bacteria barrier.

AQUACEL® Foam dressing is the only Foam that has comfort, simplicity and the healing benefits of an AQUACEL® Interface.

The comfort and simplicity of a foam plus the healing benefits of AQUACEL® dressing. Available in a range of sizes and shapes in adhesive – with a gentle silicone border designed to adhere to surrounding skin – and non-adhesive.



Medent Code	Mfr Code	Description	Unit of Sale
DWFAAQ1	420804	Convatec Aquacel Foam Adhesive Dressing 8x8 cm	Box 10
DWFAAQ2	420680	Convatec Aquacel Foam Adhesive Dressing 10x10 cm	Box 10
DWFAAQ3	420619	Convatec Aquacel Foam Adhesive Dressing 12.5x12.5 cm	Box 10
DWFAAQ4	420625	Convatec Aquacel Foam Adhesive Dressing (Heel) 8x8 cm	Box 5
DWFAAQ5	420626	Convatec Aquacel Foam Adh Dressing (Sacral) 8x8 cm	Box 5
DWFAQ1	420631	Convatec Aquacel Foam Non Adhesive Dressing 5x5 cm	Box 10
DWFAQ2	420633	Convatec Aquacel Foam Non Adhesive Dressing 10x10 cm	Box 10
DWFAQ3	420635	Convatec Aquacel Foam Non Adhesive Dressing 15x15 cm	Box 5

Convatec Kaltostat Calcium Alginate Dressing

KALTOSTAT Calcium Sodium Alginate Dressing is an alginate dressing indicated as a primary dressing for moderately to highly exuding chronic and acute wounds, and for wounds with minor bleeding. On contact with exudates the dressing gradually forms a firm, moist gel which is highly absorbent.



Medent Code	Mfr Code	Description	Unit of Sale
DWCAK1	168210	Convatec Kaltostat Calcium Alginate Dressing 5 x 5 cm	Box 10
DWCAK2	168212	Convatec Kaltostat Calcium Alginate Dressing 7.5 x12cm	Box 10
DWCARK	168117	Convatec Kaltostat Rope 2gm Calcium Alginate (Firm)	Box 5

Seasorb Calcium Alginate Dressing

SeaSorb Soft Dressing is an alginate dressing consisting of 85% calcium alginate and 15% Carboxy methylcellulose (CMC). The structure is unique because the calcium alginate and CMC molecules are intertwined. Calcium alginate and CMC molecules swell when they absorb exudate. The structure ensures that the wound exudate is absorbed vertically resulting in local absorption which helps to minimise the risk of maceration.



When the dressing comes into contact with exudate, the calcium ions from the dressing are replaced by sodium ions present in wound exudate. The ion exchange allows the dressing to form a very soft, cohesive and hydrophilic gel when managing exudate.

Medent Code	Mfr Code	Description	Unit of Sale
DWCASS	29103705	Seasorb Calcium Alginate Dressing 5 x 5 cm	Box 30
DWCRCS	291037402	Seasorb 2gm Rope Calcium Alginate 44cm	Box 6

3M Tegaderm Alginate Rope Dressing

Sterile, highly conformable, and absorbent primary wound care dressings. Available in a variety of sizes. Designed to maintain high integrity of the saturated dressing. Easily removed from most wounds in one piece with minimal gel residue. Becomes an amorphous gel-like mass when saturated. May be gently irrigated from the wound bed at dressing change.

Provides clinicians with a choice of products to meet their needs for exuding wounds. Both alginate dressings become a gel-like consistency as they absorb exudate which provides a moist wound healing environment with easy removal. Ideally suited for the management of moderately to heavily exuding partial and full thickness wounds. Both styles are compatible with the 3M family of secondary dressings.




Medent Code	Mfr Code	Description	Unit of Sale
DWCART	90120	3M Tegaderm Alginate Rope High Integrity 30.4cm	Box 5

CATEGORY: DRESSING WOUND – ANTIMICROBIAL

Actilite Manuka Honey Dressings With Activon

Actilite is a light non-adherent viscose net dressing coated with 99% Manuka honey and 1% Manuka oil. It is a non-sticky dressing ideal for infection prevention.

The dressing is designed to protect a wound, promote healing and allow the passage of exudate. The antibacterial effect of Actilite has been enhanced by combining high grade antibacterial Manuka oil with Manuka honey. The combination of Manuka honey and Manuka oil has been demonstrated in-vitro to be effective against a number of major wound infecting organisms including MRSA, VRE and Providentia stuartii.



Features & Benefits

- ✓ Antibacterial, kills harmful bacteria
- ✓ Osmotic effect, draws harmful tissue away from the wound bed
- ✓ Enhanced healing, anti-inflammatory
- ✓ Odour reduction, eliminates odours without masking them
- ✓ Conformable, non-adherent
- ✓ Maintains the ideal moist wound healing environment



Medent Code	Mfr Code	Description	Unit of Sale
DWHON1	CR3849	Actilite Manuka Honey Dressings With Activon 10 x 10cm	Box 10
DWHON2	CR4281	Actilite Manuka Honey Dressings With Activon 5 x 5cm	Box 10
DWHON3	CR3831	Algivon Alginate Dressing With Activon Manuka Honey 5 x 5cm	Box 5

Allevyn Ag Silver Smith & Nephew

Smith & Nephew Allevyn Ag - The first and only hydrocellular foam dressing containing silver sulfadiazine.

Features & Benefits

- ✓ Highly effective, proven to reduce bacterial load
- ✓ Excellent exudate management
- ✓ Waterproof, bacteria-proof top film is a barrier to pathogens and promotes healing
- ✓ Sustains antibacterial activity for up to 7 days

Allevyn Ag Non-Adhesive is an absorbent antimicrobial three layer dressing consisting of an SSD (Silver Sulfadiazine) containing absorbent hydrocellular pad sandwiched between a perforated non-adherent wound contact layer and a waterproof outer film.

The dressing provides an effective barrier to microbial contamination protecting the wound from invasive pathogenic micro-organisms, thus assisting faster healing.



Medent Code	Mfr Code	Description	Unit of Sale
DWHALAG	66800073	Smith & Nephew Allevyn Ag Silver Adhesive 7.5 x 7.5 cm	Box 10
DWHALAG1	66800078	Smith & Nephew Allevyn Ag Silver Adhesive 12.5 x 12.5 cm	Box 10
DWHALAG2	66800083	Smith & Nephew Allevyn Ag Silver Non Adhesive 5 x 5 cm	Box 10
DWHALAG3	66800086	Smith & Nephew Allevyn Ag Silver Non Adhesive 10 x 10cm	Box 10

Aquacel Convatec Silver (Ag) Dressing

AQUACEL® Ag Hydrofiber® Wound Dressing with Ionic Silver. AQUACEL® Ag dressing is available in a variety of sizes.

Does more than just provide antimicrobial activity. It is indicated for moderate to high exuding wounds which are infected or at risk of increased infection. The only antimicrobial dressing that incorporates Hydrofiber® Technology.

Features & Benefits

- ✓ Versatile primary dressing, soft and conformable
- ✓ Combines the favourable gelling characteristics of Hydrofiber Technology with the broad-spectrum antimicrobial properties of ionic silver (Ag+)
- ✓ Releases ionic silver in a controlled manner as wound exudate is absorbed into the dressing
- ✓ May help protect periwound skin by helping reduce the risk of maceration
- ✓ Supports wound healing by providing a moist wound healing environment
- ✓ Provides rapid and sustained antimicrobial activity in vitro
- ✓ Provides sustained antimicrobial activity for up to 7 days



Medent Code	Mfr Code	Description	Unit of Sale
DWSAQ1	420675	Convatec Aquacel Ag Extra Dressing 5cm x 5cm	Box 10
DWSAQ2	420676	Convatec Aquacel Ag Extra Dressing 10cm x 10cm	Box 10
DWSAQR	403771	Convatec Aquacel Ag Dressing Rope 2cm x 45cm	Box 5
DWSFAAQ1	420804	Convatec Aquacel Ag Foam Adhesive Dressing 8 x 8 cm	Box 10
DWSFAAQ2	420680	Convatec Aquacel Ag Foam Adhesive Dressing 10 x 10 cm	Box 10
DWSFAAQ3	420619	Convatec Aquacel Ag Foam Adhesive Dressing 12.5 x 12.5 cm	Box 10
DWSFAAQ4	420625	Convatec Aquacel Ag Foam Adhesive Dressing 17.5 x 17.5 cm	Box 10
DWSFAAQ5	420626	Convatec Aquacel Ag Foam Adhesive Dressing Sacral 17 x 20 cm	Box 10
DWSFAQ1	420639	Convatec Aquacel Ag Foam Non Adhesive Dressing 5 x 5 cm	Box 10
DWSFAQ2	420642	Convatec Aquacel Ag Foam Non Adhesive Dressing 10 x 10 cm	Box 10
DWSFAQ3	420645	Convatec Aquacel Ag Foam Non Adhesive Dressing 15 x 15 cm	Box 5



AMD Kendall Antimicrobial Foam Dressing

Ultra-soft foam dressing is an ultra-soft hydrophilic foam that protects and cushions wounds. Small pores decrease the chance of tissue migration into the foam dressing.



Medent Code	Mfr Code	Description	Unit of Sale
DWHKEN10	55544AMDX	AMD Kendall Antimicrobial Foam Dressing 10 x 10cm	Each

AMD Kerlix Antimicrobial Gauze Roll

Ultra-soft foam dressing is an ultra-soft hydrophilic foam that protects and cushions wounds. Small pores decrease the chance of tissue migration into the foam dressing.



Medent Code	Mfr Code	Description	Unit of Sale
DWHKER	3332	AMD Kerlix Antimicrobial Gauze Roll 11.4cm x 3.7m	Each

Inadine Adaptic Tulle Gras



INADINE PVP-I Non-Adherent Dressing consists of a low adherent knitted viscose fabric impregnated with a poly ethylene glycol (PEG) base containing 10% Povidone Iodine; equivalent to 1.0% available iodine.

Iodine is a highly effective topical antimicrobial that has been used clinically in the treatment of wounds. It has a broad spectrum of antimicrobial activity with efficacy against bacteria, microbacteria, fungi, protozoa and viruses.

Medent Code	Mfr Code	Description	Unit of Sale
DWTGIN1	PO1481	Adaptic Inadine Tulle Gras 5 x 5cm	Box 25
DWTGIN2	P01512	Adaptic Inadine Tulle Gras 9.5 x 9.5cm	Box 25

Iodosorb Dressing (Iodine)



Iodosorb removes barriers to healing by its dual action antimicrobial and desloughing properties. The broad spectrum antimicrobial action is provided by the sustained release of iodine and the desloughing action is provided by the unique cadexomer matrix. By effectively removing the barriers of bacteria, slough and debris and excess exudate, Iodosorb has been shown to accelerate wound healing in leg ulcers and is effective for the treatment of infection in chronic exuding wounds.

Medent Code	Mfr Code	Description	Unit of Sale
DWIODD	66051330	Iodosorb Dressing (Iodine) 4x6cm	Box 5

Silver Dressing Acticoat Smith & Nephew

ACTICOAT Dressings offer powerful antimicrobial barrier protection. Only ACTICOAT maintains safe bactericidal concentrations of silver with patented nanocrystalline technology.

The ACTICOAT family of dressings are flexible, absorbent and conformable and provide an efficient and effective barrier to bacterial penetration. ACTICOAT has been shown to offer significant clinical and economic results.

ACTICOAT3 has an absorbent inner core and a 3 day wear time. ACTICOAT7 has a 5 layer antimicrobial barrier dressing³ with an absorbent inner core and a longer, 7 day wear time. ACTICOAT FLEX is available in both a 3 day and 7 day variant and is suitable for use under NPWT for up to 3 days and ideal for packing, filling and conforming to difficult anatomical areas.

ACTICOAT ABSORBENT has all of the benefits of ACTICOAT in a conformable³, absorbent type 1 alginate dressing suitable for moderately exuding wounds. This dressing has a 7 day wear time and is suitable for use in cavity wounds.




Medent Code	Mfr Code	Description	Unit of Sale
DWSACT1	66000808	Smith & Nephew Silver Dressing Acticoat3 3 5cmx5cm	Box 5
DWACTAM	66000791	Smith & Nephew Silver Dressing Acticoat3 10cmx10cm	Box 12
DWSACT5	66000809	Smith & Nephew Silver Dressing Acticoat7 5cmx5cm	Box 5
DWSACT3	66800396	Smith & Nephew Silver Dressing Acticoat Flex3 5cmx5cm	Box 5
DWSACT2	66001760	Smith & Nephew Silver Dressing Acticoat Absorbent 5cmx5cm	Box 5
DWSACT	66000841	Smith & Nephew Silver Dressing Acticoat Absorbent 10cmx12.5cm	Box 5
DWSRA	66000840	Smith & Nephew Silver Dressing Acticoat Absorbent Rope 2cmx30cm	Box 5

Silver Dressing Atrauman Ag

Atrauman Ag Silver Dressing. Polyamide fabric coated with metallic silver and impregnated with a non medicated neutral ointment.

- ✓ Soft, thin support fabric drapes easily over the wound to ensure close contact
- ✓ Silver provides an antimicrobial barrier, which is useful for the management of contaminated wounds
- ✓ Help reduce the risk of bacterial colonisation
- ✓ Bacteria on the dressing surface are killed by silver ions which are activated by the moist wound environment
- ✓ Effective against Gram-negative and Gram-positive organisms
- ✓ Suitable for use in chronic wounds of various aetiologies
- ✓ Impregnated with a triglyceride based ointment
- ✓ Leaves minimal residue and reduces the risk of pain at the time of dressing removal



Medent Code	Mfr Code	Description	Unit of Sale
DWSATR	499571	Atrauman Silver Dressing Ag 5cmx5cm	Box 10
DWSATR10	499573	Atrauman Silver Dressing Ag 10cmx10cm	Box 10

3M Tegaderm Ag Mesh Silver Dressing

Versatile, easy to use 3M Tegaderm Ag Mesh Dressing with Silver

Features and Benefits

- ✓ Provides an effective antimicrobial barrier in an affordable silver dressing
- ✓ Nonwoven dressing containing silver sulfate which creates an effective antimicrobial barrier for up to 7 days. The soft, absorbent dressing is supplied sterile and may be custom cut
- ✓ The porous, non-occlusive dressing conforms to the wound base and wicks drainage into the dressing where the silver ions are available to reduce the number of bacteria and yeast
- ✓ Available in five sizes and ideal for use in wound clinics, hospitals, home care and long term care



Medent Code	Mfr Code	Description	Unit of Sale
DWSTAM	90500	Tegaderm Ag Mesh Silver Dressing 5cm x 5cm	Box 5
DWSTAM1	90501	Tegaderm Ag Mesh Silver Dressing 10cm x 12.7cm	Box 5

CATEGORY: DRESSING WOUND – CLOSURE

Leukosan Skinlink Smith & Nephew

Smith & Nephew Leukosan Adhesive is a topical skin closure adhesive formulated from a patented blend of 2-octyl and n-butyl cyanocrylate monomers, this combined with a flexible delivery system enables clinicians to confidently manage a range of topical wound closure situations.

Features and Benefits

- ✓ Choice of fine and sponge applicator tips to cater for varying wound presentations
- ✓ Aluminium tube ensures product integrity and minimises 'suck back' effect once finger pressure had been released on the tube during application
- ✓ Glue has a light violet dye to aid in visualisation during application



Medent Code	Mfr Code	Description	Unit of Sale
DWLEUB1	72542-00	Smith & Nephew Leukosan Skinlink 7mm x 60mm	Box 10
DWLEUB2	72543-00	Smith & Nephew Leukosan Skinlink 12mm x 76mm	Box 10

Premium Skin Extractors (Staple Remover)

Skin Staple Remover (Plastic with Metal Tips).

The Premium extractor features a springloaded handle that automatically returns the handle and jaw to the starting position.

Medent Code	Mfr Code	Description	Unit of Sale
EMREMS	150462	Premium Skin Extractors (Staple Remover)	Each



Royal Skin Staplers 35w

The Royal™ single-use skin stapler places one staple each time the instrument trigger is activated. The staples first penetrate the skin and are then formed, holding the tissue together.

For use in the following:

- ✓ Closure of incisions in vein stripping
- ✓ Thyroidectomy
- ✓ Mastectomy procedures
- ✓ Closure of scalp incisions
- ✓ For scalp flap hemostasis
- ✓ Skin grafts
- ✓ Plastic and reconstructive surgery



Medent Code	Mfr Code	Description	Unit of Sale
EMRSS	54887	Royal Skin Staplers 35w	Box 6

Wound Closure System Steri-Strip 3M Nexcare

3M Nexcare™ Steri-Strip™ Skin Closure. Sterile, flexible skin closure strips that are made of a porous, non-woven elastic backing coated with a pressure-sensitive, hypoallergenic adhesive and reinforced with polyester filaments for added strength. Steri-Strip elastic skin closures are ideal for use in contoured areas or areas of musculoskeletal movement.

Reinforced Steri-Strip skin closures provide general wound support for increased tensile strength and finer wound edge approximation.



Reinforced Steri-Strip

Features and Benefits:

- ✓ Popular skin closure in hospitals and is used for improved cosmetic results
- ✓ Ideal for securing, closing and supporting small cuts and wounds
- ✓ For wound support following suture or staple removal
- ✓ Breathable - for added comfort and holds firmly in place
- ✓ Hypoallergenic and made in U.S.A.



Elastic Steri-Strip



Medent Code	Mfr Code	Description	Unit of Sale
DWWCSS	W8512	Wound Closure System Steri-strip Single 3M Nexcare	Box 25
DWWCS2Y	W8516	Wound Closure System Steri-strip 5 Pack 3M Nexcare	Box 25
DWWCS5	R1547	Wound Closure System Steri-strip 3M Nexcare 12mm x 100mm	Box 50
DWWCS1	R1540	Wound Closure System Steri-strip 3M Nexcare 5 Pack 3mm x 75mm	Box 50
DWWCS4	R1546	Wound Closure System Steri-strip 3M Nexcare 10 Pack 6mm x 100mm	Box 50
DWWCS3	R1542	Wound Closure System Steri-strip 3M Nexcare 6 Pack 6mmM x 38mm	Box 50
DWWCSF1	1551	Wound Closure System Steri-strip 3M Nexcare 3 Pack 6mm x 75mm Flesh (Non-reinforced)	Box 50
DWWCS2	R1541	Wound Closure System Steri-strip 3M Nexcare 3 Pack 6mm x 75mm	Box 50
DWWCS5E	E4547	Wound Closure System Steri-strip 3M Nexcare Elastic 12mm x 100mm	Box 50
DWWCS1E	E4540	Wound Closure System Steri-strip 3M Nexcare Elastic Pack 5 3mm x 75mm	Box 50
DWWCS3E	E4542	Wound Closure System Steri-strip 3M Nexcare Elastic Pack 6 6mm x 38mm	Box 50
DWWCS2E	E4541	Wound Closure System Steri-strip 3M Nexcare Elastic Pack 3 6mm x 75mm	Box 50

Wound Closure LeukoStrip Smith & Nephew

Smith & Nephew LeukoStrip* is a skin closure strip combining a very strong hypoallergenic adhesive and a porous polyamide material which allows the free passage of exudate therefore lessening the likelihood of maceration under the strip.

Features and Benefits:

- ✓ Elastic polyamide textile
- ✓ Permeable to moisture and air
- ✓ Hypoallergenic adhesive
- ✓ Reliable adhesive strength
- ✓ Sterile and individually sealed
- ✓ Available in the traditional white or skin tone: ideal for faces and hands
- ✓ Reliably and securely adapts to the edges of the wound
- ✓ Does not traumatise the edges of the wound allowing almost pain free wound closure
- ✓ Minimal scarring
- ✓ Well tolerated by the skin



Medent Code	Mfr Code	Description	Unit of Sale
DWWCL3	66002880	Smith & Nephew Wound Closure LeukoStrip 13mm x 102mm Pack 6	Box 50
DWWCL5	66002882	Smith & Nephew Wound Closure LeukoStrip 4mm x 38mm Pack 4 Flesh	Box 50
DWWCL	66002876	Smith & Nephew Wound Closure LeukoStrip 4mm x 38mm Pack 6	Box 50
DWWCL4	66002903	Smith & Nephew Wound Closure LeukoStrip 4mm x 76mm Pack 4	Box 50
DWWCL1	66002879	Smith & Nephew Wound Closure LeukoStrip 6.4mm x 102mm Pack 5	Box 50
DWWCL6	66002878	Smith & Nephew Wound Closure LeukoStrip 6.4mm x 76mm Pack 3	Box 50
DWWCSF1	66002883	Smith & Nephew Wound Closure LeukoStrip 6.4mm x 76mm Flesh	Box 50

CATEGORY: DRESSING WOUND – CLOTH WITH PAD

Covermed Cloth Dressing With Pad Roll Smith & Nephew

Smith & Nephew Covered low allergy adhesive non-woven dressing.

Sensitive skin is liable to maceration, irritation, or allergic reactions when covered with an inappropriate adhesive dressing.

Covermed uses hypoallergenic adhesive with excellent skin tolerance, which can be gently removed from skin as well as from body hair.



Medent Code	Mfr Code	Description	Unit of Sale
DWICC1	72151-00	Smith & Nephew Covered Cloth Dressing With Pad Roll 4cm x 5m	Each
DWICC2	72153-00	Smith & Nephew Covered Cloth Dressing With Pad Roll 8cm x 5m	Each
DWICC3	72152-00	Smith & Nephew Covered Cloth Dressing With Pad Roll 6cm x 5m	Each

Cutiplast (With Pad) Island Cloth Smith & Nephew

Smith & Nephew Cutiplast* Steril is a soft conformable sterile island dressing especially useful in areas of curvature and flexion such as knees and elbows, hirsute areas and in situations of high humidity as the porous material breathes.

Consists of a non-lint low adherent wound pad covered with hypoallergenic Fixomull.



Features and Benefits:

- ✓ Low adherent
- ✓ Highly absorbent pad
- ✓ Air and water vapour permeable

Medent Code	Mfr Code	Description	Unit of Sale
DWICCUT1	66001473	Smith & Nephew Cutiplast Island Cloth With Pad 10cm x 8cm	Box 50
DWICCUT2	66001474	Smith & Nephew Cutiplast Island Cloth With Pad 15cm x 8cm	Box 50
DWICCP1	66001464	Smith & Nephew Cutiplast Island Cloth Strip With Pad 4cm x 5cm	Each
DWICCUT	66001478	Smith & Nephew Cutiplast Island Cloth With Pad 5cm x 7.2cm	Box 100
DWICCP2	66001466	Smith & Nephew Cutiplast Island Cloth Strip With Pad 6cm x 5cm	Each
DWICCUT3	66001468	Smith & Nephew Cutiplast Island Cloth Strip With Pad 8cm x 5cm	Each

Medipore Island Dressing 3M

3M Medipore +Pad soft cloth adhesive wound dressing is an all-in-one sterile dressing that features a breathable soft cloth tape backing with a unique non-adherent pad.

Features and Benefits:

- ✓ May be used as a primary or secondary dressing over light to moderately draining wounds, for securing devices, or for protecting I.V. catheter sites
- ✓ Highly conformable absorbent dressing that provides superior patient comfort as dressing flexes with body movement and contours
- ✓ Hypoallergenic, non-latex adhesive that provides excellent adhesion, yet is gentle to the skin
- ✓ Designed to promote normal wound healing processes to occur with less pain and trauma



Medent Code	Mfr Code	Description	Unit of Sale
DWIMED	3564	3M Medipore Island Dressing 6 x 10cm	Box 50

Molnlycke Mepore Dressing

Mepore[®] is the worlds first gentle, secure and long lasting self adhesive dressing. Its absorption qualities makes Mepore a suitable dressing for the treatment of low to moderately exuding wounds such as surgical incisions, minor burns, cuts and abrasions. The protective layer and rounded corners enhance its ease of application.



Mepore[®] Island Dressing: Easy-to-use self adhesive dressing with increased comfort and wear time

Mepore[®] Pro is a shower, viral and bacterial barrier (microbes larger than 25nm) for patients. It has optimal absorption capabilities and is designed with an easy to remove protection paper that enables aseptic as well as effortless application. It is also designed to minimise friction between the patient and their clothing/bed.



Mepore[®] Pro: Neutral toned comfortable dressing that protects against strike through

Medent Code	Mfr Code	Description	Unit of Sale
DWICMP1	331900	Molnlycke Mepore 5m x 4cm (With Pad) Island Cloth Roll	Each
DWICMP2	332000	Molnlycke Mepore 5m x 7cm (With Pad) Island Cloth Roll	Each
DWIDM5	671200	Molnlycke Mepore Island Dressing 9cm x 25cm	Box 30
DWIDM6	671300	Molnlycke Mepore Island Dressing 9cm x 30cm	Box 30
DWIDM1	670800	Molnlycke Mepore Island Dressing 6cm x 7cm	Box 60
DWIDM2	670900	Molnlycke Mepore Island Dressing 9cm x 10cm	Box 50
DWIDM3	671000	Molnlycke Mepore Island Dressing 9cm x 15cm	Box 50
DWIDM4	671100	Molnlycke Mepore Island Dressing 9cm x 20cm	Box 30
DWIDMP1	670820	Molnlycke Mepore Pro Sterile 6cm x 7cm	Box 60
DWIDMP2	670920	Molnlycke Mepore Pro Sterile 9cm x 10cm	Box 40
DWIDMP3	671020	Molnlycke Mepore Pro Sterile 9cm x 15cm	Box 40

Mepore Island Cloth Roll



Primapore (With Pad) Island Cloth Smith & Nephew

Primapore[®] Island Dressing is a soft conformable dressing which can be used to cover all wounds, particularly those in curved areas of the body, such as surgical wounds or general cuts and lacerations. Consists of a Melolin type pad and a water resistant low allergy Hypafix backing.



Medent Code	Mfr Code	Description	Unit or Sale
DWICP1	7135	Smith & Nephew Primapore (With Pad) Island Cloth 6cm x 8.3cm	Box 50
DWICP2	66000318	Smith & Nephew Primapore (With Pad) Island Cloth 15cm x 8cm	Box 20

- ✓ Soft breathable cover
- ✓ Highly absorbent pad
- ✓ Low Allergy Adhesive
- ✓ Conformable pad and cover

Surgin Pore Island Cloth

Exclusive to Medent Medical.

Surgin Pore/Surgin Pad is a surgical adhesive plaster with a non adherent pad. Made of a soft non woven fabric which allows air permeability, is water repellent and hypoallergenic. Suitable for surgical wounds, lacerations or general cuts. Surgin Pad has a thicker absorbent pad (10 ply cotton pad). CE Product. EOG Sterile. 5 years shelf life:

Features and Benefits:

- ✓ Surgical adhesive plaster with non adherent pad
- ✓ Made of soft non woven fabric
- ✓ Air permeable, water repellent, hypoallergenic
- ✓ Gentle on the skin
- ✓ Suitable for fragile skin, surgical wounds, lacerations, general cuts

Medent Code	Mfr Code	Description	Unit of Sale
DWICSP3	SP50	Surgin Pore Island Cloth 5 x 8 cm (Size No.50)	Box 50
DWICSP4	SP60	Surgin Pore Island Cloth 6 x 10 cm (Size No.60)	Box 50
DWICSP5	SP80	Surgin Pore Island Cloth 8 x 13 cm (Size No.80)	Box 30



CATEGORY: DRESSING WOUND – CLOTH WITHOUT PAD

Medipore Soft Cloth Surgical Tape 3M

3M Medipore- H Soft Cloth.

Features and Benefits:

- ✓ Soft Cloth, hypoallergenic tape
- ✓ More gentleness than Medipore tape
- ✓ Higher adhesion than Medipore tape
- ✓ Latex-free and hypoallergenic for sensitive patients
- ✓ Excellent cross and diagonal stretch to accommodate moderate swelling or movement
- ✓ Easy-tear perforated rolls for quick tailoring
- ✓ Water-resistant



Medent Code	Mfr Code	Description	Unit of Sale
DWSCMP1	2862	3M Medipore Soft Cloth Surgical Tape 5cmx9.1m	Each
DWSCMP2	2868	3M Medipore Soft Cloth Surgical Tape 20cmx9.1m	Each
DWSCMP3	2864	3M Medipore Soft Cloth Surgical Tape 10cmx9.1m	Each
DWSCMP4	2961	3M Medipore Soft Cloth Surgical Tape 2.5cmx9.1m	Pack 2



Fixomull Stretch Soft Cloth Smith & Nephew

Fixomull Stretch is a hypoallergenic, porous, adhesive, non-woven tape.

Features and Benefits:

- ✓ Skin-friendly and safe, quick-release backing paper
- ✓ For use on frequently mobile and highly contoured parts of the body
- ✓ Permeable to moisture and air, reducing the risk of maceration
- ✓ Radiotransparent so it does not need to be removed for examinations
- ✓ Hypo-allergenic poly acrylate adhesive adheres reliably
- ✓ Soft backing material of non-woven polyester conforms to body contours
- ✓ Can be cut to the required size for many uses
- ✓ Saves material, no unnecessary wastage



Medent Code	Mfr Code	Description	Unit of Sale
DWSCF1	02036-00	Smith & Nephew Fixomull Stretch Soft Cloth 5cm x 10m	Each
DWSCF2	02037-00	Smith & Nephew Fixomull Stretch Soft Cloth 10cm x 10m	Each
DWSCF3	02038-00	Smith & Nephew Fixomull Stretch Soft Cloth 15cm x 10m	Each
DWSCF4	02039-00	Smith & Nephew Fixomull Stretch Soft Cloth 20cm x 10m	Each

Hypafix Soft Cloth Retent Sheet Smith & Nephew

Hypafix Conformable Retention Tape is a low allergy adhesive, non-woven dressing retention sheet. Hypafix has excellent conformability that simplifies shaping around joints and awkward body contours, to allow the patient greater freedom of body movement. Complete coverage of the dressing reduces the risks of casual contamination. The broad range of sizes increases convenience and patient comfort.



Features and Benefits:

- ✓ Combines security with comfort
- ✓ Highly conformable
- ✓ Provides reliable fixation for catheters and drainage tubes
- ✓ Reduces risks of casual contamination
- ✓ Allows joints to be dressed without constraining movement



Medent Code	Mfr Code	Description	Unit of Sale
DWSCH1	71443-00	Smith & Nephew Hypafix Soft Cloth Retent Sheet 2.5cmx10m	Each
DWSCH2	71443-03	Smith & Nephew Hypafix Soft Cloth Retent Sheet 15cmx10m	Each
DWSCH3	71443-01	Smith & Nephew Hypafix Soft Cloth Retent Sheet 5cmx10m	Each
DWSCH4	71443-02	Smith & Nephew Hypafix Soft Cloth Retent Sheet 10cmx10m	Each
DWSCH5	71443-04	Smith & Nephew Hypafix Soft Cloth Retent Sheet 20cmx10m	Each

Molnlycke Mefix Soft Cloth Strip

MOLNLYCKE MEFIX SOFT CLOTH, self-adhesive fabric. Mefix® can be used in a wide variety of situations where fixation is required such as the fixation of dressings, swabs, catheters, and tubes.

Features and Benefits:

- ✓ Pre-measured protection paper for convenient measurement and cutting
- ✓ Perforated split in protection paper for easy handling
- ✓ Skin-friendly, waterbased, solvent-free polyacrylate adhesive for secure and gentle fixation
- ✓ Soft, elastic nonabsorbent nonwoven for patient comfort



Medent Code	Mfr Code	Description	Unit of Sale
DWSCMF1	310299	Molnlycke Mefix Soft Cloth 2.5cmx10m Strip	Each
DWSCMF2	310599	Molnlycke Mefix Soft Cloth 5cmx10m Strip	Each
DWSCMF3	311099	Molnlycke Mefix Soft Cloth 10cmx10m Strip	Each
DWSCMF4	311599	Molnlycke Mefix Soft Cloth 15cmx10m Strip	Each
DWSCMF5	312000	Molnlycke Mefix Soft Cloth 20cmx10m Strip	Each
DWSCMF6	313000	Molnlycke Mefix Soft Cloth 30cmx10m Strip	Each

CATEGORY: DRESSING WOUND – FILM WITH PAD

Bosui Pad Film Island Dressing

Exclusive to Medent Medical.


Bosui Pad: standard transparent dressing with wound pad. Special non-stick pad makes no pain when removing. Moisture permeability allows the skin to breathe. Waterproof and bacterial barrier allows bathing and showering. CE Product. EOG Sterile. 5 years shelf life:

Features and Benefits:

- ✓ Breathable film assures superior air permeability
- ✓ Waterproof & hypoallergenic
- ✓ Skin coloured pad blends gently with skin
- ✓ Dressing of choice for elderly skin
- ✓ Like Opsite Post-Op/Tegaderm With Pad


BOSUI PAD (30µm)

Standard transparent dressing with wound pad



Special non-stick pad

makes no pain when removing. Flesh colored pad blends with skin tone.



Moisture permeability

allows the skin to breathe.

Waterproof and bacterial barrier

allows bathing and showering.

Medent Code	Mfr Code	Description	Unit of Sale
DWIDTPB5	BP50	Bosui 5 x 7 cm (Pad) Film Island Dressing (Size No.50)	Box 50
DWIDTPB6	BP60	Bosui 6 x 9 cm (Pad) Film Island Dressing (Size No.60)	Box 50
DWIDTPB8	BP80	Bosui 8 x 12 cm (Pad) Film Island Dressing (Size No.80)	Box 30



Cutifilm Plus Island Dressing Smith & Nephew

Cutifilm Plus is a sterile dressing for the treatment of post operative wounds.

Features and Benefits:

- ✓ Sterile
- ✓ Transparent
- ✓ Hypoallergenic adhesive
- ✓ Bacteria proof and water proof
- ✓ Secure adhesion
- ✓ Low adherent pad



Medent Code	Mfr Code	Description	Unit of Sale
DWIDC1	76246	Smith & Nephew Cutifilm Plus Island Dressing 10cm x 12.5cm	Box 25
DWIDCP1	76244	Smith & Nephew Cutifilm Plus Island Dressing 8cm x 10cm	Box 25
DWIDCP2	76245	Smith & Nephew Cutifilm Plus Island Dressing 8cm x 15cm	Box 25
DWIDCP3	76247	Smith & Nephew Cutifilm Plus Island Dressing 10cm x 20cm	Box 25
DWIDCP4	76248	Smith & Nephew Cutifilm Plus Island Dressing 7.2cm x 5cm	Box 50



Molnlycke Mepore Film & Pad

Molnlycke Mepore Film & Pad, self-adhesive film dressing. Transparent fluid, viral and bacteria repellent that makes it easy to monitor the patient's wound.

The highly absorbent pad and high vapour permeable film offers excellent fluid handling. It is a skin friendly adhesive, gentle to the skin and at the same time provides secure fixation. Used to treat patients with low to moderate wound exudate levels, such as surgical wounds, cuts and abrasions.

Features and Benefits:

- ✓ Transparent dressing for easy wound monitoring that reduces risk for maceration
- ✓ Flexible film and high absorbency pad minimises risk of skin irritation
- ✓ Breathable, shower-, viral- and bacterial-proof backing film



Medent Code	Mfr Code	Description	Unit of Sale
DWMF1	275200	Molnlycke Mepore Film & Pad 5 x 7cm	Box 85
DWMF2	275400	Molnlycke Mepore Film & Pad 9 x 10cm	Box 30





Opsite Post-Op Island Dressing Smith & Nephew

Smith & Nephew OpSite Post-Op is a revolutionary post-operative dressing combining a high Moisture Vapour Transmission Rate (MVTR) film and the cushioning and absorbency of a low adherent pad together with an easy application carrier system.

The hydrophilic polyurethane film on OpSite Post-Op has Reactic properties that 'switch on' in the presence of moisture allowing any excess to transpire into the atmosphere.

OpSite Post-Op Visible is a transparent film with a unique see-through absorbent pad. The highly breathable film is waterproof and bacteria-resistant, and the unique lattice structure foam pad enables regular wound assessment without the need to lift or remove the dressing.

Ideal for moderately exuding wounds including surgical wounds, lacerations and abrasions. The low allergy adhesive reduces the risk of irritation and keeps the dressing in place, even over awkward contours.



Medent Code	Mfr Code	Description	Unit of Sale
DWIDOP1	66000710	Smith & Nephew Opsite Post-Op Island Dressing 10cm x 12cm	Box 10
DWIDOP2	66000713	Smith & Nephew Opsite Post-Op Island Dressing 10cm x 20cm	Box 20
DWIDOP3	66000712	Smith & Nephew Opsite Post-Op Island Dressing 15.5cm x 8.5cm	Box 20
DWIDOP4	66000708	Smith & Nephew Opsite Post-Op Island Dressing 6.5cm x 5cm	Box 100
DWIDOP5	66000714	Smith & Nephew Opsite Post-Op Island Dressing 25cm x 10cm	Box 20
DWIDOP6	66000715	Smith & Nephew Opsite Post-Op Island Dressing 30cm x 10cm	Box 20
DWIDOP7	66000709	Smith & Nephew Opsite Post-Op Island Dressing 9.5cm x 8.5cm	Box 20
DWIDOPV1	66800136	Smith & Nephew Opsite Post-Op Visible Dressing 8cm x 10cm	Box 20
DWIDOPV2	66800138	Smith & Nephew Opsite Post-Op Visible Dressing 20cm x 10cm	Box 20
DWIDOPV3	66800139	Smith & Nephew Opsite Post-Op Visible Dressing 25cm x 10cm	Box 20
DWIDOPV4	66800137	Smith & Nephew Opsite Post-Op Visible Dressing 15cm x 10cm	Box 20

Tegaderm (Pad) Island Dressing 3M

3M Tegaderm is a sterile composite dressing, waterproof bacterial barrier which consists of a non-adherent absorbent pad bonded to a larger thin film transparent dressing.

Tegaderm +Pad transparent dressings, with the unique 'picture frame' delivery system, make it quick and easy to achieve a precise, secure placement of the dressing every time.



Medent Code	Mfr Code	Description	Unit of Sale
DWIDTP1	3582	3M Tegaderm (Pad) Island Dressing 5 x 7 cm	Box 50
DWIDTP2	3584	3M Tegaderm (Pad) Island Dressing 6 x 10 cm	Box 50
DWIDTP3	3586	3M Tegaderm (Pad) Island Dressing 9 x 10 cm	Box 25
DWIDTP4	3587	3M Tegaderm (Pad) Oval Dressing 9 x 10 cm	Box 25
DWIDTP5	3589	3M Tegaderm (Pad) Island Dressing 9 x 15 cm	Box 25
DWIDTP6	3590	3M Tegaderm (Pad) Island Dressing 9 x 20 cm	Box 25
DWIDTP7	3591	3M Tegaderm (Pad) Island Dressing 9 x 25 cm	Box 25



CATEGORY: DRESSING WOUND – FILM WITHOUT PAD

Opsite Flexifix Film Dressing Smith & Nephew



Smith & Nephew OpSite Flexifix Transparent Film Roll.

OpSite Flexifix is a roll of transparent adhesive film, ideal for use in dressing fixation, tube fixation, dressing reinforcement and the treatment of painful peripheral neuropathy.

It combines all the advantages of **OpSite** film with a unique application system, to allow easy and secure fixation of all kinds of dressings. **OpSite**'s strong, flexible film also offers protection against skin breakdown due to friction and moisture.

Medent Code	Mfr Code	Description	Unit of Sale
DWTDOS2	66000040	Smith & Nephew OpSite Flexifix Film Dressing 5cm x 10m Roll	Each
DWTDOS3	66000041	Smith & Nephew OpSite Flexifix Film Dressing 10cm x 10m Roll	Each
DWTDOS5	66000375	Smith & Nephew OpSite Flexifix Film Dressing 15cm x 10m Roll	Each

Opsite Flexigrid Film Dressing Smith & Nephew

Smith & Nephew OpSite Flexigrid - Moisture Vapour Permeable Adhesive Film Dressing

OpSite Flexigrid is a transparent, adhesive film, with a unique wound measurement grid. The film is conformable and extensible. It is clinically proven to provide moist wound healing for superficial wounds; it is also ideal for use as a secondary dressing and to help prevent pressure sores.



Features and Benefits:

- ✓ Faster, safer, less painful wound healing
- ✓ Superior film with polyacrylate adhesive
- ✓ Waterproof film
- ✓ Simple, effective application technique
- ✓ Adaptable to awkward areas

Medent Code	Mfr Code	Description	Unit of Sale
DWTDOS1	4628	Smith & Nephew OpSite Flexigrid 6 x 7cm Film Dressing	Box 100
DWTDOS4	4630	Smith & Nephew OpSite Flexigrid 10 x 12cm Film Dressing	Box 50
DWTDOS6	4629	Smith & Nephew OpSite Flexigrid 10 x 12cm Film Dressing	Box 10



Opsite IV 3000 Dressing Smith & Nephew

Smith & Nephew OpSite Flexigrid IV3000* Moisture Responsive Cannula Dressing. **IV3000** is a moisture responsive transparent film dressing specifically designed to meet the needs of catheter fixation and to keep the catheter site dry.

Independent clinical experts have proven **IV3000** to offer safe, secure and cost effective catheter fixation and believe it to be a technological advance for IV site care.



Features and Benefits:

- ✓ Moisture responsive high MVTR 1 film
- ✓ Low allergy grid pattern adhesive
- ✓ Transparent, waterproof and conformable
- ✓ 1 Hand delivery system
- ✓ Film acts as a barrier to bacteria including MRSA2



Medent Code	Mfr Code	Description	Unit of Sale
DWTDOSIV	4925	Smith & Nephew Opsite IV 3000 Dressing 10cm x 14cm	Box 10

Tegaderm Film Dressing 3M

3M Tegaderm transparent film dressings designed for protecting skin and wound sites. Frame style allows customisation of shape and size to fit any site. Unique sacral shape is designed for challenging and difficult-to-dress areas.

Waterproof dressing provides a bacterial and viral barrier to outside contaminants and allows the patient to shower. Breathable film allows moisture vapor and oxygen exchange while providing a moist wound environment for enhanced healing.



Medent Code	Mfr Code	Description	Unit of Sale
DWTD1	1622W	3M Tegaderm 4.4 x 4.4cm Film Dressing	Box 100
DWTD3	1624W	3M Tegaderm 6 x 7cm Film Dressing	Box 100
DWTD4	1626W	3M Tegaderm 10 x 12cm Film Dressing	Box 50
DWTD5	1627	3M Tegaderm 10 x 25cm Film Dressing	Box 20
DWTD6	1628	3M Tegaderm 15 x 20cm Film Dressing	Box 10
DWTD7	1629	3M Tegaderm 20x 30cm Film Dressing	Box 10
DWTD8	9543HP	3M Tegaderm Sacral Film Dressing	Box 12

Tegaderm IV Film Dressing 3M

3M Tegaderm IV transparent film dressings designed to secure the catheter in challenging I.V. applications. Many with notches and borders reinforced with a soft cloth fabric, designed to reduce edge lift and to seal around the insertion site.

Tegaderm™ I.V. Film has a 'frame' delivery system. Sterile soft cloth tape strips are also provided with many dressings for anchoring hubs, lumens and tubing. Dressings are waterproof and latex-free; Tegaderm™ Plus Transparent Dressing packaging contains latex. Available in a variety of sizes and shapes to address multiple applications. Tegaderm™ Plus contains an iodophor in the dressing.




Medent Code	Mfr Code	Description	Unit of Sale
DWTD2	1623W	3M Tegaderm IV Film Dressing 6 x 7cm	Box 100
DWTD8	1633	3M Tegaderm IV Film Dressing 7 x 8.5cm	Box 100
DWTD9	1610	3M Tegaderm IV Paediatric Film Dressing 5 x 5.7cm	Box 100
DWTD10	1635	3M Tegaderm IV Film Dressing 8.5 x 10.5cm	Box 50

Tegaderm Film Roll 3M

3M Tegaderm Roll provides maximum versatility and the affordable choice for skin protection and waterproof cover dressings.

Waterproof backing provides a bacterial and viral barrier to outside contaminants and allows the patient to shower. The roll can be cut and customised to desired length. It is transparent to allow for easy monitoring and easy-to-use delivery reduces application waste.

Laboratory testing has proven 3M™ Tegaderm™ and Tegaderm™ HP Transparent Dressings provide a viral (HIV-1 and HBV) barrier while dressings remain intact without leakage. Dimensions of dressing listed, includes hydrocolloid and film backing.



Medent Code	Mfr Code	Description	Unit of Sale
DWTD1	16002	3M Tegaderm Film Roll 5cm x 10m	Each
DWTD2	16004	3M Tegaderm Film Roll 10cm x 10m	Each

CATEGORY: DRESSING WOUND – FOAM

Allevyn Hydro Cellular Dressing Smith & Nephew

Smith & Nephew Allevyn Non Adhesive Hydrocellular Foam Dressing.

Allevyn Heel dressing has been designed for a superior anatomical and comfortable fit to enhance management of wounds on the heel. One piece shaped construction. Protective padding

Features and Benefits:

- ✓ Unique trilaminar structure
- ✓ Non adherent wound contact layer
- ✓ Highly absorbent central hydrocellular layer
- ✓ Waterproof outer film layer
- ✓ Non adhesive, easy to use
- ✓ Helps create and maintain a moist wound environment
- ✓ Soft and cushioning



Medent Code	Mfr Code	Description	Unit of Sale
DWHAL1	66007643	Smith & Nephew Allevyn Hydrocellular Foam Non Adhesive Dressing 5 x 5 cm	Box 10
DWHAL2	66007637	Smith & Nephew Allevyn Hydrocellular Foam Non Adhesive Dressing 10 x 10 cm	Box 10
DWHALH	66007630	Smith & Nephew Allevyn Hydrocellular Foam Non Adhesive Dressing 10 x 10 cm	

Allevyn Adhesive Hydro Cellular Dressing Smith & Nephew

Smith & Nephew Allevyn Adhesive is an economically efficient dressing specifically designed to manage chronic or acute exuding wounds.

Features and Benefits:

- ✓ Unique triple action technology helps to maintain an optimal balance in fluid creating a moist wound environment for healing
- ✓ Low allergy adhesive helps to keep the dressing securely in place without adhering to the wound
- ✓ Provides an effective barrier function to exogenous bacteria as well as helping to prevent fluid and bacterial strike through
- ✓ Easy to apply and remove
- ✓ Conforms to body contours and is very comfortable for patients
- ✓ Available in a wide range of sizes which can be cut to fit awkward areas such as elbows and heels
- ✓ Sacral sizes available which combine the same high-absorbency and secure self-adhesion in an anatomically-shaped dressing




Medent Code	Mfr Code	Description	Unit of Sale
DWHALA1	66000044	Smith & Nephew Allevyn Adhesive Hydrocellular Foam Dressing 12.5 x 12.5 cm	Box 10
DWHALA2	66000043	Smith & Nephew Allevyn Adhesive Hydrocellular Foam Dressing 7.5 x 7.5 cm	Box 10
DWHASS	66000700	Smith & Nephew Allevyn Adhesive Sacrum Small Hydrocellular Foam Dressing 17 x 17 cm	Box 10

Allevyn Cavity Circular Foam Hydro Dressing Smith & Nephew

Smith & Nephew Allevyn Cavity provides Allevyn hydrocellular technology in a unique 3 dimensional structure for effective management of deep wounds.

Features and Benefits:

- ✓ Highly Absorbent
- ✓ Highly Conformable to Wound Contours
- ✓ Easy to use and apply
- ✓ Superior Non-Adherency
- ✓ Promotes Faster Healing
- ✓ Core of highly absorbent hydrocellular foam chips
- ✓ Outer non-adherent perforated film
- ✓ Range of sizes and shapes
- ✓ Soft pliable construction



Medent Code	Mfr Code	Description	Unit of Sale
DWFAC	66007326	Smith & Nephew Allevyn Cavity 5cm Circular Foam Hydro Dressing	Box 10
DWFAC10	66007327	Smith & Nephew Allevyn Cavity 10cm Circular Foam Hydro Dressing	Box 5

Allevyn Gentle Border Adhesive Smith & Nephew

Smith & Nephew Allevyn Gentle Border minimises trauma to the wound and pain to the patient during dressing changes.

Features and Benefits:

- ✓ Triple action technology - absorbs - retains - transpires
- ✓ Silicone gel adhesive flows gently into the crevices of the skin holding the dressing securely in place



ALLEVYN[®] GENTLE BORDER Silicone Gel Adhesive Dressing

Sensitive Skin & Skin Tears

Easy to apply
Gentle to remove
Comfortable
Absorbent
Sterile



2 Dressings (7.5cm x 7.5cm)

Medent Code	Mfr Code	Description	Unit of Sale
DWHALGB	66800269	Smith & Nephew Allevyn Gentle Border 7.5 x 7.5 Adhesive	Box 10
DWHALGB1	66800272	Smith & Nephew Allevyn Gentle Border 12.5 x 12.5 Adhesive	Box 10

Allevyn Thin Hydro Colloid Dressing Smith & Nephew

Smith & Nephew Allevyn Thin is a highly conformable, semi-permeable, semi-transparent wound dressing.

Features and Benefits:

- ✓ Sterile
- ✓ Self adhesive
- ✓ Flexible and highly conformable
- ✓ Maintains a moist wound environment
- ✓ Waterproof and bacteriaproof
- ✓ Permeable to oxygen and water vapour
- ✓ Does not adhere to the wound
- ✓ Sticks well to dry, intact skin
- ✓ Semi transparent film



Medent Code	Mfr Code	Description	Unit of Sale
DWHTA	66047576	Smith & Nephew Allevyn Thin Hydro Colloid Dressing 5 x 6cm	Box 10
DWHTA1	66047578	Smith & Nephew Allevyn Thin Hydro Colloid Dressing 10 x 10cm	Box 5

Convatec Versiva XC Adhesive Dressing

Convatec Versiva® XC® dressing with Hydrofiber® Technology has a wound contact layer made from Hydrofiber® Technology, which is unique to Convatec.

Features and Benefits:

- ✓ Locking in wound exudate
- ✓ Rapidly turns to a clear gel on contact with fluid
- ✓ Absorbs and retains fluid, preventing lateral wicking
- ✓ Traps harmful bacteria and enzymes contained in exudate
- ✓ Contouring to the wound bed
- ✓ Minimizes dead space where harmful bacteria can thrive
- ✓ The gel covers the wound bed
- ✓ Responding to wound conditions
- ✓ Controlled MVTR
- ✓ Maintains a moist wound healing environment



Medent Code	Mfr Code	Description	Unit of Sale
DWER1	410607	Convatec Versiva xC Non Adhesive Dressing 11cm x 11cm	Box 10
DWER2	410606	Convatec Versiva xC Non Adhesive Dressing 7.5cm x 7.5cm	Box 10
DWER3	410609	Convatec Versiva xC Adhesive Dressing 10cm x 10cm	Box 10
DWER4	410612	Convatec Versiva xC Adhesive Heel Dressing 18.5cm x 20.5cm	Box 5
DWER5	410613	Convatec Versiva xC Adhesive Sacral Dressing 21cm x 25cm	Box 5
DWER6	410610	Convatec Versiva xC Adhesive Dressing 14cm x 14cm	Box 10

Copa Hydrophilic Foam Dressings

Kendall COPA ultra-soft foam dressing is an ultra-soft hydrophilic foam that protects and cushions wounds. Small pores decrease the chance of tissue migration into the foam dressing.



Medent Code	Mfr Code	Description	Unit of Sale
DWFS1	55544	Kendall Copa Hydrophilic Foam Dressings 10cm x 10cm	Box 10

Molnlycke Lyofoam

MOLNLYCKE Lyofoam® Polyurethane Foam Dressing. Soft, conformable and comfortable Lyofoam® dressing is designed for the management of moderately to highly exuding wounds.

Its high absorbency and fluid handling capacity combined with the reduced risk of maceration and leakage provide a longer wear time.



Features and Benefits:

- ✓ High absorbency and fluid handling capacity
- ✓ Reduced risk of maceration, works under compression
- ✓ High quality and cost effective
- ✓ Waterproof backing film acts as a barrier to bacterial and viral penetration

Medent Code	Mfr Code	Description	Unit of Sale
DWFLX1	STN 603088	Molnlycke Lyofoam 10 x 10cm Extra Adhesive	Box 10
DWFLXA1	STN 603086	Molnlycke Lyofoam 15 x 15cm Extra Adhesive	Box 10
DWFL1	STN 603073	Molnlycke Lyofoam 9 x 9cm Extra Adhesive	Box 10
DWFL3	603201	Molnlycke Lyofoam Max 10 x 10cm	Box 10
DWFL5	STN603090	Molnlycke Lyofoam 20 x 15cm Extra Adhesive	Box 10
DWFL6	603014	Molnlycke Lyofoam T 9 x 6.5cm Extra Adhesive	Box 10

Molnlycke Mepilex Foam

Mepilex® Foam is a soft and conformable foam dressing. Suitable for a wide range of acute and chronic wounds like venous leg ulcers, pressure ulcers or diabetic ulcers. It minimises pain and wound or skin damage at dressing change, and can be cut to fit making it ideal for difficult to dress areas.

Mepilex® Lite is a thin foam dressing for low exuding acute and chronic wounds. Through Safetac® technology, it minimises pain and wound or skin damage at dressing change. It is very soft and conformable and clinically proven to significantly reduce the severity of skin reactions from radiotherapy treatment.



Medent Code	Mfr Code	Description	Unit of Sale
DWFM1	294100	Molnlycke Mepilex Foam 10cm x 10cm	Box 5
DWFM2	294200	Molnlycke Mepilex Foam 10cm x 20cm	Box 5
DWFM3	294300	Molnlycke Mepilex Foam 15cm x 15cm	Box 5
DWFMAG	287110	Molnlycke Mepilex Foam Ag 10cm x 10cm	Box 5
DWFML1	284000	Molnlycke Mepilex Foam Lite 6cm x 8.5cm	Box 5
DWFML2	284100	Molnlycke Mepilex Foam Lite 10cm x 10cm	Box 5
DWFML3	284300	Molnlycke Mepilex Foam Lite 15cm x 15cm	Box 5

Molnlycke Mepilex Foam Border

Mepilex® Border - all-in-one foam dressing with Safetac® technology. All-in-one foam dressing proven to minimise patient pain and trauma to the wound and surrounding skin.

Designed to create an optimal healing environment for moderate to high exuding wounds such as pressure ulcers, leg and foot ulcers, traumatic wounds, and other secondary healing wounds.

Mepilex® Border Lite - thin all-in-one foam dressing with Safetac® technology. Thin and conformable to dress wounds on feet and toes, hands and fingers. Very thin and discrete to allow patients to live a normal life. Shower proof – with a film backing providing a barrier to virus and bacteria (of size >25nm).

Mepilex® Border



Safetac®
TECHNOLOGY



Medent Code	Mfr Code	Description	Unit of Sale
DWFMB1	295200	Molnlycke Mepilex Foam Border 7.5cm x 7.5cm	Box 5
DWFMB2	295300	Molnlycke Mepilex Foam Border 10cm x 10cm	Box 5
DWFMB3	295400	Molnlycke Mepilex Foam Border 15cm x 15cm	Box 5
DWFMBL1	281000	Molnlycke Mepilex Foam Border Lite 4cm x 5cm	Box 10
DWFMBL2	281200	Molnlycke Mepilex Foam Border Lite 7.5cm x 7.5cm	Box 5
DWFMBL3	281300	Molnlycke Mepilex Foam Border Lite 10cm x 10cm	Box 5
DWFMBL4	281500	Molnlycke Mepilex Foam Border Lite 15cm x 15cm	Box 5
DWFMBL5	281100	Molnlycke Mepilex Foam Border Lite 5cm x 12.5cm	Box 5
DWFMBP1	495100	Molnlycke Mepilex Post Op Border 6cm x 8cm	Box 5

Tegaderm Foam Adhesive Dressing Oval 3M

3M Tegaderm Foam Adhesive Dressing. High Performance Foam. Full line of highly absorbent, breathable adhesive foam dressings. Provide total fluid management by a combination of fast wicking, high absorbency and breathability.

The innovative spoke delivery system allows fast, easy application for wounds over body contours. Available in a wide variety of sizes and shapes to meet your clinical needs.

Oval and square dressings are constructed of a conformable absorbent, poly urethane foam pad, an additional absorbent nonwoven layer and a border of transparent film with moist skin adhesive technology.

The mini wrap dressing is constructed of a conformable, absorbent, poly urethane foam pad with a highly breathable film backing, reinforced with soft cloth tape.



Medent Code	Mfr Code	Description	Unit of Sale
DWHFA1	90610	3M Tegaderm Foam Adhesive Dressing 8.8 x 8.8cm Oval	Box 10
DWHFA2	90611	3M Tegaderm Foam Adhesive Dressing 10 x 11cm Oval	Box 10
DWHFA3	90614	3M Tegaderm Foam Adhesive Dressing 6.9 x 7.6cm Mini Oval	Box 10
DWHFA10	90612	3M Tegaderm Foam Adhesive Dressing 14.3 x 14.3cm Oval	Box 10
DWFA3M	90619	3M Tegaderm Foam Adhesive Dressing Heel Elbow	Box 5

Features and Benefits:

- ✓ Reduced potential for skin maceration
- ✓ Unique moist skin adhesive provides adherence to dry and moist skin conditions
- ✓ Soft, conformable - patient comfort
- ✓ Oval design conforms to contours
- ✓ Mini wrap design ideal for toes, feet, noses, elbows and chins

Tegaderm Foam Non Adhesive Dressing 3M

3M Tegaderm Foam Non Adhesive Dressing. The wound dressing is a fast wicking, non-swelling poly urethane foam with a highly breathable film backing which prevents wound exudate strike-through and is a barrier to outside contamination while the dressing remains intact without leakage.

Heavily draining wounds often result in skin maceration, frequent dressing changes and patient discomfort. Tegaderm foam dressing provides total fluid management by a combination of fast wicking, high absorbency and breathability. The edges of the dressing must be bordered with 3M Tegaderm Transparent Dressing for barrier protection against bacteria and viruses.



Medent Code	Mfr Code	Description	Unit of Sale
DWFNAD1	90600	3M Tegaderm Foam Non Adhesive Dressing 5 x 5cm	Box 10
DWFNAD2	90601	3M Tegaderm Foam Non Adhesive Dressing 10 x 10cm	Box 10
DWHFRNA	90605	3M Tegaderm Foam Roll Non Adhesive Dressing 10 x 60cm	Each

Molnlycke Mepilex Foam Border

Mepilex® Border - all-in-one foam dressing with Safetac® technology. All-in-one foam dressing proven to minimise patient pain and trauma to the wound and surrounding skin.

Designed to create an optimal healing environment for moderate to high exuding wounds such as pressure ulcers, leg and foot ulcers, traumatic wounds, and other secondary healing wounds.

Mepilex® Border Lite - thin all-in-one foam dressing with Safetac® technology. Thin and conformable to dress wounds on feet and toes, hands and fingers. Very thin and discrete to allow patients to live a normal life. Shower proof – with a film backing providing a barrier to virus and bacteria (of size >25nm).

Mepilex® Border



Safetac
TECHNOLOGY

Medent Code	Description	Unit of Sale
DWFMB1	Molnlycke Mepilex Foam Border 7.5cm x 7.5cm	Box 5
DWFMB2	Molnlycke Mepilex Foam Border 10cm x 10cm	Box 5
DWFMB3	Molnlycke Mepilex Foam Border 15cm x 15cm	Box 5
DWFMBL1	Molnlycke Mepilex Foam Border Lite 4cm x 5cm	Box 10
DWFMBL2	Molnlycke Mepilex Foam Border Lite 7.5cm x 7.5cm	Box 5
DWFMBL3	Molnlycke Mepilex Foam Border Lite 10cm x 10cm	Box 5
DWFMBL4	Molnlycke Mepilex Foam Border Lite 15cm x 15cm	Box 5
DWFMBL5	Molnlycke Mepilex Foam Border Lite 5cm x 12.5cm	Box 5
DWFMBP1	Molnlycke Mepilex Post Op Border 6cm x 8cm	Box 5



CATEGORY: DRESSING WOUND – GELS/HYDROGELS/SILICONE/POWDER

Activon Manuka Honey Tube

Activon Tube contains 100% Medical Grade Manuka honey with no additives from New Zealand. The Manuka honey is filtered and sterilised, making it suitable for application to wounds. Manuka honey will de-bride and de-slough, eliminate odours and provides a moist wound healing environment.

Activon is ideal for debriding necrotic tissue, or for topping up dressings where the honey has been washed away by exudate. Great for using in cavities, just wash away with saline solution.



Features and Benefits:

- ✓ Antibacterial, kills harmful bacteria
- ✓ Osmotic effect, draws harmful tissue away from the wound bed
- ✓ Enhanced healing, anti-inflammatory
- ✓ Odour reduction, eliminates odours without masking them
- ✓ Conformable, non-adherent
- ✓ Maintains the ideal moist wound healing environment

Medent Code	Mfr Code	Description	Unit of Sale
DWHON	CR3830	Activon Manuka Honey Tube 25g	Each

Tricotex Wound Layering Dressing Smith & Nephew

Smith & Nephew Tricotex dressings are constructed from knitted viscose rayon and are designed to act as an interface between ulcerating or granulating wounds and conventional absorptive dressings.

The open knitted construction allows exudate to pass freely from the wound site. Disruption of the wound and its healing tissue is avoided by leaving Tricotex in place while the absorbent pads are changed as required.

Being non-medicated, Tricotex dressings are unlikely to elicit any tissue reaction so are ideal for use on their own but can be used in conjunction with a medicament of choice.

Tricotex dressings can be removed from the wound at any stage by gently stretching diagonal corners and lifting upwards from the edges view of the wound



Medent Code	Mfr Code	Description	Unit of Sale
DWWLT	7111	Smith & Nephew Tricotex Wound Layering Dressing 9.5 x 9.5cm	Box 50

Aqueous Cream (BP)

Features and Benefits:

- ✓ A popular cream for healing
- ✓ Soothes and moisturises skin
- ✓ Ideal as a hand and body barrier lotion
- ✓ Also often used as a soap substitute



Medent Code	Mfr Code	Description	Unit of Sale
DWAQCR	2072254	Aqueous Cream (BP) 500g	Each

Cica Care Silicone Gel Sheet Smith & Nephew

Smith & Nephew Cica-Care is a self-adhesive gel sheet that is medically proven to be up to 90% effective in the improvement of red, dark or raised scars. It has already been used by the medical profession to treat over one million people worldwide.

Cica-Care is durable, comfortable and re-usable, ideal for day and night use. It is also safe to use on children and adults making it suitable for the whole family. Cica-Care can even be used on scars up to 20 years old. Within 2-4 months, scars you thought were permanent can be treated and their appearance improved - permanently.



Medent Code	Mfr Code	Description	Unit of Sale
DWTGJCC1	66250706	Smith & Nephew Cica Care Silicone Gel Sheet 12x15cm	Each

Convatec Duoderm Gel Tube

DuoDERM Hydroactive Gel is a clear, preservative-free, viscous hydrogel. It is indicated for the management of partial and full-thickness wounds, ideally as a filler for dry cavity wounds to provide a moist healing environment. This product contains latex.



Medent Code	Mfr Code	Description	Unit of Sale
DWG2	187990	Convatec Duoderm Gel 15g Tube	Each
DWG3	187987	Convatec Duoderm Gel 30g Tube	Each
DWG3A	187987	Convatec Duoderm Gel 30g Tube	Box 3

Convatec Saf Gel Tube

SAF-Gel Hydrating Dermal Wound Dressing has an alginate containing formula, indicated for use on chronic and acute wounds. Its multiple-occasion, single-patient use package helps eliminate waste.



Medent Code	Mfr Code	Description	Unit of Sale
DWGSA	145730	Convatec Saf Gel 85g Tube	Each

Durafiber Absorbent Non Woven Gelling Fibre Smith & Nephew

Smith & Nephew Durafiber Absorbent Non Woven Gelling Fibre. A strong, gelling, fibre dressing, with the benefit of easy one piece removal. Conforms to the wound bed, has a high absorbent capacity, and locks in fluid and harmful bacteria it may contain.



Medent Code	Mfr Code	Description	Unit of Sale
DWDF5	66800559	Smith &Nephew Durafiber Absorbent Non Woven Gelling Fibre 5 x 5cm	Box 10
DWDF10	66800560	Smith &Nephew Durafiber Absorbent Non Woven Gelling Fibre 10 x 10cm	Box 10
DWDFR	66800563	Smith &Nephew Durafiber Absorbent Non Woven Gelling Rope 2 x 45cm	Box 5

Duratouch Dressings Silicon Smith & Nephew

Smith & Nephew Duratouch Silicon Dressing. An open mesh, soft silicone wound contact layer that allows exudate passage into a secondary dressing.

The soft silicone provides gentle adhesion to the surrounding skin as well as adhering to the absorbent layer. It contributes to minimal trauma on removal and helps to protect the wound bed and surrounding skin.



Medent Code	Mfr Code	Description	Unit of Sale
DWSDUR1	66801156	Smith &Nephew Duratouch Dressings Silicon 5cm x 7cm	Box 5
DWSDUR3	66801158	Smith &Nephew Duratouch Dressings Silicon 12cm x 15cm	Box 5

Flamazine Cream (Silver Sulfadiazine 1%) Pom Smith & Nephew

The mechanism of silver sulfadiazine's antibacterial action has not been fully elucidated. After exposure to the drug, structural changes in the bacterial cell membrane occur, including distortion and enlargement of the cell and a weakening of the cell wall membrane. This is accompanied by reduced viability in sensitive strains due to interference with macromolecular synthesis.

The sulfadiazine moiety also provides a bacteriostatic action against sensitive organisms. In adults, up to 10% of the sulfadiazine may be absorbed and 60-85% of the absorbed amount is excreted in the urine. In children with 13% body surface area burns, the urinary sulfadiazine concentration was 31.8 mg/L.



Medent Code	Mfr Code	Description	Unit of Sale
DWGFL	66800531	Smith & Nephew Flamazine Cream 50g (Silver Sulfadiazine 1%) Pom	Each

Gel Burnaid Tube

Burnaid Burn Gel Tube. First Aiders Choice. Gel for burns, scalds and sun burn.

Features and Benefits:

- ✓ Relieves the pain of small spot burns
- ✓ Cools
- ✓ Water soluble gel
- ✓ Can be rinsed off
- ✓ Non-toxic
- ✓ Sterile



Medent Code	Mfr Code	Description	Unit of Sale
DWGELB	PHK740	Burnaid Gel 25gram Tube	Each
DWGELB1	751588	Burnaid Gel 50 Gram Tube	Each

Gel Intrasite Wound Smith & Nephew

Smith & Nephew Gel Intrasite Wound. The unique Applipak application system containing Intrasite Gel is available in a range of different sizes. Its unique iconic shape was specifically designed to provide directable and controlled application of the gel.



Medent Code	Mfr Code	Description	Unit of Sale
DWGINT1	7308	Smith & Nephew Gel Intrasite Wound 8g	Box 10
DWGINT1E	7308	Smith & Nephew Gel Intrasite Wound 8g Singles	Each
DWGINT2	7311	Smith & Nephew Gel Intrasite Wound 15g	Box 10
DWGINT2E	7311	Smith & Nephew Gel Intrasite Wound 15g Singles	Each

Iodosorb Ointment Tubes Smith & Nephew

Smith & Nephew Iodosorb Ointment Tubes. Iodosorb removes barriers to healing by its dual action antimicrobial and desloughing properties. The broad spectrum antimicrobial action is provided by the sustained release of iodine and the desloughing action is provided by the unique cadexomer matrix.

By effectively removing the barriers of bacteria, slough and debris and excess exudate, Iodosorb has been shown to accelerate wound healing in leg ulcers and is effective for the treatment of infection in chronic exuding wounds.



Medent Code	Mfr Code	Description	Unit of Sale
DWGIO1	66051240	Smith & Nephew Iodosorb Ointment Tubes 10g	Box 4
DWGIO2	66051240	Smith & Nephew Iodosorb Ointment Tubes 10g Singles	Each

Molnlycke Mepiform

Molnlycke Mepiform. Everyday scar care dressing. Mepiform® is a thin and discreet dressing, designed for the management of both old and new hypertrophic and keloid scars. It can also be used as a prophylactic therapy on closed wounds for prevention of hypertrophic or keloid scarring. With the patented Safetac® technology, it also minimises pain during dressing changes.

Features and Benefits:

- ✓ Minimises pain and trauma during dressing changes
- ✓ Clinically proven to reduce scarring and reduce scars
- ✓ Conforms well to body contours, no need for extra fixation
- ✓ Can be worn during daily activities and be reapplied



Medent Code	Mfr Code	Description	Unit of Sale
DWTGSF1	293200	Molnlycke Mepiform 5cm x 7.5cm	Box 5

Molnlycke Mepitel Soft Silicone

Molnlycke Mepitel Soft Silicone wound contact layer is a gentle, yet effective wound contact layer with the unique Safetac® technology. The open mesh design enables good transfer of exudate to a secondary dressing and easy delivery of topical treatments. Mepitel is clinically proven to heal scalds faster and with less pain and cost.



Features and Benefits:

- ✓ Good transparency lets you instantly assess healing progress without removing the contact layer
- ✓ Does not leave residues and maintains functional qualities over time
- ✓ True undisturbed healing – can remain in place for up to 14 days
- ✓ Can cost-effectively be used in a wide variety of indications
- ✓ Safetac technology for minimal tissue trauma and patient stress and pain

Medent Code	Mfr Code	Description	Unit of Sale
DWTGSM1	290510	Molnlycke Mepitel Soft Silicone 5cm x 7.5cm	Box 10
DWTGSM2	290710	Molnlycke Mepitel Soft Silicone 7.5cm x 10cm	Box 10
DWTGSM3	291010	Molnlycke Mepitel Soft Silicone 10cm x 18cm	Box 10

Molnlycke Mesalt Debriders

Molnlycke Mesalt Debriders® supports the cleansing of moderately to heavily discharging wounds including infected wounds.

Features and Benefits:

- ✓ Absorbs exudate, bacteria and necrotic material
- ✓ Effectively supports the cleansing of the wounds
- ✓ Easy application

Medent Code	Mfr Code	Description	Unit of Sale
DWSALT	285780	Molnlycke Mesalt Debriders 7.5cm x 7.5cm	Box 30



Opsite Flexifix Gentle Roll Film Dressing with Silicone Gel Smith & Nephew

Smith & Nephew Opsite Flexifix Gentle Film Roll with Silicone Gel combines breathability, transparency and a seven day wear time to offer a versatile silicone film roll for patients with fragile or at-risk skin.

It removes cleanly, without disrupting fragile skin layers or causing patients undue pain. This dressing combines all the benefits of Opsite film, with a gentle silicone gel adhesive that enables the dressing to be lifted and repositioned without losing adherent properties.

The film is transparent, shower and bacteria proof, allowing monitoring of the site under the film without lifting the dressing. The low friction film helps to protect the skin from the formation of pressure injuries.



Medent Code	Mfr Code	Description	Unit of Sale
DWTDOSG1	66801195	Smith & Nephew Opsite Flexifix Gentle Film Roll with Silicon Gel 2.5cm x 5m	Each
DWTDOSG2	66801196	Smith & Nephew Opsite Flexifix Gentle Film Roll with Silicon Gel 5cm x 5m	Each
DWTDOSG3	66801197	Smith & Nephew Opsite Flexifix Gentle Film Roll with Silicon Gel 10cm x 5m	Each

Features and Benefits:

- ✓ Ideal for patients with fragile skin
- ✓ Minimizes the risk of pain or discomfort on removal
- ✓ Transparent, shower and bacteria proof
- ✓ Cut to required size

Purilon Hydrogel Tube

Purilon Hydrogel is a gentle hydrogel consisting of natural ingredients without any additives. Used for wounds that need effective and gentle debridement of necrotic tissue thereby paving the way for faster wound healing.

Features and Benefits:

- ✓ Unique match of hydrating and absorbing properties
- ✓ Easy to apply and remove
- ✓ Minimises leakage and maceration
- ✓ Hydrates necrotic tissue and absorbs excess exudate, slough and debris from the wound



Medent Code	Mfr Code	Description	Unit of Sale
DWGPUR	3900	Purilon Hydrogel 15g Tube	Box 10


Silflex Dressings Silicon

Silflex Dressings Silicon. A wound contact non-adherent dressing layer made from polyester mesh coated with Silflex Soft Silicone.

Silflex is designed to gently adhere to the skin surrounding a wound and not to the wound bed this minimises the pain and trauma associated with dressing change. The dressing allows the passage of exudate.

Features and Benefits:

- ✓ Non adherent (adheres only to dry surrounding skin not to a moist wound)
- ✓ Atraumatic soft silicone, reduced pain on dressing change
- ✓ Soft and conformable
- ✓ Allows passage of exudate
- ✓ Hydrophobic Silflex soft silicone coating



Medent Code	Mfr Code	Description	Unit of Sale
DWSFLEX	CR3922	Silflex Dressings Silicon 5 x 7cm	Box 10
DWSFLEX1	CR3923	Silflex Dressings Silicon 8 x 10cm	Box 10

Solosite Gel Tube Smith & Nephew

Smith & Nephew SoloSite Gel Preserved Multi-Use Hydrogel. SoloSite Gel creates and maintains an optimal moist wound environment as it desloughs and debrides. Unlike gels formulated with water miscible polymers.

SoloSite's water swellable polymer remains gel-like until saturated, so it stays in the wound longer and requires fewer dressing changes

SoloSite has been proven non-cytotoxic, non-irritating and non-sensitising. Its crystal-clear, physically stable formula doesn't change colour, enabling an unimpeded

Features and Benefits:

- ✓ Preserved for safe, multi-dose dispensing
- ✓ Convenient tube presentations, including the 20 gram
- ✓ Effective debrider and deslougher
- ✓ Promotes moist wound healing
- ✓ Cost-effective



Medent Code	Mfr Code	Description	Unit of Sale
DWGSOL20	36100614	Smith & Nephew SoloSite Gel Tube 20g	Box 10
DWGSOL20	36100614	Smith & Nephew SoloSite Gel Tube 20g Singles	Each
DWGSOL	449600	Smith & Nephew SoloSite Gel Tube 85g Singles	Each

CATEGORY: DRESSING WOUND – HYDROCOLLOIDS

Collodion Flex BP Liquid To Dry

Flexible Collodion. Clear or slightly opalescent, highly flammable. A syrupy liquid compounded of pyroxylin, ether, and alcohol, which dries to a transparent, tenacious film.

Used as a topical protectant, applied to the skin to close small wounds, abrasions, and cuts, to hold surgical dressings in place, and to keep medications in contact with the skin.




Medent Code	Mfr Code	Description	Unit of Sale
DWHCOL	PHF093	Collodion Flex BP Liquid To Dry 100ml	Box 10

Comfeel Mini Hydro Colloid Dressing

Comfeel Ulcer Hydrocolloid Dressing. Contains liquid-absorbing particles derived from natural cellulose in an elastic, self-adhesive mass, which is covered on one side by a semi-permeable, non-adhesive polyurethane film. The cellulose-like particles expand strongly when they absorb liquids.

Maintains a moist environment which softens necrotic tissue and thus promotes biological debridement. On absorbing wound exudate the dressing forms a yellowish gel which fills the entire wound cavity and protects granulation tissue and nerve endings.



Comfeel® Plus Ulcer Dressing. Faster wound healing by sealing and protection of the wound. Used as primary and secondary dressing. Can be used directly on low to moderately exuding wounds. Available in shapes to dress any part of the body. Calcium alginate has been added to the dressing to increase absorption.



Comfeel Plus Film Hydro Colloid Dressing. Transparent Dressing allows for wound inspection without dressing change. Ultra thin and flexible for sealing and healing low-to-low exuding wounds. Flexible to dress any part of the body.

Medent Code	Mfr Code	Description	Unit of Sale
DWHCU1	3213	Comfeel Ulcer 10 x 10cm Hydro Colloid Dressing	Box 10
DWHCU15	3218	Comfeel Ulcer 15 x 15cm Hydro Colloid Dressing	Box 5
DWHCU20	3223	Comfeel Ulcer 20 x 20cm Hydro Colloid Dressing	Box 5
DWHCM	3233	Comfeel Mini 6 x 4cm Hydro Colloid Dressing	Box 30
DWHTCP	29103283	Comfeel Plus Contour 9cm x 11cm Hydro Colloid Dressing	Box 5
DWHPC	29103280	Comfeel Plus Contour 6cm x 8cm Hydro Colloid Dressing	Box 5
DWHTCPT	29103533	Comfeel Plus Film 10 x 10cm Hydro Colloid Dressing	Box 10
DWHTCT	29103530	Comfeel Plus Film 5 x 7cm Hydro Colloid Dressing	Box 10
DWHTCPT2	29103536	Comfeel Plus Film 9 x 14cm Hydro Colloid Dressing	Box 10
DWHTCPTS	29103585	Comfeel Plus Film Sacral 18 x 20cm Hydro Colloid Dressing	Box 5
DWHTCPT1	29103542	Comfeel Plus Film 15 x 20cm Hydro Colloid Dressing	Box 5
DWHTCPU2	29103110	Comfeel Plus Ulcer 10 x 10cm Hydro Colloid Dressing	Box 10
DWHTCPU15	29103115	Comfeel Plus Ulcer 15 x 15cm Hydro Colloid Dressing	Box 5
DWHTCPU20	29103120	Comfeel Plus Ulcer 20 x 20cm Hydro Colloid Dressing	Box 5
DWHTCPU1	29103146	Comfeel Plus Ulcer 6 x 4cm Hydro Colloid Dressing	Box 30

Convatec Duoderm Hydrocolloid Dressing


DuoDERM CGF Sterile Dressing is a hydrocolloid, moisture-retentive wound dressing used for partial and full-thickness wounds with exudate. It incorporates a unique Convatec hydrocolloid formulation that distinguishes it from other hydrocolloid dressings.

Indicated for chronic wounds - pressure ulcers (Stage I-IV), leg ulcers, acute wounds - traumatic wounds (minor abrasions, lacerations) partial thickness burns and donor sites. dressing forms a yellowish gel which fills the entire wound cavity and protects granulation tissue and nerve endings.

DuoDERM Signal Dressing is an adhesive, tapered edge, modern hydrocolloid wound dressing with a change indicator. The adhesive hydrocolloid matrix and the outer film on DuoDERM Signal Dressing are permeable by moisture vapour and provides an occlusive barrier against bacteria, viruses, dirt and water. Indicated for use on dermal ulcers including full-thickness wounds, pressure ulcers (II-IV), leg ulcers, superficial wounds, partial thickness burns and donor sites.

DuoDERM Extra Thin CGF Dressing is a hydrocolloid dressing indicated for the management of lightly exuding wounds.

DuoDERM CombiDERM ACD Sterile Dressing is a sterile, occlusive dressing with a hydrocolloid adhesive and central absorbent pad indicated for the management of low to moderately exuding chronic and acute wounds. It combines a non-adherent wound contact layer surrounded by an adhesive hydrocolloid border and an integral pad containing granules which absorb exudate and form a gel to lock fluid away.



Medent Code	Mfr Code	Description	Unit of Sale
DWHD1	187660	Convatec Duoderm CGF 10 x 10cm Hydrocolloid Dressing	Box 5
DWHD2	187661	Convatec Duoderm CGF 14.4 x 14.4cm Hydrocolloid Dressing	Box 5
DWHD3	187643	Convatec Duoderm CGF 15 x 20cm Hydrocolloid Dressing	Box 5
DWHD4	187662	Convatec Duoderm CGF 20 x 20cm Hydrocolloid Dressing	Box 5
DWHS1	403327	Convatec Duoderm SIGNAL 14 x 14cm Hydrocolloid Dressing	Box 5
DWHS2	410501	Convatec Duoderm SIGNAL Sacral 19.8 x 22.3cm Hydrocolloid Dressing	Box 5
DWHS3	403326	Convatec Duoderm SIGNAL 10 x 10cm Hydrocolloid Dressing	Box 5
DWHDX1	187901	Convatec Duoderm x-THIN 7.5 x 7.5cm Hydrocolloid Dressing	Box 20
DWHDX2	187955	Convatec Duoderm x-THIN 10 x 10cm Hydrocolloid Dressing	Box 10
DWHDX3	187957	Convatec Duoderm x-THIN 15 x 15cm Hydrocolloid Dressing	Box 10
DWHDX4	187900	Convatec Duoderm x-THIN 15 x 20cm Hydrocolloid Dressing	Box 20
DWHDX5	187961	Convatec Duoderm x-THIN 5 x 20cm Hydrocolloid Dressing	Box 10
DWHDX6	187932	Convatec Duoderm x-THIN SPOTS 4.4 x 3.8cm Hydrocolloid Dressing	Box 20
DWHCO2	187774	Convatec Hydro Colloid COMBIDERM 7.5 x 7.5cm	Box 10



Cutinova Hydro Dressing Smith & Nephew

Smith & Nephew Cutinova Hydro Dressing is a Hydro-selective, self adhesive sterile wound dressing consisting of two layers - a poly urethane gel matrix and a poly urethane top-film which is waterproof and acts as a bacterial barrier.



Medent Code	Mfr Code	Description	Unit of Sale
DWHCUT	66047443	Smith & Nephew Cutinova Hydro Dressing 10cm x 10cm	Box 5

Replicare Thin Hydro Colloid Dressing Smith & Nephew

Smith & Nephew REPLICARE Thin is indicated for exudate absorption and the management of partial to full thickness wounds. Some typical wounds are: Ulcers (venous, arterial, diabetic) Pressure Sores Donor Sites Surgical Incisions Surgical Excisions Burns (1st and 2nd degree)

Creation and maintenance of a moist wound environment. Moist wound environments have been established as optimal environments for the management of the wound. Provides physical separation between the wound and external environments to assist in preventing bacterial contamination of the wound.

Precautions REPLICARE Thin is not indicated for use on third degree burns.



Medent Code	Mfr Code	Description	Unit of Sale
DWHCR1	59484000	Smith & Nephew Replicare Thin 5cm x 7cm Hydro Colloid Dressing	Box 10

Tegaderm Hydrocolloid Thin Dressing 3M

3M Tegaderm Hydrocolloid Thin Dressing May be worn up to 7 days. Conformable dressing and easy to apply. Dressing maintains moist wound environment for enhanced healing. Outer film is 3M - Tegaderm™ Transparent Film Dressing technology which applies easily and eliminates need for taping edges.

Easy to clean, minimizes exudate leakage and is impermeable to liquids, bacteria and viruses. For minimally draining wounds.



Medent Code	Mfr Code	Description	Unit of Sale
DWTH1	90022	3M Tegaderm Hydrocolloid Thin Dressing 10 x 10cm	Box 5

CATEGORY: DRESSING WOUND – MEDICATED

Allewyn Life Dressing Medium Smith & Nephew

Smith & Nephew ALLEVYN Life Dressing. Multi-layered design incorporating hydrocellular foam, hyper-absorber lock away core and masking layer which has been designed for people and their every day life through.

Features and Benefits:

- ✓ Unique quadrilobe shape with a wide border designed to fit the contours of the human body more securely and allows patients to shower
- ✓ The protective masking layer minimises the visual impact of strikethrough (in-vitro), providing the patient with confidence that their dressing may not attract the negative attention of others
- ✓ Excellent fluid locking under pressure (in-vitro)
- ✓ Provides the patient with peace of mind leakage prevention (in-vitro)
- ✓ Common wound odour absorption when tested (in-vitro) on representative compounds
- ✓ ALLEVYN Life's multi-layered design provides cushioning and helps to spread pressure (in-vitro)
- ✓ Silicone wound contact layer offers minimal pain and discomfort on removal while allowing the dressing to be repositioned
- ✓ Change indicator optimises dressing change frequency and minimises wastage



Medent Code	Mfr Code	Description	Unit of Sale
DWHLIFS	66801067	Smith & Nephew ALLEVYN Life Dressing 10.3cm	Box 10
DWHLIFM	66801068	Smith & Nephew ALLEVYN Life Dressing 12.9cm	Box 10

Gelocast Zinc Impregnated Bandage (Stretched) Roll Smith & Nephew

Smith & Nephew Gelocast elastic is a cellulose-based, zinc gel, semi-rigid support and compression bandage primarily used for the treatment of venous diseases.

Consists of an elastic crepe bandage coated with zinc oxide gel. The zinc gel mass has a creamy consistency which allows rapid and neat application. After quickly drying, the bandage converts into a semi-rigid state and produces a low resting pressure and a high working pressure.

Gelocast elastic is supplied in sterile sachets which give a long shelf life.



Medent Code	Mfr Code	Description	Unit of Sale
DWPA5	45090-00	Smith & Nephew Gelocast Zinc Impregnated Bandage 8cm x 7m (Stretched) Roll	Each



Ichthopaste Bandage Smith & Nephew

Smith & Nephew Ichthopaste bandage is a Zinc oxide paste and ichthammol bandage. Used to assist the management of venous leg ulcers. Where venous insufficiency exists the paste bandage should be used under graduated compression bandaging.

For the management of chronic eczema/dermatitis where occlusion is indicated. Do not use in known cases of sensitivity or allergy to any of the ingredients.



Medent Code	Mfr Code	Description	Unit of Sale
DWPA2	4959	Smith & Nephew Ichthopaste Bandage	Each

Viscopaste Zinc Paste Bandage Smith & Nephew

Smith & Nephew Viscopaste zinc paste bandage is a general-purpose first line treatment for wet and dry leg ulcers. Soothes and protect the skin, particularly where eczema or dermatitis is present.



Medent Code	Mfr Code	Description	Unit of Sale
DWPA3	4948	Smith & Nephew Viscopaste Zinc Paste Bandage	Each

Zipzoc Zinc Oxide Stockings Smith & Nephew

Smith & Nephew Zipzoc Non-Preserved Zinc Oxide Impregnated Stocking is a sterile zinc oxide impregnated stocking. It is indicated for the treatment of chronic leg ulcers where chronic venous insufficiency exists.

Zipzoc can be used as a primary contact layer under compression bandaging or hosiery. It is a preservative free dressing thus reducing sensitivity reactions, and it contains white soft paraffin and liquid paraffin so it will not dry or adhere.



Medent Code	Mfr Code	Description	Unit of Sale
DWPA4	66000748	Smith & Nephew Zipzoc Zinc Oxide Stockings	Each

CATEGORY: DRESSING WOUND – ODOUR CONTROL

Carbonet Odour Control Dressing Smith & Nephew

Smith & Nephew Carbonet dressings have been designed to facilitate easier wound management and a reduction in nursing time in the treatment of malodorous and infected wounds.

Carbonet utilise the important, well-proven benefits of an activated charcoal layer, with the additional features of absorbency and non-adherency. The dressings create a favourable environment for wound management, whilst being cost effective and easy to use.

Features and Benefits:

- ✓ Activated charcoal layer
- ✓ Wound contact layer
- ✓ Absorbent padding layer
- ✓ Absorbs bacteria and offensive odours
- ✓ Flexible and conformable



Medent Code	Mfr Code	Description	Unit of Sale
DWODCAR1	7064	Smith & Nephew Carbonet 10 x 10cm Odour Control Dressing	Box 10

Convatec Carboflex Odour Control

Convatec CarboFLEX® Odor Control Dressing is a sterile, non-adhesive dressing indicated for shallow wounds or as a secondary dressing over a wound filler for deeper wounds.

It is a five layer dressing specifically designed to address the management problems associated with malodorous wounds.



Medent Code	Mfr Code	Description	Unit of Sale
DWOC1	403202	Convatec CarboFLEX® Odor Control Dressing 10 x 10 cm	Box 10

CATEGORY: DRESSING WOUND – PLASTERS (BANDAIDS)

Bandaid Plasters

Band-Aid adhesive bandage. Heals the hurt faster. Cover and protects cuts and scrapes.

Band-Aid with quilt aid technology. Heals cuts twice as fast. Cover and protects cuts and scrapes. Effective absorption. Non-stick pad technology.



Medent Code	Mfr Code	Description	Unit of Sale
DWPLTLS	6074923	Bandaid Plasters Princess Kids	Box 15
DWPLBB	6041709	Bandaid Plasters WaterProof Transformers Kids	Box 20
DWPLMM	6080783	Bandaid Plasters Waterproof Micky Mouse Kids	Box 15
DWPLSP3	6074765	Bandaid Plastic Strip Plasters with Quilt-Aid	Box 50



Covermed Injection Plasters Smith & Nephew

Smith & Nephew Covermed Injection Plasters. Low allergy adhesive non-woven dressing. Sensitive skin is liable to maceration, irritation, or allergic reactions when covered with an inappropriate adhesive dressing.

Covermed uses hypoallergenic adhesive with excellent skin tolerance, which can be gently removed from skin as well as from body hair.

Features and Benefits:

- ✓ Hypoallergenic acrylic adhesive
- ✓ Excellent skin tolerance
- ✓ Highly permeable to air and water vapour
- ✓ Soft and conformable non-woven substrate
- ✓ Absorbent wound pad laminated with polyethylene net
- ✓ Easy, painless dressing change



Medent Code	Mfr Code	Description	Unit of Sale
DWPLD1	72222-00	Smith & Nephew Injection Plasters Covermed 2cm x 4cm	Box 500

Coverplast First Aid Dressing Smith & Nephew

Smith & Nephew Coverplast. The reliable sterile adhesive wound dressing. For flexible comfort that stretches with every movement. Robust and durable that protect and cushion the wound, providing the practical answer to every day minor injuries.

For normal skin, Coverplast is a plaster that covers the wound and sticks firmly in place, with a wide selection of suitable shapes.



Coverplast **Classic Fabric** dressings are composed of an absorbent, low-adherent pad backed by a high quality stretch fabric



Coverplast **Barrier** dressings are waterproof barrier dressings.



Coverplast **Detectable** dressings are X-Ray and metal detectable, waterproof dressings

Coverplast **Visual** dressings are visually detectable, waterproof dressings

Medent Code	Mfr Code	Description	Unit of Sale
DWPLEL2	72601-05	Smith & Nephew Coverplast Classic Finger Tip Dressing Large 7.5 x 4.4	Box 100
DWPLEL4	72601-06	Smith & Nephew Coverplast Classic Knuckle Plasters 7.2 x 3.8	Box 100
DWPLEL5	72601-20	Smith & Nephew Classic Fingertip Dressing Small 5 x 4.4cm	Box 50
DWPLAIR2	72143-05	Smith & Nephew Coverplast Barrier Plasters 3.8 x 6.3	Box 100
DWPLAD	72143-09	Smith & Nephew Coverplast Blue Detectable Metal 7.2cm x 2.2cm	Box 100
DWPLAV	71300-06	Smith & Nephew Coverplast Blue Detectable Visual 7.2cm x 2.2cm	Box 100
DWIDE1	72601-03	Smith & Nephew Coverplast Classic 7.2 x 2.2cm Fabric Plasters	Box 100
DWIDE2	72601-00	Smith & Nephew Coverplast Classic 3.8 x 2.2cm Fabric Plasters	Box 100
DWIDE3	71151-03	Smith & Nephew Coverplast Classic 6cm x 5m Island Dressing Strip Tan	Box 100
DWIDE4	71151-02	Smith & Nephew Coverplast Classic 8cm x 5m Island Dressing Strip Tan	Box 100



Coverplast **Classic Dressing Strip** is coated on one side with a porous zinc oxide adhesive. The dressing conforms well to body contours and is comfortable for the patient.

Elastoplast Hypo Allergy Transparent Spots Smith & Nephew

Smith & Nephew Elastoplast Hypoallergenic transparent Spots. Ideal for sun spots. Water repellent. Kind to sensitive skin. Fresh coloured pad blends with skin tone.

Smith & Nephew Elastoplast Kids Plaster Cars Suitable for all types of smaller wounds. They protect kids' wounds from dirt and bacteria while sealing all around the wound. The material is skin friendly and the colours are safe for children. The non-stick wound pad protects and cushions the wound. The secure adhesion ensures that the plaster stays in place but is painless to remove.



Medent Code	Mfr Code	Description	Unit of Sale
DWPLHY	76334	Smith & Nephew Elastoplast Hypo Allergy Transparent Spots 22x22mm	Box 24
DWPLCARS	-	Smith & Nephew Elastoplast Plaster Cars	Box 16

Help-it Fabric Plaster

Help-it Fabric Plaster. Flexible, durable, breathable and wash resistant and feature a non-stick island dressing.

Medent Code	Mfr Code	Description	Unit of Sale
DWIDE6	PlaFab100	Help-it Fabric Plaster 72mm x 19mm	Box 100



Healband Plastic Spots/Plasters

Healband Plastic Spots/Strips -. Cushion comfort, non-adherent pad cushions the wound and allows for gentle removal

Protective pad keeps small wounds closed and provides protection from dirt contamination. Breathable ventilation holes promote air circulation to skin and wound.



Medent Code	Mfr Code	Description	Unit of Sale
DWPLEU2	-	Healband Plastic Plasters Sterile 72mm x 20mm	Box 100
DWPLSP2	62610516	Healband Plastic Spots - Sterile	Box 50
DWPLSP	-	Healband Plastic Spots - Non-Sterile	Box 50



Injection Pad IV Pressure Pad Defries

Unique application - Allow you to apply the pressure pad without contaminating the pad or adhesive.

Hypoallergenic adhesive



Medent Code	Mfr Code	Description	Unit of Sale
DWINJP	DEF929	Defries Injection Pad IV Pressure Pad	Box 100

Nexcare Plasters 3M

3M Nexcare Tattoo Plasters Cool Kids. Unique bandage shape provides superior waterproof protection. Designed for better adhesive seal around pad to provide the optimal healing environment that

3M Nexcare Waterproof Active Plaster Strips. Flexible foam cushions, stretches and bends. Sticks to damp skin. Cushions and protects cuts, scrapes and blisters. Latex free materials.

3M Nexcare Waterproof Clear Plaster. Easy, one-hand application. Breathable, ultra-thin material bends and flexes with movement. Unique diamond shape seals the area and minimises edge lift. Soft, absorbent pad won't stick to wounds.



Medent Code	Mfr Code	Description	Unit of Sale
DWPLTC	594-20	3M Nexcare Tattoo Plasters Cool Kids	Box 20
DWPLA1	512-30PB	3M Nexcare Waterproof Active Plaster Strips 26 x 76mm	Box 30
DWPLWP	586-20PB	3M Nexcare Waterproof Clear Plaster Medium 26 x 57mm	Box 20
DWPLWP1	588-20PB	3M Nexcare Waterproof Clear Plaster Assorted 26 x 57mm	Box 20

Nexcare Opticlude Regular Plasters (3M 1539)

3M Nexcare Opticlude Orthoptic Eye Patch. Designed for comfortable vision-occlusive therapy. Fits comfortably behind glasses, yet permits free movement of the eyelid. Opticlude™ has its own extra gentle hypoallergenic adhesive so no other tape is required to keep it in place.

- ✓ Gentle to skin, hypoallergenic adhesive
- ✓ Doctor recommended
- ✓ Sticks well, yet removes easily
- ✓ Highly breathable backing, latex free materials
- ✓ Ideal for treating amblyopia (lazy eye) and strabismus (crossed eyes)
- ✓ Comfortable to wear, absorbent non-stick pad



Medent Code	Mfr Code	Description	Unit of Sale
DWPLOP	PHH626	3M Nexcare Opticlude Orthoptic Eye Patch	Box 20

Plasters Blue Metal Detect Waterproof

Help-it Blue Metal and Visual Detection Plasters - Waterproof. Help-it fabric standard plasters are flexible, durable, breathable and wash resistant and feature a non-stick island dressing.



Medent Code	Mfr Code	Description	Unit of Sale
DWPLA3	BluPla100	Plasters Blue Metal Detect Waterproof 72mm x 25mm	Box 100

Wound Dressing

Surgical and Medical Wound Dressing. 100% Cotton. Sterile, single use only.



Medent Code	Mfr Code	Description	Unit of Sale
DWW13	AX1011051	Wound Dressing #13 Small	Each
DWW14	AX10110512	Wound Dressing #14 Medium	Each
DWW15	AX10110513	Wound Dressing #15 Large	Each

CATEGORY: DRESSING WOUND – PROCEDURE PACKS

Dialysis Procedure Tray Smith & Nephew

Medent Code	Mfr Code	Description	Unit of Sale
DWWPD	2908122	Smith & Nephew Propax Dialysis Procedure Tray	Box 114



Propax Lumps and Bumps Smith and Nephew

Medent Code	Mfr Code	Description	Unit of Sale
DWLUMP	2909377	Smith & Nephew Propax Lumps and Bumps	Box 36



Theatre Pack Sterile Dressing Pack

Medent Code	Mfr Code	Description	Unit of Sale
DWTP	BAM1524	Theatre Pack Sterile Dressing Pack	Box 100



Wound Dressing Pack 14 Sagepack Gauze

Medent Code	Mfr Code	Description	Unit of Sale
DWWPG14	11001014	Wound Dressing Pack 14 Sagepack Gauze	Box 180
DWWPG14E	11001014	Wound Dressing Pack 14 Sagepack Gauze	Each



Wound Dressing Pack Propax 6 Swabs Smith & Nephew

Medent Code	Mfr Code	Description	Unit of Sale
DWWP6	2908930-02	Smith & Nephew Wound Dressing Pack Propax 6 Swabs	Box 126
DWWP6E	2908930-02	Smith & Nephew Wound Dressing Pack Propax 6 Swabs Ctn	Each



Wound Dressing Pack Propax Gauze Balls Smith & Nephew

Medent Code	Mfr Code	Description	Unit of Sale
DWGBPCTN	2906569-02	Smith & Nephew Wound Dressing Pack Propax Gauze Balls	Box 126
DWGBPE	2906569-02	Smith & Nephew Wound Dressing Pack Propax Gauze Balls	Each



CATEGORY: DRESSING WOUND – RIBBON (PACKING) GAUZE

Intrasite Conformable Ribbon (Packing) Gauze Smith & Nephew

Smith & Nephew Intrasite Conformable Ribbon (Packing) Gauze. Combines the advantages of the gel with a non-woven dressing to aid the gentle packing of deep, shallow or open undermined wounds ensuring the gel remains in close contact with the entire wound



Medent Code	Mfr Code	Description	Unit of Sale
DWINTC	66000324	Smith & Nephew Intrasite Conformable Ribbon (Packing) Gauze 10cm x 10cm	Box 10

Ribbon Gauze Roll

Fast Edge Ribbon Gauze (FERG). A wide variety of sizes available for use in packing deep wounds. Available in x-ray detectable and plain. Synergy Health offers a range of specialised theatre products, available as sterile or non sterile surface.



Medent Code	Mfr Code	Description	Unit of Sale
DWVRG1	3820	Fast Edge Ribbon Gauze 1.25cm x 5m roll	Each
DWVRG2	21501	Fast Edge Ribbon Gauze Sterile 2.5cm x 5m roll	Each
DWVRG3	3822	Fast Edge Ribbon Gauze 5cm x 5m roll	Each
DWVRG5	21396	Fast Edge Ribbon Gauze Sterile 1.25cm x 1m roll	Each
DWVRG6	21397	Fast Edge Ribbon Gauze Sterile 2.5cm x 1m roll	Each



CATEGORY: DRESSING WOUND – SPRAYS/CREAMS/OINTMENTS

Barrier Cream Dimethicone 5%

healthE Dimethicone 5% Cream 100g. Dimethicone has water repellent properties which forms a protective 'glove' against water-soluble irritants.

Recommended as a barrier to protect skin against nappy rash, bedsores and urinary incontinence as well as irritating substances such as soap, detergents, water and some chemicals.

Ideally apply before putting babies nappies on, washing dishes, gardening, working on greasy/oily equipment or apparatus and prior to some medical procedures. May also be used to help restore dry skin chafing and dermatitis. Non-perfumed and suitable for use on sensitive skin.



Medent Code	Mfr Code	Description	Unit of Sale
DWCRBAR	PHZ847	healthE Dimethicone 5% Cream 100g	Each

Cadesorb Pro Modulating Tube

Cadesorb is a white, starch based sterile ointment Designed specifically to control the pH of its local environment.

By controlling pH levels and modulating protease activity, in addition to managing slough and exudate, Cadesorb help correct the natural balance in chronic wounds and stimulate healing.



Medent Code	Mfr Code	Description	Unit of Sale
DWOCA1	66000056	Cadesorb Pro Modulating 10g Tube	Box 5

Cavilon Barrier Spray No Sting 3M

A sting-free, alcohol-free liquid barrier film that dries quickly to form a breathable, transparent coating on the skin. Designed to protect intact or damaged skin from urine, feces, other body fluids, tape trauma, and friction.

A hypoallergenic, non-cytotoxic liquid film which will not sting even when applied to damaged or denuded skin.

Medent Code	Mfr Code	Description	Unit of Sale
DWSPC	3346	3M Cavilon Barrier Spray No Sting 28ml	Bottle



Cavilon Durable Barrier Cream 3M

A concentrated, therapeutic cream that provides durable, long-lasting protection from body fluids while it moisturizes the skin. Resists wash off and therefore requires less frequent applications.

The cream vanishes into the skin and does not transfer off to briefs or diapers. It contains dimethicone for skin protection.

Medent Code	Mfr Code	Description	Unit of Sale
DWCRB3	3355	3M Cavilon Durable Barrier Cream Tube 92g 3355	Each
DWCRB5	3354	3M Cavilon Durable Barrier Cream Tube 28g 3354	Each
DWCRB6	3353	3M Cavilon Durable Barrier Cream 2g Single Use Sachet 3353	Box 20



Cavilon Extra Dry Skin Cream 3M

A lipid-rich cream specially formulated to moisturize dry skin. Softens and replenishes rough, cracked or scaly skin. Ideal for individuals with diabetes or other dry skin conditions.



Medent Code	Mfr Code	Description	Unit of Sale
DWCRB4	3386	3M Cavilon Extra Dry Skin Cream 4oz Tube	Each

Cetaphil Moisturizing Lotion (Pump)

Cetaphil Moisturizing Lotion. Specifically formulated to soothe your chronically dry skin. Once applied, your skin feels soft and smooth. Fragrance free. Protects and enhances skin natural moisture

It is non-greasy, non-irritating and non-comedogenic (will not clog pores), leaving you feeling more comfortable in your own skin. Fast-acting, long-lasting, light-weight formula helps maintain healthy skin.



Medent Code	Mfr Code	Description	Unit of Sale
DWLOTCE	2144336	Cetaphil Moisturizing Lotion 250ml Pump	Each

Convatec Duoderm Paste

Convatec DuoDERM Paste is a hydrocolloid paste which can be used to manage chronic wounds - pressure ulcers (Stage I-IV), leg ulcers. Acute wounds - traumatic wounds (minor abrasions, lacerations) partial thickness burns and donor sites.

It incorporates a unique Convatec hydrocolloid formulation that distinguishes it from other hydrocolloid dressings.



Medent Code	Mfr Code	Description	Unit of Sale
DWPA1	187930	Convatec Duoderm Paste 30g	Each

Cream Anthisan

Relief from pain, itching and inflammation from insect bites, stings and nettle rash. Brand is Sanofi Aventis.



Medent Code	Mfr Code	Description	Unit of Sale
DWANTH	212245	Cream Anthisan 25g	Each

Cream Anti-flamme Tube

Natures Kiss Anti-flame Herbal Relief Crème. For bumps, bruises, aches and pains.

Soothes strained muscles, joint aches, bumps, bruises and common sports injuries. It contains natural ingredients including pure herb extracts of arnica, hypericum, calendula and peppermint oil.

Used and recommended by New Zealand physiotherapists. Not recommended for use during pregnancy.



Medent Code	Mfr Code	Description	Unit of Sale
DWCRANT	14886	Natures Kiss Anti-flame Herbal Relief Crème Tube 450g	Each
DWCRANT1	2128535	Natures Kiss Anti-flame Herbal Relief Crème Tube 100g	Each

Cream Aquacare

Aquacare HP Cream 100g - 10% Urea Cream. Aquacare is a topical cream formulated to moisturise, soften and smooth dry skin. For the treatment of eczema, itchiness and other hyperkeratotic skin conditions.

Apply 2-3 times daily to the affected area. For external use only.




Medent Code	Mfr Code	Description	Unit of Sale
DWCAQ	PHF185	Aquacare HP Cream 100g - 10% Urea Cream	Each

Cream Arnica Tube Weleda

Weleda Arnica Cream. Made from natural and organic ingredients. A traditional herbal medicine for the family First Aid Kit and Sports Bag.

Features and Benefits:

- ✓ Fast and effective relief for bruises, sprains, backache muscle aches and pains
- ✓ Contains pure essential oils of silver fir, lemon, orange, rosemary and sage
- ✓ Reduces pain and inflammation from injuries and over-exertion
- ✓ Supports the healing process
- ✓ Use before and after sport to relieve muscle ache
- ✓ Effective for the temporary relief of the pain of arthritis or rheumatism



Medent Code	Mfr Code	Description	Unit of Sale
DWCRARN	742988	Weleda Arnica Cream 15% 36ml Tube	Each

Cream Arnica Natures Kiss

Nature's Kiss Arnica Cream! Contains more Arnica for more care.

Soothes every day bumps, bruises, aches and pains of children and adults alike. It contains 15% Arnica oil and has no Peppermint, so it's safe to use for infants and even pregnant women.

Medent Code	Mfr Code	Description	Unit of Sale
DWCRARN1	2383616	Natures Kiss Cream Arnica 90g	Each



Cream Calendula Healing Tube Weleda

Weleda Cream Calendula Healing 36ml Tube A traditional herbal medicine for the family First Aid Kit.

For healing cuts, abrasions and wounds. Anti-inflammatory, supports the formation of healthy skin tissue. Promotes skin repair. Effective for slow healing wounds.

Medent Code	Mfr Code	Description	Unit of Sale
DWCRCAL	PHK567	Weleda Cream Calendula Healing 36ml Tube	Each



Cream Nutraplus 10% Tube

Nutraplus Cream 10% Urea - 100g. A moisturiser that contains urea as the active ingredient. Urea is a substance that is found naturally in the skin, which by its nature is very attracted to water. It is used to help rehydrate dry, scaly skin.

Fragrance free. Lanolin free. Skin treatment for dermatitis and dry scaly skin conditions like eczema.



Medent Code	Mfr Code	Description	Unit of Sale
DWCRN	N1038 255424	Nutraplus Cream 10% Urea - 100g	Each

Cream Secura (Extra Protection) Tube Smith & Nephew

Smith & Nephew Cream Secura (Extra Protection) tube. Helps to protect against skin irritation associated with nappy rash or exposure to faeces and urine. It provides durable, long-wearing protection that won't wash away, even in the most severe cases of incontinence.

Forms a durable barrier against faeces, urine, and wetness to help reduce the risk of skin irritation. pH balanced, and deodorising (Clove oil), also contains Aloe and Allantoin to moisturise and protect the skin.



Medent Code	Mfr Code	Description	Unit of Sale
DWCRTC1	59450425	Smith & Nephew Cream Secura (Extra Protection) 92g Tube	Each
DWCRS	59447525	Smith & Nephew Cream Secura (Extra Protection) 78g Tube	Each

Emulsifying Ointment

Emulsifying Ointment. A moisturising ointment that can also be used as a soap substitute. To use as a soap substitute, soften a small amount in a cup of warm water and use as you would a liquid soap. Does not lather but is an effective cleanser, without drying the skin.



Medent Code	Mfr Code	Description	Unit of Sale
DWOE	PHN135	Emulsifying Ointment 500g	Each

Flamazine Cream (Silver Sulfadiazine 1%) PM Smith & Nephew

The mechanism of silver sulfadiazine's antibacterial action has not been fully elucidated. After exposure to the drug, structural changes in the bacterial cell membrane occur, including distortion and enlargement of the cell and a weakening of the cell wall membrane. This is accompanied by reduced viability in sensitive strains due to interference with macromolecular synthesis.

The sulfadiazine moiety also provides a bacteriostatic action against sensitive organisms. In adults, up to 10% of the sulfadiazine may be absorbed and 60-85 % of the absorbed amount is excreted in the urine. In children with 13% body surface area burns, the urinary sulfadiazine concentration was 31.8 mg/L.



Medent Code	Mfr Code	Description	Unit of Sale
DWGFL5	66800532	Smith & Nephew Flamazine Cream 500g Jar (Silver Sulfadiazine 1%) PM	Each

Health E Fatty Cream

healthE Fatty Cream 500g. healthE Fatty Cream is recommended as a topical emollient for dry skin conditions such as eczema and dermatitis. This specially formulated dermatological cream makes it suitable as supportive therapy in conditions that have previously been treated with more potent corticosteroids. Apply 1-2 times daily to the affected area, or as directed by your health professional.

Features and Benefits:

- ✓ Highly moisturising, unique emollient formula
- ✓ Richer than aqueous cream and not as greasy as emulsifying ointment
- ✓ Fragrance, steroid and colour free
- ✓ Suitable for all skin types including sensitive skin
- ✓ Does not contain Sodium Lauryl Sulfate (SLS)



Medent Code	Mfr Code	Description	Unit of Sale
DWCF	2329468	Health E Fatty Cream 500g	Each

Pro Block Sunscreen, SPF

Pro Block Sunscreen. SPF 50+ with moisturising Aloe Vera & Vitamin E. Filters up to 98% of UV Rays. Broad Spectrum Protection from both UVA & UVB. 4 hours water resistant. Sensitive for the family with its Aloe Vera and Vitamin E content. Non-greasy makes it ideal for workers. Nanoparticle Free. Approved to AS/NZS 2604:2012. Sold Individually. Available in a 1L pump pack container, 500ml flip top bottle as well as a handy 125ml tube.



Medent Code	Mfr Code	Description	Unit of Sale
DWLOSUN1	SS1-50	Help-It Ultra Block Sunscreen 1L Pump	Each
DWLOSUN5	SS250-50	Help-It Ultra Block Sunscreen 250ml Flip Top Bottle	Each

Hirudoid Cream

Contain heparinoid, which is a similar to heparin, and is used to improve circulation after bruising or soft tissue injuries. Heparinoid works by dissolving small blood clots and improving the blood supply to the skin.

Hirudoid cream can be used to help relieve pain and inflammation associated with superficial thrombophlebitis.



Medent Code	Mfr Code	Description	Unit of Sale
DWHIRC	272380	Hirudoid Cream 40gm	Each

Liquid Paraffin Sterile

Sterile Liquid Paraffin BP Tube 5g. Liquid Paraffin is a transparent, colourless, almost odourless, oily liquid with good lubricating properties. Use as an aseptic dressing or lubricant for catheters and surgical instruments.

Medicinal liquid paraffin, also known as paraffinum liquidum, is a very highly refined mineral oil used for medical purposes.



Medent Code	Mfr Code	Description	Unit of Sale
DWCLP	204056	Liquid Paraffin Sterile 5g	Each

Lotion Alpha Keri

Alpha Keri Lotion The Dry and Problem Skin Solution. Alpha Keri Lotion is specifically formulated to protect against the elements, lock and seal in moisture and relieve itching from dry and scaly skin.

Alpha Keri Lotion is an oil-in-water soothing emollient lotion which relieves dryness and itching. Alpha Keri Lotion softens and smoothes dry skin. It increases the moisture content, lubricates and softens the skin, and locks in moisture to prevent further dryness.



Medent Code	Mfr Code	Description	Unit of Sale
DWLOTAK	2327457	Lotion Alpha Keri 300ml	Each

DP Lotion-HC 1%

DP Lotion-HC 1%. Hydrocortisone 1% w/w, lanolin 0.6% w/v and mineral oil 15.9% w/v lotion. A uniformly white, glossy, moderately viscous lotion.

Hydrocortisone is a corticosteroid that, when applied topically, often produces dramatic suppression of lesions having inflammation as a prominent feature. Such lesions include: eczema, various forms of dermatitis, psoriasis and intertrigo.

Lanolin is a substance obtained from sheep's wool that, when mixed with mineral oil, gives an emollient lotion that facilitates the absorption of hydrocortisone.



Medent Code	Mfr Code	Description	Unit of Sale
DWLOT	619787	DP Lotion-HC 1%	Each

Lyderm Cream

Lyderm Cream is a topical treatment for scabies, a contagious parasitic infection, caused by parasitic mite. It contains permethrin, a synthetic insecticide belonging to the family of chemicals known as pyrethroids.

Lyderm Cream kills the scabies parasite, preventing spread of infection and eventually relieving symptoms of the infection.



Medent Code	Mfr Code	Description	Unit of Sale
DWCRLY	479233	Lyderm Cream 30g	Each

Mebo Anti Itch Ointment

Mebo Anti-Itch Ointment. MEBO's unique formulations stimulate stem cell regeneration. The ointments promote the natural regeneration of stem cells located at the bottom of the hair follicle. These cells are capable of producing not only the hair shaft but also the daughter cells that create the dermis and epidermis.

MEBO anti itch works by assisting the natural repair of cracked or chafed skin and by aiding the natural repair of nerve endings in the skin tissue.



Medent Code	Mfr Code	Description	Unit of Sale
DWCMEBO1	20120000	Mebo Anti Itch Ointment 30g	Each

Mebo Burn Repair Ointment

MEBO Burn Repair Ointment. MEBO burn repair works by helping to cool the affected area, by helping the body's own immune defences, by moisturising and protecting the area, by encouraging stem cell regeneration for the growth of new skin, by helping to minimise scarring.



Medent Code	Mfr Code	Description	Unit of Sale
DWCMEBO2	20120002	Mebo Burn Repair Ointment 40g	Each

Mebo Wound Repair Ointment

MEBO Wound Repair Ointment. MEBO wound repair works by helping the body's own immune defences, by assisting the affected area to liquefy dead tissue by encouraging stem cell regeneration for the growth of new skin, by helping to minimise scarring.



Medent Code	Mfr Code	Description	Unit of Sale
DWCMEBO	20120001	Mebo Wound Repair Ointment 40g	Each

NU-GEL® Hydrogel with Alginate

NU-GEL® Hydrogel with Alginate. Transparent hydro-active amorphous gel containing sodium alginate that gently and effectively rehydrates necrotic tissue and fibrinous slough, and is designed to create an ideal moist wound healing environment for granulating and epithelising wounds.

Gently and effectively provides rehydration of the wound for autolytic debridement, significantly decreasing the amount of necrotic tissue and slough. The alginate increases the absorptive capacity of NU-GEL® Hydrogel and improves the ability to control wound exudate for longer wear time. Can be applied to the wound with a simple one-handed technique.



Medent Code	Mfr Code	Description	Unit of Sale
DWGNU	2242842	NU-GEL® Hydrogel with Alginate	Each

Oil Alpha Keri

Alpha Keri Oil. Specifically formulated as a total body bath/shower treatment that works to revive and relieve dry and sensitive skin associated with a variety of skin problems. Alpha Keri Oil moisturises dry skin by improving the hydrating properties of skin.

Soaking in an Alpha Keri Oil bath, or applying to wet skin after showering, helps seal and lock the moisture the stratum corneum takes up during the bath or shower. It is available in both fragrance free and fragrance fresh varieties.



Medent Code	Mfr Code	Description	Unit of Sale
DWOILAK	2327449	Oil Alpha Keri 300ml	Each
DWOILAK1	268003	Oil Alpha Keri 1L	Each
DWOILAK2	2149168	Oil Alpha Keri 500ml	Each

Ointment Bepanthen

Bepanthen® Cream. Antiseptic. Antiseptic cream with vitamin B5. Soothes, heals & protects. Used as a Baby Nappy Care Ointment, effective for the relief of nappy rash & mothers use in nipple care. Used to treat minor burns, cuts, abrasions, scalds, bites, stings minor skin irritations & sunburn.



Medent Code	Mfr Code	Description	Unit of Sale
DWBEP	2194929	Bepanthen® Cream Ointment 100g	Each

Ointment Ungvita

Ungvita Vitamin A Ointment (Roche). Contains Vitamin A in a petroleum jelly base to help repair damaged skin.

Provides effective treatment for dry skin conditions such as chafed skin, cracked heels and cracked nipples, and minor skin disorders such as wounds, scratches, eczema, dermatitis, cuts and abrasions

Each gram contains Retinyl palmitate (Vitamin A) 1800 IU



Medent Code	Mfr Code	Description	Unit of Sale
DWUVA	268321	Ungvita Vitamin A Ointment (Roche) 50g Tube	Each

Ointment Vitamin A

Home Essentials Vitamin A Protective Ointment.

Vitamin A is used to protect dry skin and used by beauticians to help smooth out wrinkles and help protect the skin from signs of ageing.

Protective ointment on skin may assist pimple, blemishes and fine wrinkles, help repair damaged skin, help treat dry skin conditions such as chafed skin and cracked heels, help treat minor skin conditions.



Medent Code	Mfr Code	Description	Unit of Sale
DWCVA	2350157	Home Essentials Vitamin A Protective Ointment 25g	Each
DWCVA1	2326809	Home Essentials Vitamin A Protective Ointment 50g	Each

Opsite Spray Dressing

OpSite Spray Transparent Film Dressing Quick drying film which is easy and convenient to apply. The film helps to provide protection for dry minor surgical and surface wounds.

Permeable to moisture vapour and air. Water resistant. Conformable and elastic. Quick and easy to apply. This product contains no Chlorofluorocarbon (CFC) propellants.



Medent Code	Mfr Code	Description	Unit of Sale
DWSPO	66004978	Opsite Spray Dressing 100ml	Each

Paste Magnesium Sulphate

Pharmacare Magnesium Sulphate Paste BP. A drawing paste to help treat minor skin infections including boils and reduce swelling. Contains Glycerol 61% w/w, Magnesium Sulphate 38.5% w/w, Phenol 0.5% w/w.

Mix well before use and apply paste to piece of gauze or lint. Place on skin and bandage. Change after 12-24 hours. For external use only. Prolonged or repeated use may damage surrounding skin.



Medent Code	Mfr Code	Description	Unit of Sale
DWPMS	2350009	Paste Magnesium Sulphate 80g	Each

Sunblock SPF 50 Pro Bloc

Sunscreen 250ml Bottle - SS250.

- ✓ SPF 30+ with moisturising Vitamin E
- ✓ 4 hours water resistant
- ✓ Flip top bottle

Medent Code	Mfr Code	Description	Unit of Sale
DWLOSUN5	SS250	Sunblock SPF 50 Pro Bloc 250ml	Each



Vaseline Jelly Pot

Vaseline Jelly 50G Pot. Vaseline is used to reduce friction and prevent abrasions in the sporting field. It forms a seal around burns, grazes or other wounds to form a waterproof barrier to help prevent contamination or further infection.

It can help soothe sore lips, minor burns, cuts, scrapes and sunburn, and chafing.

Medent Code	Mfr Code	Description	Unit of Sale
DWWAS	303267	Sunblock SPF 50 Pro Bloc 250ml	Each



Waxsol Ear Drop

Waxsol Ear Drop 10ml. For easy removal of ear wax. Hardened ear wax can lead to hearing problems and ear aches.

Waxsol helps soften the wax in your ears and makes ear wax removal, pain-free and effortless. Contains Docusate Sodium 0.5%. Do not use Waxsol if you have a perforated ear drum or a sore ear. Waxsol should only be applied to the ear canal.

Medent Code	Mfr Code	Description	Unit of Sale
DWWXSOL	619477	Waxsol Ear Drop 10ml	Each



Zinc & Castor Ointment

Zinc and Castor Oil Ointment. For nappy rash treatment and rash prevention. As an emollient to soothe and moisturise the skin. Also protective with mild astringent properties.

Commonly used for mild dermatological conditions such as eczema and dermatitis.

Apply 2 - 3 times a day or at each nappy change as needed. Contains Zinc Oxide, white Beeswax, Castor Oil, Cetostearyl Alcohol, Arachis Oil. For external use only. Avoid contact with eyes.



Medent Code	Mfr Code	Description	Unit of Sale
DWCRZC	544345	Zinc & Castor Ointment 500g Tub	Each

CATEGORY: DRESSING WOUND – TULLE GRAS

Actisorb Plus 25 Tulle Gras

Actisorb Plus 25. An activated charcoal dressing with silver, enclosed in a non-adherent nylon sleeve. Triple action that not only masks odour but also addresses its cause.

Activated charcoal traps odour in the dressing. It binds bacteria and toxins that impair the healing process. Silver kills the bacteria in the dressing.



Medent Code	Mfr Code	Description	Unit of Sale
DWTGACT	MAP105	Actisorb Plus 25 10.5 x10.5 (10's) Tulle Gras	Box 10

Activon Manuka Honey Tulle Dressings

Activon Tulle Dressing. A knitted viscose mesh primary dressing impregnated with 100% Medical Grade Manuka honey. Activon creates a moist healing environment and effectively eliminates wound odour while providing antibacterial action.



Medent Code	Mfr Code	Description	Unit of Sale
DWTGACT1	CR3658	Activon Manuka Honey Tulle Dressings 10cm x 10cm	Box 5

Adaptic Digit Finger Dressing

Adaptic® Digit Non-Adhering Dressing. Made of unique tacky silicone coated open-knit cellulose fabric with a double-coiled tubular bandage that is specifically designed for dressing fingers and toes.

Features and Benefits:

- ✓ Comfortable, convenient to apply and easy to remove for mild to severe digit injuries
- ✓ No applicator or scissors are needed making it easy and simple to use
- ✓ Conformable and non-bulky, allowing maximum digit flexibility
- ✓ Removes cleanly from the wound, minimising trauma and pain




Medent Code	Mfr Code	Description	Unit of Sale
DWTGADDL	MAD023	Adaptic Digit Finger Dressing Large 2.8cm	Box 10
DWTGADDM	MAD013	Adaptic Digit Finger Dressing Medium 2.4cm	Box 10

Adaptic Non Adherent Tulle Gras

Adaptic® Non-Adhering Dressing. Protects the fragile tissue in wounds by its unique structure with small mesh size, preventing tissue adherence to either the ADAPTIC® dressing or the secondary dressing on top of it.

Features and Benefits:

- ✓ Designed to help protect the wound
- ✓ Prevents dressing from adhering to the wound
- ✓ Primary dressing made of knitted cellulose acetate fabric allows ADAPTIC® to be cut to wound size without unravelling or shredding
- ✓ Impregnated with a specially formulated petrolatum emulsion
- ✓ Dressing removes cleanly from the wound, minimizing trauma and pain



Medent Code	Mfr Code	Description	Unit of Sale
DWTGAD1	2012	Adaptic Non Adher 7.6 x 7.6cm Tulle Gras	Box 50
DWTGAD2	2015	Adaptic Non Adher 7.6 x 24.6cm Tulle Gras	Box 24
DWTGAD3	2014	Adaptic Non Adher 7.6 x 40.6cm Tulle Gras	Box 36

Atrauman Dg Tulle Gras

Atrauman Tulle Gras Dressings. Low adherent, fine-weave tulle dressing impregnated with a non-medicated, non-petroleum based triglyceride for atraumatic wound treatment.

Features and Benefits:

- ✓ Soft thin support material of Atrauman is made of polyester
- ✓ Ensures close contact with the whole surface of the wound, allowing an easy passage for exudate and reducing the risk of maceration
- ✓ Fine weave and hydrophobic surface counteracts adhesion to the wound by helping to prevent new tissue from penetrating the dressing
- ✓ Enhanced by triglyceride, which keeps the wound edges soft and supple
- ✓ Can be combined with numerous secondary dressings



Medent Code	Mfr Code	Description	Unit of Sale
DWTGAT1	499510	Atrauman Dg 5 x 5 cm Tulle Gras	Box 10
DWTGAT2	499550	Atrauman Dg 5 x 5 cm Tulle Gras	Box 50
DWTGAT3	499513	Atrauman Dg 7.5 x 10cm Tulle Gras	Box 10
DWTGAT4	499553	Atrauman Dg 7.5 x 10 cm Tulle Gras	Box 50

Bactigras (Antiseptic) Tulle Gras Smith & Nephew

Smith & Nephew Bactigras* is an antiseptic, soft paraffin dressing which soothes and protects the wound whilst helping to reduce wound infection and inflammation. It has low adherence and allows the wound to drain freely into an absorbent secondary dressing. Suitable for covering wounds such as minor burns, lacerations, abrasions, graft sites and leg ulcers.

Features and Benefits:

- ✓ Soft thin support material of Atrauman is made of polyester
- ✓ Ensures close contact with the whole surface of the wound, allowing an easy passage of exudate and reducing the risk of maceration
- ✓ Fine weave and hydrophobic surface counteracts adhesion to the wound by helping to prevent new tissue from penetrating the dressing
- ✓ Enhanced by triglyceride, which keeps the wound edges soft and supple
- ✓ Can be combined with numerous secondary dressings



Medent Code	Mfr Code	Description	Unit of Sale
DWTGBAC	7456	Smith & Nephew Bactigras (Antiseptic) 5 x 5cm Tulle Gras	Box 50
DWTGBAC1	7457	Smith & Nephew Bactigras (Antiseptic) 10 x 10cm Tulle Gras	Box 10

Cuticerin Tulle Gras Smith & Nephew

Smith & Nephew Cuticerin is an impregnated dressing made of a smooth acetate fabric, with a defined mesh size, coated with a water repellent ointment containing petrolatum, paraffin and Eucerite. Cuticerin is a low adherent wound contact layer.

Medent Code	Mfr Code	Description	Unit of Sale
DWTGCUT	66045560	Smith & Nephew Cuticerin 7.5 x 7.5 cm Tulle Gras	Box 50
DWTGCUT1	66045561	Smith & Nephew Cuticerin 7.5 x 20 cm Tulle Gras	Box 50



Jelonet Tulle Gras Sterile Smith & Nephew

Smith & Nephew Jelonet Tulle Gras Sterile. Low Adherent, Impregnated Gauze Dressing. Non-Medicated Paraffin Gauze Dressing. Soft paraffin base. Sterile leno weave presentation.



Medent Code	Mfr Code	Description	Unit of Sale
DWTGJEL	7404	Smith & Nephew Jelonet 10cm x 10cm Tulle Gras Sterile	Box 10
DWTGJ1	7403	Smith & Nephew Jelonet 5cm x 5cm Tulle Gras Sterile	Box 50
DWTGJ2	7409	Smith & Nephew Jelonet 10cm x 10cm Tulle Gras Sterile	Box 100
DWTGJ4	7459	Smith & Nephew Jelonet 10cm x 40cm Tulle Gras Sterile	Box 10

Paranet Tulle Gras

Paranet™ Sterile Gauze Dressing BP. For use on superficial wounds, burns, skin grafts and traumatic injuries. Recommended as a primary wound contact layer. Paranet™ should be used in conjunction with an appropriate secondary dressing.

Made from 100% cotton open-mesh leno weave, impregnated with yellow soft paraffin. Presented in gauze strips and multi-piece dressings, Paranet™ is packed in polypropylene containers and sealed with a tamper-evident label.



Medent Code	Mfr Code	Description	Unit of Sale
DWTGP3	3475	Paranet 5 cm x 5 cm Sterile Tulle Gras	Box 36
DWTGP4	3481	Paranet 10 cm x 7 m Strip Tulle Gras	Each
DWTGP5	3479	Paranet 10 cm x 10cm Tulle Gras	Each
DWTGP6	3485	Paranet 10 cm x 10cm Sterile Tulle Gras	Box 10

CATEGORY: DRESSING WOUND – WIPES

Bzk Antiseptic Towelette Alcohol Free Wipes

Help-IT BZK Antiseptic Towelette Alcohol Free Wipes.

These antiseptic towelettes, with the active ingredient Benzalkonium Chloride, are ideal for cleansing wounds when alcohol is not advisable.

Single-use design eliminates the worry of cross contamination. Full-size, individually wrapped. Alcohol-free, sterile towelettes. For the cleaning of cuts, abrasions, and other minor wounds. Towelettes measure 19.5 x 10.5cm.

Medent Code	Mfr Code	Description	Unit of Sale
DWWIP4	WIPEAMT50S	BZK Antiseptic Towelette Alcohol Free Wipes	Box 50



Cavilon No Sting Barrier Film Wipes 3M

3M Cavilon No Sting Barrier Film Wipes. A sting-free, alcohol-free liquid barrier film that dries quickly to form a breathable, transparent coating on the skin. Designed to protect intact or damaged skin from urine, feces, other body fluids, tape trauma, and friction. A hypoallergenic, non-cytotoxic liquid film which will not sting even when applied to damaged or denuded skin.

Medent Code	Mfr Code	Description	Unit of Sale
DWCRB1	3344	3M Cavilon No Sting Barrier Film 1ml Wipes	Box 25



Convatec Remove Adhesive Wipes Convacare

Convatec Remove Adhesive Wipes Convacare. Adhesive Remover Wipes gently remove all adhesives including skin barriers, tapes and hydrocolloid dressings

Medent Code	Mfr Code	Description	Unit of Sale
DWWIP1	37443	Convatec Remove Adhesive Wipes Convacare	Box 100



Wipes Remove Smith & Nephew

Smith & Nephew Wipes Remove Solvent Wipes. Remove Solvent formulated to dissolve adhesives and assist in removing acrylic-based, rubber-based, and hydrocolloid-based residues from the skin.

Medent Code	Mfr Code	Description	Unit of Sale
DWWIP2	403100	Wipes Remove Smith & Nephew	Box 50




CATEGORY: DRESSING WOUND – OTHERS

Kendall Argyle Penrose Tubing

ARGYLE™ Penrose Tubing. Latex in construction, radiopaque and sterile. The staples are available in two sizes, wide (clip span 13.0 mm before closure) and regular (clip span 9.9 mm before closure).

Features and Benefits:

- ✓ Latex construction
- ✓ Durable
- ✓ Radiopaque
- ✓ Allows placement verification by x-ray
- ✓ Sterile
- ✓ Ease of use



Medent Code	Mfr Code	Description	Unit of Sale
EMPENTU	8888-515809	Kendall Argyle Penrose Tubing	Box 50

Puritan Wound Depth Indicator

Puritan Wound Depth Indicator. Disposable depth indicator, sterile and convenient allowing for easy measurement of wound depth.

Medent Code	Mfr Code	Description	Unit of Sale
EMWAVD	25-15061PF	Puritan Wound Depth Indicator	Box 50



Seroma Wound Drainage System & Suc Res

Seroma Cath 2 inch with Catheter - Wound drainage system 14G 2 inch Teflon catheter/stylet needle/plug with two side ports, a 2.75 inch blunt metal catheter with two side ports and a connecting drainage tube for a bulb suction reservoir.


Medent Code	Mfr Code	Description	Unit of Sale
EMCWD	1002/2001	Seroma Wound Drainage System & Suc Res 2 Inch	Each



Wound Mapping Grid 3M

3M Wound Mapping Grid. Clear, adhesive wound mapping grid for monitoring and recording the progress of wounds. Hygienic - backing material is thrown away. Simple to use.

Medent Code	Mfr Code	Description	Unit of Sale
DWMG	GRID1	3M Wound Mapping Grid 10 x 10	Box 50



Wound 3M Measuring Guide Pad

3M Wound Measuring Guide Pad. Discard after each use.



3M™ Wound Measuring Guide

Discard after single use

Date _____ / _____ / _____ Patient Initials/Number _____



Medent Code	Mfr Code	Description	Unit of Sale
DWWMG	XH0007 00498	3M Wound 3M Measuring Guide Pad	25 Sheets

To Order Contact Us at:

Medent Medical

Unit 4, 110 Whakatiki St
PO Box 40356
Upper Hutt
New Zealand

Free Phone: 0800 263 3368

Phone: 64 4 527 7702

Free Fax: 0800 527 7499

Fax: 64 4 5277499

Email Enquires: enquiry@medent.co.nz

Finding a product and ordering it – you choose an option that suits you.
Either talk with one of our team or order right here online!

- Order from our sales representative when they visit you.
To arrange a visit please phone us on 0800 263 3368
- You can phone us on 0800 263 3368
or Fax us on 0800 527 7499

All products and more available at www.medent.co.nz



Dressing Wound



PRODUCT CATALOGUE

2015 | Edition 1



[Online Sales & Ordering](#)



HEAD OFFICE:
PH 04 527 7702 FAX 04 527 7499
PH OUTSIDE WELLINGTON: 0800 26 333 68
TOLL FREE FAX: 0800 52 774 99

[Email Enquiries](#)