



**HEINE GAMMA® Series**

Shockproof and Latex-free Sphygmomanometers  
and Stethoscopes



## HEINE GAMMA® STETHOSCOPES

- ✦ **Solid chest piece** and ideal weight avoid stray noise from manipulation and ensure excellent resonance. Hygienic matt chrome finish.
- ✦ **Chest piece with non-chill rims** for the comfort of the patients.
- ✦ Ergonomic, chrome-plated binaurals. Individually adjustable.
- ✦ **High-quality Y-piece tubing** for perfect acoustics and sound reproduction. Length 56 cm.
- ✦ **Latex-free** materials.

## [ HEINE GAMMA® 3.1 Pulse Stethoscope ]

GAMMA 3.1 | ORDER NO.:

M-000.09.941

**For fine acoustic performance.**

- ✦ **Stethoscope with solid flat chest piece.** The chest piece can be easily inserted under the sphygmomanometer cuff.
- ✦ **High-quality membrane,** 44 mm dia.

- ✦ **Single-tube construction** with excellent sound transmission.
- ✦ **Includes ear olives** (large/soft) for sound insulation.
- ✦ **Optional available:** Ear olives, (small/hard) (M-000.09.946).



## [ HEINE GAMMA® 3.2 Acoustic Stethoscope ]

GAMMA 3.2 | ORDER NO.:

M-000.09.942

**For excellent acoustic performance.**

Technical data as GAMMA 3.1, but with a double chest piece.

- ✦ **Combined double chest piece** with two different frequency ranges to choose from: A high-quality membrane on one side and a conical acoustic head on the other side.

- ✦ **High-quality membrane,** 44 mm diameter for excellent acoustic performance and sound transmission.
- ✦ **Cone diameter 30 mm.**



## [ HEINE GAMMA® 3.3 Acoustic Stethoscope ]

GAMMA 3.3 | ORDER NO.:

M-000.09.943

**For the examination of children.**

Technical data as GAMMA 3.2.

- ✦ **Reduced contact area** of 31 mm makes GAMMA 3.3 ideal for examining children.



## [ HEINE GAMMA® C3 Cardio Stethoscope ]

GAMMA C3 | ORDER NO.:

M-000.09.944

**Specially-developed for cardiology.**

- ✦ **Exceptional acoustic performance,** particularly in the cardiac frequency range.
- ✦ **Combined double chest-piece** with two different frequency ranges. Instant selection of a double membrane on one side or a conical acoustic head on the other.

- ✦ **High-performance double membrane,** 45 mm diameter gives ideal sound transmission and perfect skin contact for exceptional sound reproduction.
- ✦ **Two-in-one tube system** for enhanced performance. Very supple, yet kink-resistant with exceptional acoustic reproduction.
- ✦ **Includes ear olives** (large/soft) and one pair of replacements (small/hard) for sound insulation.

