

## Easy Access To Product Codes

This catalog contains a condensed listing of DAVIS & GECK and USS Suture product codes. For a complete list of all available suture codes please contact our customer service representative at: 1-800-962-9888, Option #4.

Your catalog has been organized to provide easy access to specific suture product codes and descriptive information.

	PRODUCT/BRAND	MATERIAL	PAGE
	<b>DEXON II, POLYSORB</b>	BRAIDED SYNTHETIC ABSORBABLE	21-22
	<b>BIOSYN</b>	MONOFILAMENT SYNTHETIC ABSORBABLE	23
	<b>MAXON</b>	MONFILAMENT SYNTHETIC ABSORBABLE	24
	<b>Chromic Gut</b>	NATURAL ABSORBABLE, CHROMIC	25
	<b>Plain Gut</b>	NATURAL ABSORBABLE, PLAIN	25
	<b>MONOSOF, DERMALON</b>	MONOFILAMENT NYLON	26
	<b>NOVAFIL</b>	SYNTHETIC NONABSORBABLE MONOFILAMENT POLYBUTESTER	27
	<b>Surgical Silk, SOFSILK</b>	NONABSORBABLE COATED BRAIDED SILK	28

# Absorbable Sutures

## DEXON II


	D&G		Box Qty	Size								
	Order #	Color		Length	6-0	5-0	4-0	3-0	2-0	0	1	2
no needle	8886900351	Green	60"	3dz						x		
pre-cut strands	8886900741	Green	3x18"	3dz				x				
	8886900751	Green	3x18"	3dz					x			
ligating reel	8886900852	Beige	60"	3dz					x			

### REVERSE CUTTING D&G USS ETHICON

D&G	USS	ETHICON
CE-2		
CE-2		
CE-4	C-13	FS-2
CE-4	C-13	FS-2
CE-4	C-13	FS-2
CE-4	C-13	FS-2
CE-4	C-13	FS-2
CE-6	C-14	FS-1
CE-6	C-14	FS-1
C-6	GS-10	CP-2
C-6	GS-10	CP-2
PRE-1	P-10	P-1
PRE-2	P-13	P-3
PRE-4	P-12	PS-2
PRE-4	P-12	PS-2

### D&G Box Size

Order #	Color	Length	Box Qty	Size							
				6-0	5-0	4-0	3-0	2-0	0	1	2
8886953821	Beige	18"	3dz		x						
8886953831	Beige	18"	3dz			x					
8886960821	Beige	18"	3dz		x						
8886960831	Beige	18"	3dz			x					
8886961221	Beige	27"	3dz		x						
8886961231	Beige	27"	3dz			x					
8886961241	Beige	27"	3dz				x				
8886961831	Beige	30"	3dz			x					
8886964147	Violet	30"	3dz				x				
8886974547	Violet	27"	3dz				x				
8886974541	Green	27"	3dz				x				
8886953211	Beige	18"	3dz	x							
8886953621	Beige	18"	3dz		x						
8886972821	Beige	18"	3dz		x						
8886972831	Beige	18"	3dz			x					




### TAPER POINT D&G USS ETHICON

D&G	USS	ETHICON
T-31	EV-23	SH-2
T-31	EV-23	SH-2
T-16	EV-25	SH-1
T-16	EV-25	SH-1
T-16	EV-25	SH-1
T-3	GS-23	FN-2,CT-3
T-5	EV-20	SH
T-5	EV-20	SH
T-5	EV-20	SH
T-5	EV-20	SH
T-19	GS-22	FN,CT-2
T-19	GS-22	FN,CT-2


### D&G Box Size


Order #	Color	Length	Box Qty	Size							
				5-0	4-0	3-0	2-0	0	1	2	
8886910727	Violet	27"	3dz	x							
8886910731	Green	27"	3dz		x						
8886911737	Violet	27"	3dz		x						
8886911747	Violet	27"	3dz			x					
8886911621	Beige	27"	3dz	x							
8886922641	Beige	27"	3dz			x					
8886923347	Violet	27"	3dz			x					
8886923337	Violet	27"	3dz		x						
8886923347	Violet	27"	3dz			x					
8886923357	Violet	27"	3dz				x				
8886924541	Green	27"	3dz			x					
8886924557	Violet	27"	3dz				x				




# Absorbable Sutures

## POLYSORB

	REVERSE CUTTING			USS		Box		Size							
	D&G	USS	ETHICON	Order #	Color	Length	Qty	6-0	5-0	4-0	3-0	2-0	0	1	2
	CE-3	C-12	FS-3	SL-630	Violet	30"	3dz		x						
	CE-3	C-12	FS-3	SL-607	Violet	30"	3dz			x					
	CE-3	C-12	FS-3	SL-608	Violet	30"	3dz				x				
	CE-4	C-13	FS-2	SL-635	Violet	30"	3dz			x					
	CE-4	C-13	FS-2	SL-636	Violet	30"	3dz				x				
	CE-4	C-13	FS-2	SL-822	Undyed	30"	3dz					x			
	CE-6	C-14	FS-1	SL-694	Violet	30"	3dz						x		
	CE-6	C-14	FS-1	SL-692	Undyed	30"	3dz			x					
	C-30	C-23	CP	SL-632	Violet	30"	3dz						x		
	PRE-1	P-10	P-1	SL-1625	Undyed	18"	3dz	x							
	PRE-2	P-13	P-3	SL-3687	Undyed	18"	1dz		x						
	PRE-4	P-12	PS-2	SL-1628	Undyed	18"	3dz						x		

	TAPER POINT			USS		Box		Size						
	D&G	USS	ETHICON	Order #	Color	Length	Qty	5-0	4-0	3-0	2-0	0	1	2
	TE-31	CV-23	RB-1	UL-202	Violet	30"	3dz	x						
	TE-31	CV-23	RB-1	UL-203	Violet	30"	3dz		x					
	TE-31	CV-23	RB-1	UL-204	Violet	30"	3dz			x				
	T-16	CV-25	SH-1	GL-181	Violet	30"	3dz		x					
	T-16	CV-25	SH-1	GL-182	Violet	30"	3dz				x			
	T-5	V-20	SH	GL-120	Violet	30"	3dz	x						
	T-5	V-20	SH	GL-121	Violet	30"	3dz		x					
	T-5	V-20	SH	GL-122	Violet	30"	3dz				x			
	T-8	V-30		GL-223	Violet	30"	3dz						x	

## POLYSORB - Ophthalmic

	POLYSORB		USS		Box		Size	
	Premium Point Spatula	USS	Order #	Color	Length	Qty	8-0	7-0
	SE-140-8	dbl-arm	L-2747	Violet	8"	12	x	
	SE-140-8	dbl-arm	L-2748	Violet	12"	12	x	
	SE-140-8	dbl-arm	L-1796	Violet	12"	12		x
	SE-140-8	single	L-2750	Violet	12"	12	x	
	SE-140-8	single	L-2751	Violet	12"	12		x

# Absorbable Sutures

## BIOSYN

### REVERSE CUTTING


D&G    USS    ETHICON

PRE-2	<b>P-13</b>	P-3
PRE-2	<b>P-13</b>	P-3
CE-4	<b>C-13</b>	FS-2
C-30	<b>C-23</b>	X-1
C-30	<b>C-23</b>	X-1
CE-4	<b>C-13</b>	FS-2
CE-6	<b>C-14</b>	FS-1
CE-6	<b>C-14</b>	FS-1
C-9	<b>GS-11</b>	CP-1
C-9	<b>GS-11</b>	CP-1

### USS

Order #    Color    Length    Box Qty    4-0    3-0    2-0    0    1    Size

<b>SM-3687</b>	Clear	18"	12dz	x					
<b>SM-3690</b>	Clear	18"	12dz		x				
<b>SM-691</b>	Clear	27"	3dz	x					
<b>SM-632</b>	Violet	30"	3dz		x				
<b>SM-633</b>	Violet	30"	3dz	x					
<b>SM-822</b>	Clear	27"	3dz		x				
<b>SM-693</b>	Clear	27"	3dz		x				
<b>SM-823</b>	Clear	27"	3dz			x			
<b>CM-816</b>	Violet	27"	3dz			x			
<b>CM-817</b>	Violet	27"	3dz				x		



### TAPER POINT


D&G    USS    ETHICON

T-31	<b>CV-23</b>	SH-2
T-31	<b>CV-23</b>	SH-2
T-31	<b>CV-23</b>	SH-2
T-5	<b>V-20</b>	SH-2
T-5	<b>V-20</b>	SH-2

### USS


Order #    Color    Length    Box Qty    5-0    4-0    3-0    2-0    0    1    Size

<b>UM-202</b>	Violet	30"	3dz	x					
<b>UM-203</b>	Violet	27"	3dz		x				
<b>UM-204</b>	Violet	27"	3dz			x			
<b>GM-121</b>	Violet	27"	3dz		x				
<b>GM-122</b>	Violet	27"	3dz			x			




# Absorbable Sutures

## MAXON



REVERSE CUTTING			D&G		Box		Size							
D&G	USS	ETHICON	Order #	Color	Length	Qty	7-0	6-0	5-0	4-0	3-0	2-0	0	1
CE-4	C-12	FS-3	8886661221	Green	30"	3dz			x					
CE-4	C-12	FS-3	8886663931	Green	30"	3dz				x				
CE-4	C-12	FS-3	8886663941	Green	30"	3dz					x			
CE-6	C-14	FS-1	8886664141	Green	30"	3dz					x			
CE-6	C-14	FS-1	8886664151	Green	30"	3dz						x		
CE-10	C-16	FSL	8886662641	Clear	27"	3dz						x		
C-6	GS-10	CP-2	8886674441	Clear	27"	3dz							x	
C-6	GS-10	CP-2	8886674551	Green	27"	3dz								x
PRE-1	P-10	P-1	8886653201	Clear	18"	3dz	x							
PRE-1	P-10	P-1	8886653311	Green	18"	3dz		x						
PRE-2	P-13	P-3	8886653621	Clear	18"	3dz			x					
PRE-2	P-13	P-3	8886653631	Clear	18"	3dz				x				
PRE-2	P-13	P-3	8886653921	Green	18"	3dz			x					
PRE-4	P-12	PS-2	8886660331	Green	18"	3dz					x			
PRE-4	P-12	PS-2	8886660341	Green	18"	3dz						x		
PRE-4	P-12	PS-2	8886660441	Clear	18"	3dz							x	




TAPER POINT			D&G		Box		Size						
D&G	USS	ETHICON	Order #	Color	Length	Qty	6-0	5-0	4-0	3-0	2-0	0	1
TE-9	CV-CVF	FS-3	8886611511	Green	18"	3dz	x						
TE-30	CV-22	RB-2	8886620511	Green	30"	3dz	x						
TE-30	CV-22	RB-2	8886620521	Green	30"	3dz		x					
T-31	EV-23	RB-1	8886620721	Green	30"	3dz			x				
T-31	EV-23	RB-1	8886620731	Green	30"	3dz				x			
T-31	EV-23	RB-1	8886620741	Green	27"	3dz					x		
T-16	CV-25	SH-1	8886621731	Green	27"	3dz					x		
T-16	CV-25	SH-1	8886621741	Green	27"	3dz						x	
T-5	EV-20	SH-1	8886623331	Green	27"	3dz						x	
T-5	EV-20	SH-1	8886623341	Green	27"	3dz							x
T-5	EV-20	SH-1	8886623241	Clear	27"	3dz							x
T-5*	V-20	SH-1	8886623541	Green	36"	3dz							x
T-19	GS-22	FN,CT-2	8886624541	Green	27"	3dz							x

\*double arm

# Absorbable Sutures

## Chromic Gut




### REVERSE CUTTING

D&G	USS	ETHICON
<b>CE-2</b>		C-3
CE-4	<b>C-13</b>	FS-2
CE-4	<b>C-13</b>	FS-2
CE-10	<b>C-16</b>	FSL
<b>C-6</b>	GS-10	CP-2
<b>C-6</b>	GS-10	CP-2

### D&G/USSC

Order #	Length	Box Qty	Size						
			5-0	4-0	3-0	2-0	0	1	
8886056421	18"	3dz	x						
SG635	18"	3dz		x					
SG636	18"	3dz			x				
SG637	36"	3dz				x			
8886056041	27"	3dz			x				
8886056051	27"	3dz				x			



### TAPER POINT

D&G	USS	ETHICON
T-16	<b>CV-25</b>	SH-1
T-16	<b>CV-25</b>	SH-1
TE-5	<b>V-20</b>	SH
TE-5	<b>V-20</b>	SH
TE-5	<b>V-20</b>	SH
TE-5	<b>V-20</b>	SH
T-19	<b>GS-22</b>	CT-2
T-19	<b>GS-22</b>	CT-2
T-19	<b>GS-22</b>	CT-2

### USS

Order #	Length	Box Qty	Size				
			4-0	3-0	2-0	0	1
GG-181	30"	3dz	x				
GG-182	30"	3dz		x			
GG-121	30"	3dz	x				
GG-122	30"	3dz		x			
GG-123	30"	3dz			x		
GG-124	30"	3dz				x	
CG-882	30"	3dz		x			
CG-883	30"	3dz			x		
CG-884	30"	3dz				x	


### No Needle

Standard length
Ligating reel
Standard length
Standard length
Standard length
Standard length
Standard length
Ligating reel
Ligating reel

### D&G

Order #	Length	Box Qty	Size				
			3-0	2-0	0	1	2
G-83	18"	2dz			x		
LG-113	60"	2dz			x		
G-214	60"	2dz				x	
8886050641	54"	3dz	x				
8886050651	54"	3dz		x			
8886050661	54"	3dz			x		
8886050681	54"	3dz					x
8886051852	54"	1dz		x			
8886051862	54"	1dz			x		

## Plain Gut



### D&G

D&G	USS	ETHICON
<b>no needle</b>		
<b>TE-10</b>	CV-1,CVF	BV-1
<b>CE-2</b>		C-3
CE-4	<b>C-13</b>	FS-2
CE-4	<b>C-13</b>	FS-2
CE-6	<b>C-14</b>	FS-1

### D&G

Order #	Length	Box Qty	Size						
			6-0	5-0	4-0	3-0	2-0	0	1
8886000571	54"	3dz							x
8886011013	18"	1dz	x						
8886032316	18"	1dz	x						
SG821	30"	3dz			x				
SG822	30"	3dz				x			
SG810	30"	3dz				x			

# Non-Absorbable Sutures

## MONOSOF/DERMALON


	REVERSE CUTTING			USS		Box Qty	Size							
	D&G	USS	ETHICON	Order #	Length		6-0	5-0	4-0	3-0	2-0	0	1	2
	CE-3	C-12	FS-3	SN-660	18"	3dz	x							
	CE-3	C-12	FS-3	SN-659	18"	3dz		x						
	CE-4	C-13	FS-2	SN-661	18"	3dz		x						
	CE-4	C-13	FS-2	SN-662	18"	3dz			x					
	CE-4	C-13	FS-2	SN-641	30"	3dz				x				
	CE-4	C-13	FS-2	SN-644	36"	3dz			x					
	CE-4	C-13	FS-2	SN-665	36"	3dz				x				
	CE-4	C-14	FS-2	SN-658	18"	3dz			x					
	CE-4	C-14	FS-2	SN-663	18"	3dz				x				
	CE-4	C-14	FS-2	SN-669	30"	3dz				x				
	CE-4	C-14	FS-2	SN-649	36"	3dz			x					
	CE-4	C-14	FS-2	SN-648	36"	3dz				x				
	CE-8	C-15	FS	SN-764	18"	3dz					x			
	CE-8	C-15	FS	SN-664	30"	3dz					x			
	CE-10	C-16	FSLX	SN-670	30"	3dz			x					
CE-10	C-16	FSLX	SN-671	30"	3dz				x					
CE-10	C-16	FSLX	SN-688	30"	3dz					x				
CE-10	C-16	FSLX	SN-788	36"	3dz					x				
PRE-1	P-10	P-1	SN-1697	18"	3dz	x								
PRE-2	P-13	P-3	SN-1698	18"	3dz		x							
PRE-2	P-13	P-3	SN-1699	18"	3dz			x						
PRE-4	P-12	PS-2	SN-3667	18"	1dz			x						
PRE-6	P-14	PS-1	SN-1663	18"	3dz				x					
C-4	C-22	J-1	SN-630	30"	3dz			x						
C-4	C-22	J-1	SN-631	30"	3dz				x					

### STRAIGHT REVERSE CUTTING

	D&G	USS	ETHICON	USS		Box Qty	Size 3-0
				Order #	Length		
	CS-1	SC-2	KS	SN-627	30"	3dz	x
		SC-1	TS	SN-616	36"	3dz	x

# Non-Absorbable Sutures

## NOVAFIL




### REVERSE CUTTING

D&G	USS	ETHICON
CE-4	C-13	FS-2
CE-4	C-13	FS-2
CE-6	C-14	FS-1
CE-6	C-14	FS-1
CS-1	SC-2	KS
PRE-2	P-13	P-3
PRE-2	P-13	P-3
PRE-2	P-13	P-3
PRE-3	P-11	PS-3
PRE-4	P-12	PS-2
PRE-4	P-12	PS-2

### D&G

Order #	Color	Length	Box Qty	Size								
				6-0	5-0	4-0	3-0	2-0	0	1	2	
8886442023	Blue	18"	1dz		x							
8886442033	Blue	18"	1dz			x						
8886442243	Blue	18"	1dz				x					
8886442253	Blue	18"	1dz					x				
8886449041	Blue	30"	3dz					x				
8886441213	Blue	18"	1dz	x								
8886441223	Blue	18"	1dz		x							
8886441233	Blue	18"	1dz			x						
8886441423	Blue	18"	1dz		x							
8886441633	Blue	18"	1dz			x						
8886441643	Blue	18"	1dz				x					



### TAPER POINT

D&G	USS	ETHICON
T-31	EV-23	RB-1
T-5	EV-20	SH
T-5	EV-20	SH
T-19	GS-22	FN,CT-2

### D&G

Order #	Color	Length	Box Qty	Size								
				5-0	4-0	3-0	2-0	0	1	2		
8886452021	Blue	36"	3dz	x								
8886445241	Blue	30"	3dz			x						
8886445251	Blue	30"	3dz				x					
8886445951	Blue	30"	3dz				x					




# Non-Absorbable Sutures


## Silk

D&G	D&G Order #	Color	Length	Box Qty	Size
pre-cut strands	8886103541	Black	7x30"	3dz	x
ligating reel	8886102342	Black	144"	2dz	x

REVERSE CUTTING			D&G	Color	Length	Box Qty	Size				
D&G	USS	ETHICON	Order #				6-0	5-0	4-0	3-0	2-0
CE-21	C-1	C-2	8886129613	Black	18"	1dz	x				
CE-2			8886107311	Black	18"	3dz	x				
CE-2			8886107321	Black	18"	3dz		x			
CE-2			8886107331	Black	18"	3dz			x		
CE-2 double armed			8886123813	Black	18"	1dz	x				
CE-2 double armed			8886123833	Black	18"	1dz			x		
CE-4	C-13	FS-2	8886107421	Black	18"	3dz		x			
CE-4	C-13	FS-2	8886107431	Black	18"	3dz			x		
CE-4	C-13	FS-2	8886107441	Black	18"	3dz				x	
CE-4	C-13	FS-2	8886107433	Black	18"	1dz			x		
CE-6	C-14	FS-1	8886107641	Black	18"	3dz				x	
CE-6	C-14	FS-1	8886107651	Black	18"	3dz					x



TAPER POINT			D&G	Color	Length	Box Qty	Size						
D&G	USS	ETHICON	Order #				6-0	5-0	4-0	3-0	2-0	0	1
T-16	CV-25	SH-1	8886112421	Black	30"	3dz		x					
T-16	CV-25	SH-1	8886112451	Black	30"	3dz					x		
T-5	EV-20	SH	8886112341	Black	30"	3dz				x			
T-5	EV-20	SH	8886112351	Black	30"	3dz					x		
T-12	GS-21	CT-1	8886113261	Black	30"	3dz						x	




## SOFSILK

REVERSE CUTTING			D&G	Color	Length	Box Qty	Size						
D&G	USS	ETHICON	Order #				5-0	4-0	3-0	2-0	0	1	
PRE-6	P-14	PS-1	SS-1684	Black	18"	3dz			x				
CE-23*	HE-3	G-3	S-1783	Black	18"	1dz		x					

TAPER POINT			D&G	Color	Length	Box Qty	Size						
D&G	USS	ETHICON	Order #				5-0	4-0	3-0	2-0	0	1	
TE-5	V-20	SH	GS-63M	Black	5x18"	2dz			x				
TE-5	V-20	SH	GS-62M	Black	5x18"	2dz				x			
TE-1	CV-11	C-1	VS-890	Black	30"	3dz	x						
T-37	CV-20	RB-3	GS-44M	Black	18"	2dz		x					
precut strands			S-183	Black	18"	3dz				x			
precut strands			S-185	Black	18"	3dz					x		
precut strands			S-307	Black	6x30"	3dz						x	

\*double arm






# Suture Needles

## General Surgery




### Taper Point Needles

**D&G** **USS** Ethicon






1/2 Circle

T-3	GS-23	FN-2, CT-3	_____	
T-19	GS-22	FN, CT-2	_____	
T-12	GS-21	CT-1	_____	

5/8 Circle

TT-2	GU-46	UR-6	_____	
TT-3	GU-45	UR-5	_____	
TT-4	GU-44	UR-4	_____	

1/2 Circle

T-3	HGS-23	MO-7	_____	
T-40	HGS-22	MO-6	_____	
—	HGS-20	MO-5	_____	
T-56	HGS-21	MO-4	_____	
T-25	HGS-24	MO-2	_____	














### Straight Taper Point Needles

TS-9	ST-1	ST	_____	51 mm
------	------	----	-------	-------

### Taper Point Needles Cardiovascular, Gastro-Intestinal, Pediatric






**D&G** **USS** Ethicon

3/8 Circle

TE-10, TE-11, CV-310	CV-1, CVF-1, CV, CVF	BV-1, BV	_____	
TE-1, CV-301	CV-11, CVF-11,	C-1	_____	
TE-7, CV-307	CV-15	BB	_____	
T-37, CV-337	CV-20	RB-3	_____	
—	CVF-21	TF	_____	
TE-30, CV-330	CV-22, CVF-22	RB-2	_____	
TE-31, CV-331	CV-23, CVF-23	RB-1	_____	
T-31	EV-23	RB-1	_____	
CV-316	CV-24	SH-2	_____	
T-16	CV-25	SH-1	_____	
T-16	EV-25	SH-1	_____	
TE-5	V-20, VF-20	SH	_____	
T-5	EV-20	SH	_____	

### Tapercutting Needles Vascular & Pediatric

3/8 Circle

DTE-10	KV-1, KVF-1	CC	_____	
DTE-1, DTE-2	KV-11, KVF-11	CC-1	_____	
DT-31	KV-15	V-4	_____	
DT-31	KV-5	MO-7	_____	
DT-5	KV-7	MO-6	_____	

# Suture Needles

## Plastic/Cosmetic

### Premium Reverse Cutting Needles

3/8 Circle

D&G	USS	Ethicon	
PRE-20	P-16	P-6	7 mm
PRE-1, SBE-1	P-10	P-1	11 mm
PRE-2, SBE-2	P-13	P-3	13 mm
PRE-3, SBE-3	P-11	PS-3	16 mm
PRE-4, SBE-4	P-12	PS-2	19 mm
PRE-6, SBE-6	P-14	PS-1	24 mm
			13 mm
PR-2	P-22	PS-6	16 mm
PR-4	P-242	PS-4	

### Premium Conventional Cutting Needles

3/8 Circle

SCE-2	PC-13	PC-1	11 mm
—	PC-11	PC-3	16 mm
SCE-4	PC-12	PC-5	19 mm

## Cuticular

### Reverse Cutting Needles

D&G USS Ethicon

3/8 Circle

CE-21	C-1	C-2	11 mm
CE-3	C-12	FS-3	16 mm
CE-4	C-13	FS-2	19 mm
CE-6	C-14	FS-1	24 mm
CE-8	C-15	FS	26 mm
C-2	C-21	M-1	13 mm
C-4	C-22	J-1	20 mm
C-30	C-23	X-1	24 mm
C-41	C-25	—	30 mm

### Reverse Cutting Needles

—	SC-11	—	11 mm
—	SC-4	SC-1	24 mm

## Microsurgery

### Taper Point Needles

3/8 Circle

TE-50	MV-50-3	BV-50-3/2	4 mm
TE-70	MV-70-3	BV-75-3	4 mm
TE-100	MV-100-3	BV-130-3	4 mm
TE-100	MV-135-3	BV-130-3	5 mm
—	MV-70-4	BV-75-4	5 mm
TE-143	MV-100-4	BV-100-4	5 mm
CV-345	MV-135-4	BV-130-4/5	5 mm
CV-351	MV-135-4	BV-130-5	6 mm
CV-351	MV-175-8	BV-175-6	8 mm
CV-351	MVF-175-8	BV-175-6	8 mm
—	MV-175-9	BV-175-8	10 mm

### Taper Point Needles

3/8 Circle

—	MVK-70-3	V-75-3	4 mm
—	MVK-100-4	V-100-3	5 mm
		V-130-3	

# Suture Needles

## Ophthalmic

### Premium Pointed Spatula Needles

D&G	USS	Ethicon	
1/4 Circle			
LO-90-8	SE-100-8	CS/TG-100-8	6 mm
3/8 Circle			
LE-140-6	SE-140-6	CS/TG-140-6	6 mm
LE-140-8	SE-140-8	CS/TG-140-8	6 mm
1/2 Circle			
LA-160-4	SE-160-4	CS/TG-160-4	6 mm
LA-160-6	SE-60-6	CS/TS-160-6	6 mm
LA-160-8	SE-160-8	CS/TG-160-8	6 mm
L-175-8	SE-175-8	CS/TG-175-8	7 mm
Compound Curve			
LB-90/50-6	SE-CC-6	S-29	5 mm

### Spatula Needles

1/4 Circle			
SLO-110-11	SS-29	S-29	7 mm
DO-5	SS-24	S-24	8 mm
—	SS-14	S-14	10 mm
1/2 Circle			
D-1	SS-2	KS	8 mm

### Premium Point Reverse Cutting Needles

3/8 Circle			
CE-30	HE-6	G-6	7 mm
CE-21	HE-1	G-1	11 mm
CE-23	HE-3	G-3	13 mm
1/2 Circle			
C-21	HE-7	G-7	7 mm
C-22	HE-7	G-7	9 mm

### Taper Point Needles

D&G	USS	Ethicon	
—	ESK	EN-S	23 mm
TS-20	EST	ST-4	19 mm
1/2 Circle			
T-31	EV-23	SH-2	17 mm
T-16	EV-25	SH-1	22 mm
T-5	EV-20	SH	26 mm
T-19	EGS-22	CT-2	27 mm

### Straight Tapercutting Needles

EKS	ST-1		24 mm
-----	------	--	-------